

R-TRAK WITH TB FLEX SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

£5,495

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 3" x 9" Flat Free Front Castors



- SRS Suspension System
- NE+ Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

15"W.	FOC	19"W.	FOC
16"W.	FOC	20"W.	FOC
17"W.	FOC	21"W.	FOC
18"W.	Standard	22"W.	FOC

SEAT PACKAGES - SEAT DEPTH

15"D.	FOC	19"D.	FOC
16"D.	FOC	20"D.	FOC
17"D.	FOC		
18"D.	Standard		

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>				Standard
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>				£450
Powered Seat Lift & Tilt. <i>110" lift & 50° tilt. Uses 2 power options.</i>				£1895
Manual Seat to Floor Height. <i>From seat pan.</i>		16.5"	17.5"	
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18"	Changes Tilt range to -5° to 25°	18.5" Changes Tilt range to -10° to 20°
Powered Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18"	18.5"	

STANDARD BACK SUPPORT POSTIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
Attendant Controlled Manual Recline. <i>Adj. range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adj. range 95° to 125°. Uses 1 power option.</i>	£345

STANDARD BACK SUPPORT

Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	Standard
TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	£135
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Enter ADI Back Support part number here for quotation.	

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour</i>				£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour</i> seat pan.				£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>				£559
Back Height	50cm. <i>Actual height.</i>	Standard	57cm. <i>Actual height.</i>	FOC

The back support width will match seat width selected in seat dimensions.

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>				Standard
TB Flex Removable, Quick Release. <i>Adj. range 9" to 14" from seat pan.</i>				£85
Single Post Quick Height Adjustable Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14"</i>				£100
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16"</i>				£279
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14"</i>					Standard	
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>					FOC	
Power Elevating Articulating Centre Mount Foot Platform. <i>Minimum seat width 14". Uses 1 power option.</i>					£700	
Power Elevating Articulating & Lower To Floor Centre Mount Foot Platform. <i>Minimum seat width 14". Lowers to the floor to aid standing transfers. Uses 2 power options.</i>					£1112	
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". Manually Adj. & powered elevating only	FOC	13"-19". Manually Adj. & powered elevating only	Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>				FOC	
	10"W x 8"D Split Footplates.		FOC	10"W x 11"D Split Footplates.		Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	}	Calf Panel. <i>18"W x 8"H.</i>		£42
	Reverse Dartex Calf Pad Covers.	£42		Reverse Dartex Calf Panel Cover.		£62
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>				£16	
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>				£16	
	Footplate Connector.		£10	Outer Edge Supports.		£10
	Wheel Bumpers.		£10	Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.		

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Adjustment range 12" to 15"</i>	FOC
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5"</i>	£330
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20"</i>	£250
Power Elevating Swing-away Legrests. <i>Adj. range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650

Enter seat pan to heel measurement.

If no measurement is supplied the standard factory length option will be fitted.

Stealth Residual Limb Support (each). *Includes a 60-90° Swing-away leg support hanger for mounting.*

LEFT

£300

Limb Support

£300

RIGHT

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M/ ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

Thoracic Lateral Supports. <i>Supplied as a pair.</i>	LEFT	£150	Lateral Thoracic Swing-away Hardware.				£150	RIGHT			
		£35	3" H x 5" L Curved Lateral Pad.				£35				
		£35	4" H x 6" L Curved Lateral Pad.				£35				
Stealth 4-Point Harness.	Extra Small	£137	}	Classic	STD*	}	Dynamic	STD*	}	Vinyl Face	FOC
	Small	£137		Contoured	FOC		Structured	FOC		Zipper Face**	FOC
	Medium	£137									
	Large	£137									
	Extra Large	£137									
Chest Strap.	Extra Small	£87	Small	£87	}	Side Release		Standard			
	Medium	£87	Large	£87		Velcro Release		£19			
	Extra Large	£87									

*STD = Standard option **Classic Only

PELVIC SOLUTIONS

Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.						Standard	
	Retractable Lap Belt, Push Button Release.						£156	
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	}	Single Pull	Standard	}	Side Release	Standard
	Small	£78		Rear Pull	FOC		Push Button	FOC
	Medium	£78		Dual Pull	FOC			
	Large	£78						
		Extra Large	£78					
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£120	Small	£120	}	Push Button	Standard	
	Medium	£120	Large	£120		Side Release	£20	

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>		Standard
NE + Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>		£225
Q-Logic 3e Controls. <i>2 power option ready.</i>		£300
Q-Logic 3e Controls. <i>5 power option ready.</i>		£500
Q-Logic 3 EX Controls. <i>No power options.</i>		£275
Q-Logic 3 EX Controls. <i>2 power option ready.</i>		£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>		£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount.</i>		£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>		£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>		£316
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>	£300	<div> <div> Push Button iAccess Toggle Lever iAccess </div> <div> If iAccess is selected choose either Push Button or Toggle lever option. </div> </div>
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						£75
Retract 4 Swing-away. <i>Full 180° Movement.</i>						£150
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
Accu-Trac 4-Pole Motors. <i>Requires Q-Logic 3 EX. Delivers enhanced traction and tracking, especially with switched drive controls. Only available with 6 mph speed program and lights package.</i>	£450
LED Front Castor and Rear Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Front Castor and Rear Fender Lights.	FOC

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Back-up Camera.	£300

NOTES & SPECIAL REQUESTS

PRICE

**R-TRAK WITH
TB FLEX SEATING**



BASE PRICE	£5,495
------------	---------------

COST OPTIONS	
--------------	--

TOTAL	
-------	--

6



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 31/08/23 V1



**Crash
Tested**



V1: 14/06/2022

