

**RETAIL**

**EPIC SEATING**

**ORGANISATION DETAILS**

Organisation Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

**DELIVERY DETAILS**


Organisation Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓞ Clients Weight \_\_\_\_\_

**CLIENT SEAT AREA MEASUREMENTS**



**MOBILITY BASE INFORMATION**

Make/Model\*.      Leggero Enzo™      Quickie® Iris®      ki Mobility Focus CR      Invacare® Solara® 3G

Frame Width\*.      Frame Depth\*.

Best when paired with foot-operated tilt-release option and non-reclining canes.

## EPIC SEATING INFORMATION

Posterior  
Pelvic  
Support  
Sizes.

### Width Options

**Small:** Fits 10"-16" frame width  
**Medium:** Fits 16"-21" frame width  
**Large:** Fits 21"-26" frame width

### Height

The posterior pelvic support is 6 inches tall.  
The pad is adjustable along the cane, allowing the pad to rest between the iliac crest and the PSIS.

## RESISTANCE LEVELS > LOW

### Low to Normal Tone

Ideal for users with limited active movement who may also experience posterior pelvic tilt. Used to allow relaxation of postural muscles while providing an articulating posterior thoracic support or to allow for manual adjustment and accommodation of limited active movement or kyphosis without the use of tools.

## HIGH < RESISTANCE LEVELS

### Low to High Tone

Ideal for users with normal to highly active movement, such as dystonic users who may exhibit an extension pattern of spasticity. Provides an articulating posterior thoracic support which can absorb energy exerted by the increased tone.

## LEGGERO ENZO

LOW

FORCE

EPiC Seating ASM Small 50N. <i>ESAT040401050.</i>	£1600
EPiC Seating REM ASM Small 50N. <i>ESAT040404050.</i>	£1600
EPiC Seating ASM Medium 50N. <i>ESAT040402050.</i>	£1600
EPiC Seating REM ASM Med 50N. <i>ESAT040405050.</i>	£1600
EPiC Seating ASM Large 50N. <i>ESAT040403050.</i>	£1600
EPiC Seating REM ASM Large 50N. <i>ESAT040406050.</i>	£1600

## HIGH

EPiC Seating ASM Small 100N. <i>ESAT040401100.</i>	£1600
EPiC Seating REM ASM Small 100N. <i>ESAT040404100.</i>	£1600
EPiC Seating ASM Medium 100N. <i>ESAT040402100.</i>	£1600
EPiC Seating REM ASM Med 100N. <i>ESAT040405100.</i>	£1600
EPiC Seating ASM Large 100N. <i>ESAT040403100.</i>	£1600
EPiC Seating REM ASM Large 100N. <i>ESAT040406100.</i>	£1600

## QUICKIE IRIS

LOW

FORCE

EPiC Seating ASM Small 50N. <i>ESAT090901050.</i>	£1600
EPiC Seating ASM Medium 50N. <i>ESAT090902050.</i>	£1600
EPiC Seating ASM Large 50N. <i>ESAT090903050.</i>	£1600

## HIGH

EPiC Seating ASM Small 100N. <i>ESAT090901100.</i>	£1600
EPiC Seating ASM Medium 100N. <i>ESAT090902100.</i>	£1600
EPiC Seating ASM Large 100N. <i>ESAT090903100.</i>	£1600

## KI FOCUS CR

LOW

FORCE

EPiC Seating ASM Small 50N. <i>ESAT030501050.</i>	£1600
EPiC Seating ASM Medium 50N. <i>ESAT030502050.</i>	£1600
EPiC Seating ASM Large 50N. <i>ESAT030503050.</i>	£1600

## HIGH

EPiC Seating ASM Small 100N. <i>ESAT030501100.</i>	£1600
EPiC Seating ASM Medium 100N. <i>ESAT030502100.</i>	£1600
EPiC Seating ASM Large 100N. <i>ESAT030503100.</i>	£1600

## INVACARE SOLARA 3G

LOW

FORCE

EPiC Seating ASM Small 50N. <i>ESAT021601050.</i>	£1600
EPiC Seating ASM Medium 50N. <i>ESAT021602050.</i>	£1600
EPiC Seating ASM Large 50N. <i>ESAT021603050.</i>	£1600

## HIGH

EPiC Seating ASM Small 100N. <i>ESAT021601100.</i>	£1600
EPiC Seating ASM Medium 100N. <i>ESAT021602100.</i>	£1600
EPiC Seating ASM Large 100N. <i>ESAT021603100.</i>	£1600

## SEAT CUSHION SELECTION PORT

Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

Omit Cushion. *For use with an aftermarket cushion.*

### TRU-COMFORT 2 SPP

	ADD WIDTH	ADD DEPTH	
TRU-Comfort 2 SPP. 14"-21"W.	ST-TC2SPP	<input type="text"/>	£165
TRU-Comfort 2 SPP. 22"+	ST-TC2SPP	<input type="text"/>	£185



Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

	Width										
	14	15	16	17	18	19	20	21	22	24	
14											
15											
16											
17											
18											
19											
20											
21											
22											
24											

### ESSENCE SPP

	ADD WIDTH	ADD DEPTH	
Essence SPP, No Wedge. 15"-20"W.	ST-ESSSPP	<input type="text"/>	£210
Essence SPP, No Wedge. 22"W.	ST-ESSSPP	<input type="text"/>	£265
Essence SPP, Small Wedge. 15"-20"W.	ST-ESSSPP	<input type="text"/>	£265
Essence SPP, Small Wedge. 22"W.	ST-ESSSPP	<input type="text"/>	£295
Essence SPP, Med Wedge. 15"-20"W.	ST-ESSSPP	<input type="text"/>	£265
Essence SPP, Medium Wedge. 22"W.	ST-ESSSPP	<input type="text"/>	£295
Essence SPP, Large Wedge. 15"-20"W.	ST-ESSSPP	<input type="text"/>	£265
Essence SPP, Large Wedge. 22"W.	ST-ESSSPP	<input type="text"/>	£295




Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

	Width										
	10	12	14	15	16	17	18	19	20	21	22
10											
12											
14											
15											
16											
17											
18											
19											
20											
21											
22											
24											

### SOLUTION SPP

	ADD WIDTH	ADD DEPTH	
Solution SPP. 10"-20"W.	ST-SOLSPP	<input type="text"/>	£190
Solution SPP. 22"-24"W	ST-SOLSPP	<input type="text"/>	£210



Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

	Width										
	10	12	14	15	16	17	18	19	20	22	24
10											
12											
14											
15											
16											
17											
18											
19											
20											
22											
24											

## SPECTRUM FOAM SPP

Spectrum Foam SPP. 16"-20"W.	ST-SPCFSP	<input type="text"/>	<input type="text"/>	£265
Spectrum Foam SPP. 22"-24"W	ST-SPCFSP	<input type="text"/>	<input type="text"/>	£310
Spectrum Foam SPP. 26"-32"W	ST-SPCFSP	<input type="text"/>	<input type="text"/>	£310




Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

	Width											
	14	16	17	18	19	20	22	24	26	28	30	32
15												
16												
17												
18												
19												
20												
21												
22												
24												
26												
28												

## SPECTRUM GEL SPP

Spectrum Gel SPP. 16"-20"W.	ST-SPCGSPP	<input type="text"/>	<input type="text"/>	£265
Spectrum Gel SPP. 22"-24"W	ST-SPCGSPP	<input type="text"/>	<input type="text"/>	£310




Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

	Width							
	14	16	17	18	19	20	22	24
15								
16								
17								
18								
19								
20								
21								
22								
24								
26								

## GLACIAL SP

Glacial SP. 10"-20"W.	ST-GLSP	<input type="text"/>	<input type="text"/>	£190
Glacial SP. 22"-24"W	ST-GLSP	<input type="text"/>	<input type="text"/>	£210



Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

	Width											
	10	12	14	15	16	17	18	19	20	21	22	24
10												
12												
14												
15												
16												
17												
18												
19												
20												
21												
22												
23												
24												

## ZEN SP

Zen SP. 10"-21"W.	ST-ZENSP <input type="text"/> <input type="text"/>	£150
Zen SP. 22"-24"W	ST-ZENSP <input type="text"/> <input type="text"/>	£170



Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

		Width											
		10	12	14	15	16	17	18	19	20	21	22	24
Depth	10												
	12												
	14												
	15												
	16												
	17												
	18												
	19												
	20												
	21												
	22												
	23												
	24												

## PREMIERE P

Premiere P. 14"-20"W.	ST-PREMP <input type="text"/> <input type="text"/>	£210
Premiere P. 22"W.	ST-PREMP <input type="text"/> <input type="text"/>	£230



Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

		Width							
		14	15	16	17	18	19	20	22
Depth	13								
	14								
	15								
	16								
	17								
	18								
	19								
	20								
22									

## SIMPLICITY G

Zen SP. 10"-20"W.	ST-SIMG <input type="text"/> <input type="text"/>	£75
Zen SP. 22"-24"W	ST-SIMG <input type="text"/> <input type="text"/>	£85



Select your size on the table to the right, then complete the Part Number by adding the Width first then Depth.

		Width											
		10	12	14	15	16	17	18	19	20	21	22	24
Depth	10												
	12												
	14												
	15												
	16												
	17												
	18												
	19												
	20												
	21												
	22												
	23												
	24												

## BACK SUPPORT

Choosing the right height back support: posterior support is divided into two segments, Pelvic and Thoracic. The 6" pelvic support is a standard component, please factor in a 2" gap between the two segments when selecting your thoracic back height.

Omit Backrest. <i>For use with an aftermarket backrest.</i>		FOC?
<b>Back Support Width.</b>		
ADI Active Back Support.	ADI AL, Low Height, 2Pt Pro Hardware 10".	£520
	ADI AL, Medium Height, 2Pt Pro Hardware 13".	£520
	ADI AL, Low Height, 4Pt Pro Hardware 10".	£520
	ADI AL, Medium Height, 4Pt Pro Hardware 13".	£520
	ADI AL, Low Height, Width Adjustable Elite Hardware 11".	£650
	ADI AL, Medium Height, Width Adjustable Elite Hardware 13".	£650
ADI Deep Back Support.	ADI AL, Low Height, 2Pt Pro Hardware 10".	£650
	ADI AL, Medium Height, 2Pt Pro Hardware 13".	£650
	ADI AL, Low Height, 4Pt Pro Hardware 10".	£650
	ADI AL, Medium Height, 4Pt Pro Hardware 13".	£650
	ADI AL, Low Height, Width Adjustable Elite Hardware 11".	£725
	ADI AL, Medium Height, Width Adjustable Elite Hardware 13".	£725
	ADI AFT, Low Height, 2Pt Pro Hardware 10".	£650
	ADI AFT, Medium Height, 2Pt Pro Hardware 13".	£650
	ADI AFT, Low Height, 4Pt Pro Hardware 10".	£650
	ADI AFT, Medium Height, 4Pt Pro Hardware 13".	£650
	ADI AFT, Low Height, Width Adjustable Elite Hardware 11".	£725
	ADI AFT, Medium Height, Width Adjustable Elite Hardware 13".	£725
ADI Extra Deep Back Support.	ADI AL, Medium Height, 2Pt Pro Hardware 13".	£650
	ADI AL, Medium Height, 4Pt Pro Hardware 13".	£650
	ADI AL, Medium Height, Width Adjustable Elite Hardware 13".	£725
Custom Back Selection.	Custom Back Selection. <i>E2617. See TFB Flex, TFB Contour or TFB Linear Order Forms.</i>	

## THORACIC LATERALS

Curved Lateral Pads.	UniLink 3x5 Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP535050320.</i>	£40
	UniLink 4x6 Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP546050320.</i>	£40
	UniLink 5x7 Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP557050320.</i>	£40
Curved Lateral Pads Gel.	UniLink 3x5 Gel Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP535057120.</i>	£66
	UniLink 4x6 Gel Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP546057120.</i>	£66
	UniLink 5x7 Gel Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP557057120.</i>	£66
Flat Lateral Pads.	UniLink 3x5 Flat Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP835050320.</i>	£40
	UniLink 3x5 Flat Lateral Pad, CoolCore, Metric, 1" Spacing. <i>Part: ULP835060320.</i>	£40
	UniLink 3x5 Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing. <i>Part: ULP835080320.</i>	£40
	UniLink 4x6 Flat Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP846050320.</i>	£40
	UniLink 4x6 Flat Lateral Pad, CoolCore, Metric, 1" Spacing. <i>Part: ULP846060320.</i>	£40
	UniLink 4x6 Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing. <i>Part: ULP846080320.</i>	£40
	UniLink 5x7 Flat Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP857050320.</i>	£40
	UniLink 5x7 Flat Lateral Pad, CoolCore, Metric, 1" Spacing. <i>Part: ULP857060320.</i>	£40
Flat Lateral Pads Gel.	UniLink 3x5 Gel Flat Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP835057120.</i>	£66
	UniLink 3x5 Gel Flat Lateral Pad, CoolCore, Metric, 1" Spacing. <i>Part: ULP835067120.</i>	£66
	UniLink 3x5 Gel Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing. <i>Part: ULP835087120.</i>	£66
	UniLink 4x6 Gel Flat Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP846057120.</i>	£66
	UniLink 4x6 Gel Flat Lateral Pad, CoolCore, Metric, 1" Spacing. <i>Part: ULP846067120.</i>	£66
	UniLink 4x6 Gel Flat Lateral Pad, Stealth Tek, Metric, 1" Spacing. <i>Part: ULP846087120.</i>	£66
	UniLink 5x7 Gel Flat Lateral Pad, Neoprene, Metric, 1" Spacing. <i>Part: ULP857057120.</i>	£66
	UniLink 5x7 Gel Flat Lateral Pad, CoolCore, Metric, 1" Spacing. <i>Part: ULP857067120.</i>	£66
UniLink Hardware.	<b>LEFT</b> £175 <i>ULS25100-L.</i> UniLink Swing-Away Lateral. <i>ULS25100-R.</i> £175 <b>RIGHT</b>	
	£175 <i>ULS25000-L.</i> UniLink Swing-Away Lateral Short. <i>ULS25000-R.</i> £175	
UniLink Mount Options.	1 inch Cane Clamp. <i>Part: UIMCNI100.</i>	£36
	1 inch Crescent Clamp. <i>Part: UIMCRI100.</i>	£36
	90 Degree Lateral Bracket Mount (Planar, TFB I, Flex Mod.). <i>Part: UIMBK1090.</i>	£36
	105 Degree Lateral Bracket Mount (TFB Rectangle). <i>Part: UIMBK1105.</i>	£36
	135 Degree Lateral Bracket Mount (TFB Contour, TFB Flex Active, ADI, Q-Back). <i>Part: UIMBK1135.</i>	£36
TWB Thin Series Hardware.	<b>LEFT</b> £179 Lateral Thoracic Thin Swing-Away Hardware <i>Part: TWSTL-L</i> £179 <b>RIGHT</b>	
	£179 Lateral Thoracic Thin Swing-Away Hardware (No Link) <i>Part: TWBTL-L</i> £179 <i>Part: TWBTL-R</i>	
Thin Series Mount Options.	Lateral Thoracic Cane Mount, 1 inch. <i>Part: LTHW500.</i>	£36



## HEAD SUPPORT

Headrest Mounting Hardware.	Standard Straight Mounting Hardware, 17mm Ball, Metric. <i>Part: HMO400-17M.</i>		£86	
	Multi-Axis Hardware, 17mm Ball, Metric. <i>Part: HMO475-17M.</i>		£86	
	UniLink Headrest HW, 17mm Ball, 1/2" SQ 7"VB. <i>Part: UHB2D33200 HCPCS: E1028.</i>		£111	
	UniLink Headrest HW, 17mm Ball, 15MM SQ 12"VB. <i>Part: UHB3D33200 HCPCS: E1028.</i>		£111	
	UniLink Headrest HW, 17mm Ball, 1/2" SQ 12"VB. <i>Part: UHB4D33200 HCPCS: E1028.</i>		£111	
Headrest Hardware Mount.	Universal Headrest Adapter Plate. <i>Part: SBM300U.</i>		£35	
Headrest Tone Deflector.	Headrest Tone Deflector, for 17MM Ball, Metric. <i>Part: TD100M HCPCS: K0108.</i>		£35	
Comfort Plus Nino Head Support Pads.	6" Headrest, Flat Plate, Metric. <i>Part: CP06CCFPM.</i>	£108	8" Headrest, Flat Plate, Metric. <i>Part: CP08CCFPM.</i>	£108
	10" Headrest, Flat Plate, Metric. <i>Part: CP10CCFPM.</i>	£108	14" Headrest, Flat Plate, Metric. <i>Part: CP14CCFPM.</i>	£108
	19" Headrest, Flat Plate, Metric. <i>Part: CP19CCFPM.</i>		£108	
Adjustable Comfort Plus Head Support Pads.	Adjustable Comfort Plus 14" Headrest, Metric. <i>Part: ACP150M.</i>		?	
	Adjustable Comfort Plus 19" Headrest, Metric. <i>Part: ACP550M.</i>		£268	

## EPIC SEATING REPLACEMENT PARTS

Replacement Posterior Pelvic Block.	10 wide, 6in Tall. <i>Part: ESR0008010.</i>	£122	11 wide, 6in Tall. <i>Part: ESR0008011.</i>	£122	12 wide, 6in Tall. <i>Part: ESR0008012.</i>	£122
	13 wide, 6in Tall. <i>Part: ESR0008013.</i>	£122	14 wide, 6in Tall. <i>Part: ESR0008014.</i>	£122	15 wide, 6in Tall. <i>Part: ESR0008015.</i>	£122
	16 wide, 6in Tall. <i>Part: ESR0008016.</i>	£122	17 wide, 6in Tall. <i>Part: ESR0008017.</i>	£122	18 wide, 6in Tall. <i>Part: ESR0008018.</i>	£122
	19 wide, 6in Tall. <i>Part: ESR0008019.</i>	£122	20 wide, 6in Tall. <i>Part: ESR0008020.</i>	£122	21 wide, 6in Tall. <i>Part: ESR0008021.</i>	£122
	22 wide, 6in Tall. <i>Part: ESR0008022.</i>			£122		
Posterior Pelvic Support.	Lateral Pelvic Wedge Kit. <i>Part: ESR0006000.</i>					£79
Replacement Struts.	Single Gas Strut, 50N. <i>Part: ESR0003050.</i>				£130	x Qty
	Single Gas Strut, 100N. <i>Part: ESR0003100.</i>				£130	x Qty
	Cable Assembly. <i>Part: ESR0004000.</i>				£60	x Qty
Pivot Parts.	Replacement Release Pin. <i>Part: ESR0001000.</i>				£31	x Qty
	Replacement Bearings. <i>Part: ESR0002000.</i>				£69	x Qty
Attachment Tube (pair).	Invacare Solara 3G, Attachment Tube. <i>Part: ESR020000.</i>					£47
	Ki Focus CR, Attachment Tube. <i>Part: ESR0300000.</i>					£47
	Leggero Enzo, Attachment Tube. <i>Part: ESR0400000.</i>					£47
	Quickie Iris, Attachment Tube. <i>Part: ESR0900000.</i>					£47
Foam Grip Handles.	Replacement Foam Grip Handles. <i>Part: ESR0005000.</i>				£22	x Qty

## NOTES & SPECIAL REQUESTS

## PRICE



**Stealth Products**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**EPIC SEATING -20/11/23 V1**



V1: 14/06/2022