

MOTABILITY ID 714948

## ASPEN WITH TB FLEX SEAT / POWER TILT / 6 MPH

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Clients Weight \_\_\_\_\_

### BASE MODEL

- 4 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 3" x 9" Flat Free Front Castors
- ATX Suspension System



- Powered Tilt
- NE+ Control System
- 2 x 62 ah Batteries
- 6 Mph Motors
- Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

### FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
3" x 9" Pneumatic Black Tyres.	FOC

### COLOURS

Back in Black	Standard
Electric Blue	FOC
Silver Bullet	FOC
Red Corvette	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

## SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified.

## SEAT POSITIONING

Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	Standard
Seat to Floor Height. <i>From seat pan.</i>	17.5"      18" Changes tilt range to -5° to 25°      18.5" Changes tilt range to -5° to 25°

## STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°.</i>	Standard
---	----------

## STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
---	----------

## ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16"</i>	FOC
Arm Support Pads	10" Standard      Standard      10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	14" Standard      FOC      14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC

## LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard	
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	FOC	
Manual Elevating Swing-away Legrests. <i>Adjustment range 15" to 20".</i>	FOC	
Seat pan to heel measurement	(Leave blank for standard factory)	Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.

Standard (STD)

## NON-EXPANDABLE CONTROLS

NE + Controls. *2 power option ready.*

Standard

Control Method. Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

## HAND CONTROL MOUNTS

Fixed.

Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.*

FOC

Mounting Position.      Left                      FOC                      Right                      Standard

## SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.

Standard

## ACCESSORIES

Kerb Climber.

FOC

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**MOTABILITY ID 714948 - 31/08/23 V1**



**Crash Tested**



V1: 14/06/2022

