

EDGE 3 STRETTO - CONTOURED SEAT / STANDARD

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight (required for suspension settings) _____

BASE MODEL

- 14 Colour Options
- Flip-up Armrests
- 70° Swing-Away Legrests
- Contoured Seat with Fold-Forward Backrest
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



CONTOURED SEAT PACKAGES

16"W x 16"D	FOC	18"W x 18"D	Standard	20"W x 20"D	FOC
-------------	-----	-------------	----------	-------------	-----

Seat package includes standard headrest and auto buckle lap belt.

SEAT CUSHIONS

Contoured Seat Cushion. *136 Kg weight limit.*

Standard

Seat cushion dimensions will match seat package selected above.

SEAT POSITIONING

Manual Settable Tilt. *0° to 30° tool adjustable with a telescopic tube.*

Standard

Manual Seat to Floor Height. *From seat pan.*

16.5"

17.5"

18.5"

ARM SUPPORTS

Std Flip-up Height Adjustable. *Includes contour arm pads. Do not select arm pads below.*

Standard

LEG SUPPORTS

70° Swing-away Legrests. *Adjustment range 12" to 15".*

Standard

NON-EXPANDABLE CONTROLS

NE + Controls. *With no power options.*

Standard

VR2 Controls. *With no power options.*

FOC

HAND CONTROL MOUNTS

Fixed.

Standard

Mounting Position.

Left

FOC

Right

Standard

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. *Lights and indicators not included.*

Standard

NOTES & SPECIAL REQUESTS

2



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
Motability ID 706811 - 17/03/23 V1



V1: 14/06/2022



Crash Tested

