

MOTABILITY ID 707223

EDGE 3 STRETTO TB FLEX SEAT / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name _____
 Order Date _____
 Address _____
 _____ Postcode _____
 Phone _____ Email _____
 PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____
 _____ Postcode _____
 Quantum Specialist _____
 Ⓢ Clients Weight (required for suspension settings) _____

BASE MODEL

- 17 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors
- 2" x 6" Flat Free Rear Castors



- SRS Suspension System
- Powered Tilt
- Powered Recline
- NE+ Control System
- 2 x 40 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Military Green	FOC
Night Train	FOC	Lemon Crush	FOC
Midnight Blue	FOC	Cherry Bomb	FOC
Electric Blue	FOC	Evolution Orange	FOC
Blue Jean	FOC	White Satin	FOC
Deep Purple	FOC	Grey Street	FOC
Sugar Plum	FOC	Green Machine	FOC
Red Corvette	FOC	Raspberry Beret	FOC
Pretty in Pink	FOC		

Please download our free PQ-365 app to view colours and other product features in more detail.



STANDARD SEAT DIMENSIONS

15" Width	FOC	16" Width	FOC	15" Depth	FOC	16" Depth	FOC
17" Width	FOC	18" Width	Standard	17" Depth	FOC	18" Depth	Standard
19" Width	FOC	20" Width	FOC	19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	Standard
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5" 18.5" 19.5"

STANDARD BACK SUPPORT POSITIONING

Powered Back Recline. <i>Adj. range 95° to 125°. Uses 1 power option.</i>	Standard
---	----------

STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
---	----------

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard	
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC	
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i>	FOC	
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC	
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>	FOC	
Arm Support Pads	10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
90° Paediatric Swing-away Legrests. <i>Adjustment range 7" to 11".</i>	FOC
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	FOC
Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 10" to 16". 5"W x 11"L Footplates.</i>	FOC
Seat pan to heel measurement	(Leave blank for standard factory)

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M/CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M/CP14CCFPM.</i>	FOC

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

NON-EXPANDABLE CONTROLS

NE + Controls. <i>2 power option ready.</i>	Standard	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

HAND CONTROL MOUNTS

Fixed.	Standard				
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>	FOC				
Mounting Position.	<table border="0"> <tr> <td>Left</td> <td>FOC</td> <td>Right</td> <td>Standard</td> </tr> </table>	Left	FOC	Right	Standard
Left	FOC	Right	Standard		

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
---	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 707223 - 18/01/23 V1



Crash Tested



V1: 14/06/2022

