

MOTABILITY ID 707382

EDGE 3 STRETTO WITH TB3 SEATING / iLEVEL / POWER TILT

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight (required for suspension settings) _____

BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- iLevel Seat Lift and Powered Tilt
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-16"W
14"-16"W
16"-20"W
Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-16"D
14"-16"D
16"-20"D
Set depth at:

Minimum seat size of 14"W x 14"D required for manually adjustable and powered centre mount foot platform when iLevel is selected.

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

iLevel Powered Seat Lift & Tilt. <i>12" lift & 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options.</i>	Standard
Powered iLevel Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	18" 19" 20"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a pair of plastic lateral wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Lateral Wedges. <i>Select size in Back Support Dimensions</i>	Standard
---	----------

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports or powered recline.</i>	FOC

BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Lateral Wedges. <i>Select size in Back Support Dimensions</i>	Standard
Back Cane Only Support Positioning. Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
Back Cane Only Support Positioning. Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>	FOC
Back Cane Height.	16"H 17"H 18"H 19"H 20"H
<i>From seat pan.</i>	21"H 22"H 23"H 24"H 25"H

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Ignore Back Support Dimensions if omit back support is selected in TB3 or Back Canes sections. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.		12"-14"W	14"-16"W	16"-20"W
		18"-22"W	20"-24"W	
Back Support Width.	12"W	14"W	16"W	17"W
	18"W	19"W	20"W	21"W
	22"W	24"W		
Back Support Height. <i>From seat pan.</i>	14"H. <i>12"W back only.</i>	15"H. <i>12"W back only.</i>	16"H.	17"H.
	18"H.	19"H.	20"H.	21"H.
	22"H.	23"H.	24"H.	25"H. <i>Not available with 24"W back.</i>
	26"H. <i>Not available with 24"W back.</i>	27"H. <i>Not available with 24"W back.</i>	28"H. <i>Not available with 24"W back.</i>	29"H. <i>Not available with 24"W back.</i>

19" minimum backrest height when powered recline is selected.

ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				Standard
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>				FOC
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>				FOC
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				FOC
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16"</i>	FOC		

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>				Standard
Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>				FOC
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". <i>Manually adjustable & powered elevating only</i>	Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>	FOC	10"W x 8"D Split Footplates.	FOC
	10"W x 11"D Split Footplates.	Standard		
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard		

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-away Legrests. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>				FOC
90° Paediatric Swing-away Legrests. <i>4"W x 6"D angle & depth adj. footplates. Adjustment range 7" to 11".</i>				FOC
Manual Elevating Swing-Away Legrests. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>				FOC

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-away Legrests. <i>Select lowers below. Supplied with 8.5"W x 6"D extruded footplates. Minimum seat width 15"</i>						FOC
60° to 90° Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	13.5"-16.5"	Standard	15.5"-19.5"	FOC

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>2 power option ready.</i>	Standard						
<table border="0"> <tr> <td rowspan="2"> iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified. Supplied FOC when ilevel is selected.</i> </td> <td rowspan="2">STD</td> <td rowspan="2">}</td> <td>Push Button iAccess</td> <td rowspan="2">If iAccess is selected choose either Push Button or Toggle lever option.</td> </tr> <tr> <td>Toggle Lever iAccess</td> </tr> </table>	iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified. Supplied FOC when ilevel is selected.</i>	STD	}	Push Button iAccess	If iAccess is selected choose either Push Button or Toggle lever option.	Toggle Lever iAccess	
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified. Supplied FOC when ilevel is selected.</i>				STD		}	Push Button iAccess
	Toggle Lever iAccess						
Omit iAccess OLED 4 Button Switch Module.							
<table border="0"> <tr> <td rowspan="2">Control Method.</td> <td>Hand Control. <i>Select position in Hand Control Mounts below.</i></td> <td>Standard</td> </tr> <tr> <td>Attendant Control. <i>Hand control mounted at rear.</i></td> <td>FOC</td> </tr> </table>	Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard	Attendant Control. <i>Hand control mounted at rear.</i>	FOC		
Control Method.		Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard				
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC					

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.	Standard				
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>	FOC				
Mounting Position.	<table border="0"> <tr> <td>Left</td> <td>FOC</td> <td>Right</td> <td>Standard</td> </tr> </table>	Left	FOC	Right	Standard
Left	FOC	Right	Standard		

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	Standard
6 Mph Program with Lights and indicators.	FOC
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 707382 - 17/03/23 V1



Crash Tested



V1: 14/06/2022

