

EDGE 3 STRETTO WITH TB3 SEAT / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight (required for suspension settings) _____

BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- Powered Tilt and Powered Recline
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|--|----------|
| 3" x 12.5" Flat Free Black Tyres. | Standard |
| 3" x 12.5" Pneumatic Black Tyres. | FOC |
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | FOC |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC |

COLOURS

| | | | |
|---------------|----------|-------------------|-----|
| Back in Black | Standard | Pretty in Pink | FOC |
| Green Machine | FOC | Lemon Crush* | FOC |
| Midnight Blue | FOC | Night Train* | FOC |
| Electric Blue | FOC | Evolution Orange* | FOC |
| Blue Jean | FOC | Military Green* | FOC |
| Deep Purple | FOC | Grey Street* | FOC |
| Red Corvette | FOC | Raspberry Beret* | FOC |

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-16"W

14"-16"W

16"-20"W

Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-16"D

14"-16"D

16"-20"D

Set depth at:

Minimum seat size of 14"W x 14"D required for manually adjustable and powered centre mount foot platform when iLevel is selected.

STEALTH SEAT CUSHIONS

No Seat Cushion Required

FOC

Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. *158.7 Kg weight limit.*

Standard

Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. *160 Kg weight limit.*

FOC

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.*

FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. *0° to 50° COG tilt. Uses 1 power option.*

Standard

Powered Tilt Seat to Floor Height. *From seat pan.*

17.5"

18.5"

19.5"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. *Select TB3 back options below. Ignore all 'Back Cane' options.*

Standard

Includes a metal active contour back support shell, TRU-Comfort 2 back cushion, and a pair of plastic lateral wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Lateral Wedges.

Standard

Select size in Back Support Dimensions

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Powered Recline. *175° with sheer reduction. Uses 1 power option.*

Standard

BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB3 Back Support System.

| Backrest Support Width Package. | | 12"-14"W. | 14"-16"W. | 16"-20"W. |
|--|-------|-----------|-----------|-----------|
| | | 18"-22"W. | 20"-24"W. | |
| Back Support Width. | 12"W. | 14"W. | 16"W. | 17"W. |
| | 18"W. | 19"W. | 20"W. | 21"W. |
| | 22"W. | 24"W. | | |
| Back Support Height. <i>From seat pan.</i> | 19"H. | 20"H. | 21"H. | 22"H. |
| | 23"H. | 24"H. | 25"H. | 26"H. |
| | 27"H. | 28"H. | 29"H. | |

ARM SUPPORTS

| | | | | |
|---|---|----------|--|----------|
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i> | | | | Standard |
| Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i> | | | | FOC |
| Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i> | | | | FOC |
| TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i> | | | | FOC |
| Arm Support Pads | 10" Standard | Standard | 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC |
| | 14" Standard | FOC | 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC |
| | 17" Standard. <i>TB3 only. Minimum seat depth 16"</i> | FOC | | |

CENTRE MOUNTED LEG SUPPORTS

| | | | | |
|--|---|----------|--|----------|
| Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i> | | | | Standard |
| Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i> | | | | FOC |
| Lower Leg Extensions. <i>From seat pan.</i> | 8"-12". <i>Paediatric centre mount only.</i> | FOC | 10"-16". <i>Manually adjustable & powered elevating only</i> | Standard |
| Footplates. | 10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i> | FOC | 10"W x 8"D Split Footplates. | FOC |
| | 10"W x 11"D Split Footplates. | Standard | | |
| Calf Support. | Individual Calf Pads. <i>6"W x 6"H.</i> | Standard | | |

INDIVIDUALLY MOUNTED LEG SUPPORTS

| | | | | |
|---|--|--|--|-----|
| 70° Swing-away Legrests. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i> | | | | FOC |
| 90° Paediatric Swing-away Legrests. <i>4"W x 6"D angle & depth adj, footplates. Adjustment range 7" to 11".</i> | | | | FOC |
| Manual Elevating Swing-Away Legrests. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i> | | | | FOC |

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

| | | | | | | |
|---|---------|-----|-------------|----------|-------------|-----|
| 60° to 90° Tool Adjustable Swing-away Legrests. <i>Select lowers below. Supplied with 8.5"W x 6"D extruded footplates. Minimum seat width 15"</i> | | | | | FOC | |
| 60° to 90° Lower Leg Extensions. <i>From seat pan.</i> | 11"-13" | FOC | 13.5"-16.5" | Standard | 15.5"-19.5" | FOC |

STEALTH HEAD SUPPORTS

| | | | | |
|--|--|--|--|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | | | | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | | | | FOC |

PELVIC BELTS

| | | | | |
|--|--|--|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | | | | Standard |
|--|--|--|--|----------|

EXPANDABLE CONTROLS

| | | | | |
|---|--|--|--|----------|
| Q-Logic 3e Controls. <i>2 power option ready.</i> | | | | Standard |
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i> | | | Standard |
| | Attendant Control. <i>Hand control mounted at rear.</i> | | | FOC |

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

| | | | | |
|---|------|-----|-------|----------|
| Fixed. | | | | Standard |
| Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i> | | | | FOC |
| Mounting Position. | Left | FOC | Right | Standard |

SPEED PROGRAM & LIGHTS PACKAGE

| | |
|---|----------|
| 4 Mph Program with Lights and indicators. | Standard |
| 6 Mph Program with Lights and indicators. | FOC |
| LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i> | Standard |
| Omit LED Fender Lights. | FOC |

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 707383 - 18/04/23 V1



Crash Tested



V1: 14/06/2022

