

## EDGE 3 STRETTO WITH TB3 SEAT / iLEVEL / POWER TILT / POWER RECLINE

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓢ Clients Weight (required for suspension settings) \_\_\_\_\_

### BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- iLevel Seat Lift
- Powered Tilt
- Powered Recline
- 12.5" Flat Free Drive Tyres



- 2" x 5" Flat Free Front Castors
- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

12"-16"W
14"-16"W
16"-20"W
Set width at:

## SEAT PACKAGES - SEAT DEPTH

12"-16"D
14"-16"D
16"-20"D
Set depth at:

Minimum seat size of 14"W x 14"D required for manually adjustable and powered centre mount foot platform when iLevel is selected.

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

iLevel Powered Seat Lift & Tilt. <i>12" lift &amp; 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options.</i>	Standard
Powered iLevel Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	18"      19"      20"

## BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below.</i>	Standard
--	----------

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Lateral Wedges. <i>Select size in Back Support Dimensions</i>	Standard
---	----------

## TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	Standard
---	----------

## BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB3 Back Support System.

Backrest Support Width Package.		12"-14"W	14"-16"W	16"-20"W
		18"-22"W	20"-24"W	
Back Support Width.	12"W	14"W	16"W	17"W
	18 "W	19"W	20"W	21"W
	22"W	24"W		
Back Support Height.	19"H.	20"H.	21"H.	22"H.
	23"H.	24"H.	25"H. <i>Not available with 24"W back.</i>	26"H. <i>Not available with 24"W back.</i>
	<i>From seat pan.</i>	27"H. <i>Not available with 24"W back.</i>	28"H. <i>Not available with 24"W back.</i>	29"H. <i>Not available with 24"W back.</i>

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16"</i>				Standard
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14"</i>				FOC
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>				FOC
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				FOC
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16"</i>	FOC		

## CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14"</i>				Standard
Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>				FOC
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". <i>Manually adjustable &amp; powered elevating only</i>	Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>	FOC	10"W x 8"D Split Footplates.	FOC
	10"W x 11"D Split Footplates.	Standard		
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard		

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-away Legrests. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15"</i>				FOC
90° Paediatric Swing-away Legrests. <i>4"W x 6"D angle &amp; depth adj. footplates. Adjustment range 7" to 11"</i>				FOC
Manual Elevating Swing-Away Legrests. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15"</i>				FOC

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-away Legrests. <i>Select lowers below. Supplied with 8.5"W x 6"D extruded footplates. Minimum seat width 15"</i>						FOC
60° to 90° Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	13.5"-16.5"	Standard	15.5"-19.5"	FOC

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M/CP10CCFPM.</i>				Standard
UniLink 14" Comfort Plus. <i>HMO475-17M/CP14CCFPM.</i>				FOC

## PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.				Standard
--	--	--	--	----------

## EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>2 power option ready.</i>		Standard			
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified. Supplied FOC when ilevel is selected.</i>	STD	<table border="0"> <tr> <td>Push Button iAccess</td> <td rowspan="2">} If iAccess is selected choose either Push Button or Toggle lever option.</td> </tr> <tr> <td>Toggle Lever iAccess</td> </tr> </table>	Push Button iAccess	} If iAccess is selected choose either Push Button or Toggle lever option.	Toggle Lever iAccess
Push Button iAccess	} If iAccess is selected choose either Push Button or Toggle lever option.				
Toggle Lever iAccess					
Omit iAccess OLED 4 Button Switch Module.					
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard			
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC			

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.	Standard				
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>	FOC				
Mounting Position.	<table border="0"> <tr> <td>Left</td> <td>FOC</td> <td>Right</td> <td>Standard</td> </tr> </table>	Left	FOC	Right	Standard
Left	FOC	Right	Standard		

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	Standard
6 Mph Program with Lights and indicators.	FOC
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
 Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**Motability ID 706909 - 17/03/23 V1**



**Crash Tested**

