

FUSION E ACTIVE CONTOUR 4 MPH

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Adjustable Seat Width & Depth Package
- TRU-Comfort 2 Seat Cushion
- Stealth Headrest
- Fold-forward Back with Stealth Pro Back Support
- 14" Solid Drive Wheel
- 10" Solid Castors



- Rear Suspension System
- Black Frame
- Yellow Wheel Accent Inserts
- VR2 Control System
- 2 x 40 ah Batteries
- Max User Weight:
23 st 6lbs. (150 kg)
- Stealth Pro Back Support

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

SEAT PACKAGE

18"W x 20"D. *Seat width can be adjusted between 16-23". Seat depth can be adjusted between 16-21".* Standard

SEAT CUSHIONS

2" Foam Seat Cushions. *Supplied with 1 x 16"W x 18"D & 1 x 20" x 22".* Standard

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.* Standard

Enter TC2 cushion size here **WIDTH** **DEPTH** *Refer to cushion order forms for size availability.*

BACK SUPPORT POSITIONING

Angle Adjustable, Fold-Forward. *3° - 33° angle adjustment in 4 steps of 10°.* Standard

STEALTH PRO BACK SUPPORTS

| | | | | |
|-----------------------------------------------------------|------------------------------|----------|------------------------------|----------|
| Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i> | | | | Standard |
| Back Height | 50 cm. <i>Actual height.</i> | Standard | 57 cm. <i>Actual height.</i> | FOC |
| Back Width | 16"W. | 16.5"W. | 17"W. | |
| | 18"W. | 19"W. | 20"W. | |
| | 21"W. | 22"W. | | |

SEAT POSITIONING

| | |
|------------------------------------------------------------------------------------------------|----------|
| Static Seat. <i>18.5" front seat to floor and 17" rear seat to floor height from seat pan.</i> | Standard |
|------------------------------------------------------------------------------------------------|----------|

ARM SUPPORTS

| | |
|-------------------------------------------------------------------------------------------------------------------|----------|
| Height Adjustable, Removable. <i>Includes 15"L x 3"W arm pads. Adjustment range 10.5" to 14.5" from seat pan.</i> | Standard |
|-------------------------------------------------------------------------------------------------------------------|----------|

LEG SUPPORTS

| | |
|----------------------------------------------------------------------------|----------|
| 70° Swing-away Legrests. <i>Adjustment range 13" to 16" from seat pan.</i> | Standard |
|----------------------------------------------------------------------------|----------|

PELVIC SOLUTIONS

| | |
|--------------------------------------------|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--------------------------------------------|----------|

NON-EXPANDABLE CONTROLS

| | |
|-------------------|----------|
| VR2 Hand Control. | Standard |
|-------------------|----------|

HAND CONTROL MOUNTS

| | |
|-------------------------------------------------------------------------|----------|
| Right Side with Fixed Position Mount. <i>Can be moved to left side.</i> | Standard |
|-------------------------------------------------------------------------|----------|

STEALTH HEAD SUPPORTS

| | |
|----------------------------------------------------------|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | FOC |

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 706356 - 25/11/22 V1



Crash Tested

