

## FUSION E ACTIVE CONTOUR 6 MPH

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓢ Clients Weight \_\_\_\_\_

### BASE MODEL

- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Adjustable Seat Width & Depth Package
- TRU-Comfort 2 Seat Cushion
- Stealth Headrest
- Fold-forward Back with Stealth Pro Back Support
- 14" Pneumatic Drive Wheel
- 10" Pneumatic Castors



- Rear Suspension System
- Black Frame
- Yellow Wheel
- Accent Inserts
- VR2 Control System
- 2 x 40 ah Batteries
- Max User Weight:  
23 st 6lbs. (150 kg)
- Stealth Pro Back Support

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### SEAT PACKAGE

18"W x 20"D. *Seat width can be adjusted between 16-23". Seat depth can be adjusted between 16-21".* Standard

### SEAT CUSHIONS

2" Foam Seat Cushions. *Supplied with 1 x 16"W x 18"D & 1 x 20" x 22".* Standard

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.* Standard

Enter TC2 cushion size here

**WIDTH**

**DEPTH**

*Refer to cushion order forms for size availability.*

### BACK SUPPORT POSITIONING

Angle Adjustable, Fold-Forward. *3° - 33° angle adjustment in 4 steps of 10°.* Standard

## STEALTH PRO BACK SUPPORTS

Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>				Standard
Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC
Back Width	16"W.	16.5"W.	17"W.	
	18"W.	19"W.	20"W.	
	21"W.	22"W.		

## SEAT POSITIONING

Static Seat. <i>18.5" front seat to floor and 17" rear seat to floor height from seat pan.</i>	Standard
--	----------

## ARM SUPPORTS

Height Adjustable, Removable. <i>Includes 15"L x 3"W arm pads. Adjustment range 10.5" to 14.5" from seat pan.</i>	Standard
---	----------

## LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 13" to 16" from seat pan.</i>	Standard
--	----------

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

## NON-EXPANDABLE CONTROLS

VR2 Hand Control.	Standard
-------------------	----------

## HAND CONTROL MOUNTS

Right Side with Fixed Position Mount. <i>Can be moved to left side.</i>	Standard
---	----------

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

## SPEED PACKAGE

6 Mph with Lights and Indicators.	Standard
-----------------------------------	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**Motability ID 706355 - 25/11/22 V1**



**Crash Tested**

