

FUSION E FLEX CONTOUR 6 MPH

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Adjustable Seat Width & Depth Package
- TRU-Comfort 2 Seat Cushion
- Stealth Headrest
- Fold-forward Back with Stealth Pro Back Support
- 14" Pneumatic Drive Wheel
- 10" Pneumatic Castors



- Rear Suspension System
- Black Frame
- Yellow Wheel
- Accent Inserts
- VR2 Control System
- 2 x 40 ah Batteries
- Max User Weight:
23 st 6lbs. (150 kg)
- Stealth Pro back Support

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

SEAT PACKAGE

18"W x 20"D. *Seat width can be adjusted between 16-23". Seat depth can be adjusted between 16-21".* Standard

SEAT CUSHIONS

2" Foam Seat Cushions. *Supplied with 1 x 16"W x 18"D & 1 x 20" x 22".* Standard

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.* Standard

Enter TC2 cushion size here WIDTH DEPTH *Refer to cushion order forms for size availability.*

BACK SUPPORT POSITIONING

Angle Adjustable, Fold-Forward. *3° - 33° angle adjustment in 4 steps of 10°.* Standard

STEALTH PRO BACK SUPPORTS

Active Flex, Coolcore Cover. <i>Tension adjustable, 12 cm contour.</i>					Standard
Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC	
Back Width	16"W.		16.5"W.		17"W.
	18"W.		19"W.		20"W.
	21"W.		22"W.		

SEAT POSITIONING

Static Seat. <i>18.5" front seat to floor and 17" rear seat to floor height from seat pan.</i>	Standard
------------------------------------------------------------------------------------------------	----------

ARM SUPPORTS

Height Adjustable, Removable. <i>Includes 15"L x 3"W arm pads. Adjustment range 10.5" to 14.5" from seat pan.</i>	Standard
-------------------------------------------------------------------------------------------------------------------	----------

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 13" to 16" from seat pan.</i>	Standard
----------------------------------------------------------------------------	----------

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--------------------------------------------	----------

NON-EXPANDABLE CONTROLS

VR2 Hand Control.	Standard
-------------------	----------

HAND CONTROL MOUNTS

Right Side with Fixed Position Mount. <i>Can be moved to left side.</i>	Standard
-------------------------------------------------------------------------	----------

STANDARD HEAD SUPPORTS

Fusion 14" Adjustable Comfort. <i>FRMASMB16245.</i>	Standard
-----------------------------------------------------	----------

SPEED PACKAGE

6 Mph with Lights and Indicators.	Standard
-----------------------------------	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 706353 - 18/08/23 V1



Crash Tested



V1: 14/06/2022

