

## Q1450 WITH CONTOURED SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_  
 Order Date \_\_\_\_\_  
 Address \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Phone \_\_\_\_\_ Email \_\_\_\_\_  
 PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Quantum Specialist \_\_\_\_\_  
 Ⓢ Clients Weight \_\_\_\_\_

### BASE MODEL

- 5 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 3" x 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Rear Castors



- 1.5" Anti-tip Wheels
- Swing-arm Suspension System
- VR2 Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h)
- MotorsMax User Weight:  
42 st 8lbs. (272 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Grey Tyres.	Standard
3" x 14" Pneumatic Grey Tyres.	FOC

### REAR CASTORS

2" x 9" Flat Free Grey Tyres.	Standard
2" x 9" Pneumatic Grey Tyres.	FOC

### COLOURS

Back in Black	Standard
Electric Blue	FOC
Red Corvette	FOC
Evolution Orange	FOC
Pretty in Pink	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEATING OPTIONS

Pinchless Hinge High Back Captain's Seat. <i>Includes a head support and contoured seat cushion.</i>	22"W-20"D.	FOC	22"W-22"D.	FOC	22"W-24"D.	Standard
	24"W-20"D.	FOC	24"W-22"D.	FOC	24"W-24"D.	FOC
HD High Back Captain's Seat. <i>Includes a head support, lumbar cut-out, static hinges, contoured seat cushion.</i>			26"W-20"D.	FOC	26"W-22"D	FOC
Pinchless Hinge High Back Contoured Seat. <i>Includes a head support and solid seat pan.</i>					22"W-20"D	FOC
Special Order Contoured Seat.	Seat type		Width		Depth	POA

## STEALTH SEAT CUSHIONS

Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	STD
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may only ordered with the Contoured Pinchless Hinge High Back seat and the solid seat pan option.

## SEAT POSITIONING

Static Seat.				Standard
Static Seat to Floor Height. <i>From seat pan.</i>	17.5"	18.5"	19.5"	

## ARM SUPPORTS

Bariatric Single Post Tool Height Adjustable Removable. <i>Includes Bariatric full length arm pads.</i>				Standard
Single Post HD Height Adjustable Removable.				FOC
HD & Bariatric Height Package. <i>From seat pan.</i>	8"-12"	FOC	10"-14"	Standard
Single Post HD Arm Support Pads.	10" Standard. <i>Not available with Bariatric.</i>			FOC
	14" Standard. <i>Not available with Bariatric.</i>			Standard
	10" Waterfall. <i>Not available with Bariatric.</i>			FOC
	14" Waterfall. <i>Not available with Bariatric.</i>			FOC
Flip-up Arm Supports <i>Includes arm pads.</i>	Full Length, Flip-up, Height Adjustable, Removable. <i>Required with pinchless hinge high back. Not available with HD seat</i>			FOC
	Desk Length, Flip-up, Height Adjustable. <i>Removable Required with pinchless hinge high back. Not available with HD seat.</i>			FOC

## CENTRE MOUNTED LEG SUPPORTS

17"W x 13"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.</i>	Standard
17"W x 16"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.</i>	FOC

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away Legrests.						FOC
70° HD Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
70° HD Footplates.	6"W x 8"D.	Standard	6"W x 11"D.	FOC	8"W x 11"D.	FOC

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	<b>LEFT</b>	£300	Limb Support	£300	<b>RIGHT</b>
--	-------------	------	--------------	------	--------------

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard
Auto Buckle Lap Belt Length.	38".	FOC	50".	FOC	60".	Standard
	70".	FOC	80".	FOC	90".	FOC

## NON-EXPANDABLE CONTROLS

VR2 Controls. <i>With no power options.</i>						Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>					Standard

## HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>						FOC
Mounting Position.	Left		Right			

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.						Standard
---	--	--	--	--	--	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**Motability ID 713664 - 24/06/23 V1**

