

## Q1450 WITH SYNERGY SEATING

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

### BASE MODEL

- 5 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 3" x 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Rear Castors



- 1.5" Anti-tip Wheels
- Swing-arm Suspension System
- VR2 Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
42 st 8lbs. (272 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Grey Tyres.	Standard
3" x 14" Pneumatic Grey Tyres.	FOC

### REAR CASTORS

2" x 9" Flat Free Grey Tyres.	Standard
2" x 9" Pneumatic Grey Tyres.	FOC

### COLOURS

Back in Black	Standard
Electric Blue	FOC
Red Corvette	FOC
Evolution Orange	FOC
Pretty in Pink	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

20"-24"W.
24"-26"W
26"-28"W
28"-30"W
30"-32"W
Set width at:

## SEAT PACKAGES - SEAT DEPTH

20"-24"D.
24"-28"D
Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	STD
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Static Seat.	Standard
Static Seat to Floor Height. <i>From seat pan.</i>	16.5" 17.5" 18.5"

## SYNERGY BACK SUPPORTS

Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i>	Standard
HD Sling Back. <i>Height range 14" to 25" tall from seat pan.</i>	FOC
Enter height preference here.	

## BACK SUPPORT POSITIONING

Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
1" Endo Cane Mounts. <i>Reduces backrest width by 2" compared to seat width.</i>	FOC
2" Endo Cane Mounts. <i>Reduces backrest width by 4" compared to seat width.</i>	FOC
Back Cane Height.	16"H. 17"H. 18"H. 19"H. 20"H.
<i>From seat pan.</i>	21"H. 22"H. 23"H. 24"H. 25"H.

## ARM SUPPORTS

HD Single Post Tool Height Adjustable Removable.	Standard	
Bariatric Single Post Tool Height Adjustable Removable.	FOC	
HD & Bariatric Height Package. <i>From seat pan.</i>	8"-12" FOC 10"-14" Standard	
Cantilever Flip-Up Armrests. <i>Adj. range 11"-17". Not available with HD sling back support. Increases overall width by 7".</i>	FOC	
Arm Support Pads.	10" Standard. <i>Not available with Bariatric.</i>	FOC
	14" Standard. <i>Not available with Bariatric.</i>	Standard
	10" Waterfall. <i>Not available with Bariatric.</i>	FOC
	14" Waterfall. <i>Not available with Bariatric.</i>	FOC
	Bariatric. <i>Bariatric arm supports only.</i>	FOC

## CENTRE MOUNTED LEG SUPPORTS

17"W x 13"D Base Mounted Flip-up Foot Platform. *Mounted to chassis. Will not move with power tilt.* Standard  
*Ignore other centre mount options if selected.*

17"W x 16"D Base Mounted Flip-up Foot Platform. *Mounted to chassis. Will not move with power tilt.* FOC  
*Ignore other centre mount options if selected.*

## NON-EXPANDABLE CONTROLS

VR2 Controls. *With no power options.* Standard

Control Method. Hand Control. *Select position in Hand Control Mounts below.* Standard

Attendant Control. *Hand control mounted at rear.* FOC

## HAND CONTROL MOUNTS

Fixed. Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.* FOC

Mounting Position. Left Right

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators. Standard

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**Motability ID 713665 - 24/06/23 V1**

