

Q6 EDGE 2.0 WITH TB3 SEAT / iLEVEL / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 13 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|--|----------|
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | Standard |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC |

COLOURS

| | | | |
|----------------|----------|-------------------|-----|
| Back in Black | Standard | White Satin* | FOC |
| Green Machine | FOC | Lemon Crush* | FOC |
| Midnight Blue | FOC | Night Train* | FOC |
| Electric Blue | FOC | Evolution Orange* | FOC |
| Pretty in Pink | FOC | Military Green* | FOC |
| Deep Purple | FOC | Grey Street* | FOC |
| Red Corvette | FOC | | |

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

14"-16"W.

16"-20"W

20"-22"W

Set width at:

SEAT PACKAGES - SEAT DEPTH

14"-16"D.

16"-20"D

20"-22"D

Set depth at:

Minimum seat size of 14"W x 14"D required for powered centre mount foot platform when iLevel is selected.

STEALTH SEAT CUSHIONS

No Seat Cushion Required

FOC

Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. *158.7 Kg weight limit.*

Standard

Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. *160 Kg weight limit.*

FOC

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.*

FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

iLevel Powered Seat Lift & Tilt. *12" lift & 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options.*

Standard

Powered iLevel Lift & Tilt Seat to Floor Height. *From seat pan.*

17.5"

18.5"

19.5"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. *Select TB3 back options below.*

Standard

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts.

Standard

Select size in Back Support Dimensions.

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Powered Recline. *175° with sheer reduction. Uses 1 power option.*

Standard

BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB3 Back Support System.

| Backrest Support Width Package. | | 12"-14"W. | 14"-16"W. | 16"-20"W. |
|---|-------|-----------|-----------|-----------|
| | | 18"-22"W. | 20"-24"W. | |
| Back Support Width. | 12"W. | 14"W. | 16"W. | 17"W. |
| | 18"W. | 19"W. | 20"W. | 21"W. |
| | 22"W. | 24"W. | | |
| Back Support Height. <i>From seat pan.</i> | 19"H. | 20"H. | 21"H. | 22"H. |
| | 23"H. | 24"H. | 25"H. | 26"H. |
| | 27"H. | 28"H. | 29"H. | 30"H. |
| | 31"H. | 32"H. | | |

ARM SUPPORTS

| | | | | | | |
|--|--------------|----------|--|--|--|----------|
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i> | | | | | | Standard |
| Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i> | | | | | | FOC |
| Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i> | | | | | | FOC |
| TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i> | | | | | | FOC |
| Arm Support Pads. | 10" Standard | Standard | 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | | | FOC |
| | 14" Standard | FOC | 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | | | FOC |

CENTRE MOUNTED LEG SUPPORTS

| | | | | | | |
|--|--|-----|----------|----------|----------|----------|
| Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i> | | | | | | Standard |
| Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i> | | | | | | FOC |
| Lower Leg Extensions. <i>From seat pan.</i> | 8"-12". <i>Paediatric centre mount only.</i> | FOC | 10"-16". | Standard | 13"-19". | FOC |
| Footplates. | 10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i> | | | | | FOC |
| | 10"W x 8"D Split Footplates. | | | | | FOC |
| | 10"W x 11"D Split Footplates. | | | | | Standard |
| Calf Support. | Individual Calf Pads. <i>6"W x 6"H.</i> | | | | | Standard |
| | Calf Panel. <i>18"W x 8"H. Not available with Paediatric centre mount.</i> | | | | | Standard |

INDIVIDUALLY MOUNTED LEG SUPPORTS

| | | | | | | |
|---|--|--|--|--|--|-----|
| 70° Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i> | | | | | | FOC |
| 90° Paediatric Swing-Away Leg Supports. <i>4"W x 6"D angle & depth adj. footplates. Adjustment range 7" to 11".</i> | | | | | | FOC |
| Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i> | | | | | | FOC |

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

| | | | | | | |
|--|---------|-----|-------------|----------|-------------|-----|
| 60° to 90° Tool Adjustable Swing-Away Leg Supports. <i>Flip-up extruded footplates 8.5"W x 6"D. Select lowers below. Minimum seat width 15".</i> | | | | | | FOC |
| Lower Leg Extensions. <i>From seat pan.</i> | 11"-13" | FOC | 13.5"-16.5" | Standard | 15.5"-19.5" | FOC |

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

STEALTH HEAD SUPPORTS

| | | | | | | |
|--|--|--|--|--|--|-----|
| UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i> | | | | | | FOC |
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | | | | | | FOC |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | | | | | | FOC |

Please visit stealthproducts.com for alternative head supports.

PELVIC BELTS

| | | | | | | |
|--|--|--|--|--|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | | | | | | Standard |
|--|--|--|--|--|--|----------|

EXPANDABLE CONTROLS

Q-Logic 3e Controls. *5 power option ready.*

Standard

iAccess OLED 4 Button Switch Module.
*Mounted to opposite armrest unless specified.
 Supplied FOC when ilevel is selected.*

FOC

Push Button iAccess
 Toggle Lever iAccess

If iAccess is selected choose either
 Push Button or Toggle lever option.

Control Method.

Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

Q-Logic 3e hand controls are supplied with built-in Bluetooth.

HAND CONTROL MOUNTS

Fixed.

Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.*

FOC

Mounting Position.

Left

Right

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.

Standard

6 Mph Program with Lights and indicators.

FOC

ACCESSORIES

Prepare Chassis for Dahl Docking System. *Dahl system not included.*

FOC

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 718709 - 23/04/23 V1



Crash Tested



V1: 14/06/2022

