

MOTABILITY ID 706911

Q6 EDGE 2.0 WITH TB3 SEAT / iLEVEL / POWER TILT / POWER RECLINE / POWER FOOTBOARD

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 13 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	White Satin*	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Pretty in Pink	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC		

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

14"-16"W.

16"-20"W

20"-22"W

Set width at:

SEAT PACKAGES - SEAT DEPTH

14"-16"D.

16"-20"D

20"-22"D

Set depth at:

Minimum seat size of 14"W x 14"D required for powered centre mount foot platform when iLevel is selected.

STEALTH SEAT CUSHIONS

No Seat Cushion Required

FOC

Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. *158.7 Kg weight limit.*

Standard

Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. *160 Kg weight limit.*

FOC

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.*

FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

iLevel Powered Seat Lift & Tilt. *12" lift & 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options.*

Standard

Powered iLevel Lift & Tilt Seat to Floor Height. *From seat pan.*

17.5"

18.5"

19.5"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. *Select TB3 back options below.*

Standard

Include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts.

Standard

Select size in Back Support Dimensions.

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Powered Recline. *175° with sheer reduction. Uses 1 power option.*

Standard

BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB3 Back Support System.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.
		18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.		
Back Support Height. <i>From seat pan.</i>	19"H.	20"H.	21"H.	22"H.
	23"H.	24"H.	25"H.	26"H.
	27"H.	28"H.	29"H.	30"H.
	31"H.	32"H.		

ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				Standard
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>				FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>				FOC
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				FOC
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

CENTRE MOUNTED LEG SUPPORTS

Power Elevating Articulating Centre Mount Foot Platform. <i>Minimum seat width 14". Uses 1 power option.</i>				Standard
Lower Leg Extensions.	10"-16".	FOC	13"-19".	Standard
<i>From seat pan.</i>				
Footplates.	10"W x 8"D Split Footplates.			FOC
	10"W x 11"D Split Footplates.			Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>			Standard
	Calf Panel. <i>18"W x 8"H.</i>			Standard

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>				FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>				FOC
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>				FOC

Please visit stealthproducts.com for alternative head supports.

PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.				Standard
--	--	--	--	----------

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>5 power option ready.</i>				Standard
iAccess OLED 4 Button Switch Module.	FOC	} ————	Push Button iAccess	If iAccess is selected choose either Push Button or Toggle lever option.
<i>Mounted to opposite armrest unless specified. Supplied FOC when ilevel is selected.</i>			Toggle Lever iAccess	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>			Standard
	Attendant Control. <i>Hand control mounted at rear.</i>			FOC

Q-Logic 3e hand controls are supplied with built-in Bluetooth.

HAND CONTROL MOUNTS

Fixed.				Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>				FOC
Mounting Position.	Left		Right	

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	Standard
6 Mph Program with Lights and indicators.	FOC

ACCESSORIES

Prepare Chassis for Dahl Docking System. <i>Dahl system not included.</i>	FOC
---	-----

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
MOTABILITY ID 706911 - 23/04/23 V1



Crash Tested



V1: 14/06/2022

