

Q6 EDGE HD WITH SYNERGY SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Clients Weight _____

BASE MODEL

- 10 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 35 st 7lbs. (226 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

COLOURS

Back in Black	Standard	Red Corvette	FOC
Night Train	FOC	Pretty in Pink	FOC
Electric Blue	FOC	Evolution Orange	FOC
Midnight Blue	FOC	Grey Street	FOC
Deep Purple	FOC	Green Machine	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

20"-24"W
24"-26"W
26"-28"W
Set width at:

SEAT PACKAGES - SEAT DEPTH

15"-20"D
20"-24"D
24"-28"D
Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	STD
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.	Standard
Static Seat to Floor Height. <i>From seat pan.</i>	16.75" 17.75" 18.75"

SYNERGY BACK SUPPORTS

Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i>	Standard
HD Sling Back. <i>Height range 14" to 25" tall from seat pan.</i>	FOC
Enter height preference here.	
Back Cane Height.	14"H. 15"H. 16"H. 17"H. 18"H. 19"H.
<i>From seat pan.</i>	20"H. 21"H. 22"H. 23"H. 24"H. 25"H.

STANDARD BACK SUPPORT POSITIONING

Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
---	----------

ARM SUPPORTS

HD Single Post Tool Height Adjustable Removable.	Standard
Bariatric Single Post Tool Height Adjustable Removable.	FOC
HD & Bariatric Height Package. <i>From seat pan.</i>	8"-12" FOC 10"-14" Standard
Arm Support Pads	10" Standard. <i>Not available with Bariatric.</i> FOC 10" Waterfall. <i>Not available with Bariatric.</i> FOC
	14" Standard. <i>Not available with Bariatric.</i> Standard 14" Waterfall. <i>Not available with Bariatric.</i> FOC
	Bariatric. <i>Bariatric arm supports only.</i> FOC

CENTRE MOUNT LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis.</i>	Standard
--	----------

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away Legrests.						FOC
70° HD Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
70° HD Footplates.	6"W x 8"D.	Standard	6"W x 11"D.	FOC	8"W x 11"D.	FOC
HD Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20". Min. seat width 15"</i>						FOC
HD Manual Elevating Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
HD Manual Elevating Footplates.	Flip-up Extruded. <i>7.25"W x 6"D</i>					FOC
	Flip-up Extruded. <i>8.5"W x 6"D</i>					Standard
	6"W x 8"D.	FOC	6"W x 11"D.	FOC	8"W x 11"D.	FOC

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M/CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M/CP14CCFPM.</i>	FOC
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M/ACP150M.</i>	FOC

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard
Auto Buckle Lap Belt Length.	38".	FOC	50".	FOC	60".	Standard
	70".	FOC	80".	FOC	90".	FOC

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>						Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>					Standard
	Attendant Control. <i>Hand control mounted at rear.</i>					FOC

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						FOC
Mounting Position.	Left		Right			

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.						Standard
---	--	--	--	--	--	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 721035 - 21/06/23 V1



Crash Tested

