



MOTABILITY ID 721025

Q6 EDGE HD WITH SYNERGY SEATING / POWER TILT

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Clients Weight _____

BASE MODEL

- 10 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
35 st 7lbs. (226 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|---------------------------------|----------|
| 3" x 14" Flat Free Black Tyres. | Standard |
| 3" x 14" Pneumatic Black Tyres. | FOC |

COLOURS

| | | | |
|---------------|----------|------------------|-----|
| Back in Black | Standard | Red Corvette | FOC |
| Night Train | FOC | Pretty in Pink | FOC |
| Electric Blue | FOC | Evolution Orange | FOC |
| Midnight Blue | FOC | Grey Street | FOC |
| Deep Purple | FOC | Green Machine | FOC |

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

| |
|---------------|
| 20"-24"W |
| 24"-26"W |
| 26"-28"W |
| Set width at: |

SEAT PACKAGES - SEAT DEPTH

| |
|---------------|
| 15"-20"D |
| 20"-24"D |
| 24"-28"D |
| Set depth at: |

STEALTH SEAT CUSHIONS

| | |
|---|-----|
| No Seat Cushion Required | FOC |
| Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i> | STD |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i> | FOC |
| Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i> | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

| | |
|--|-----------------------------|
| HD Powered Seat Tilt. <i>0° to 45°. Uses 1 power option.</i> | Standard |
| HD Tilt Seat to Floor Height. <i>From seat pan.</i> | 18.5" 19.5" 20.5" |

SYNERGY BACK SUPPORTS

| | |
|--|----------|
| Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i> | Standard |
| HD Sling Back. <i>Not available with HD powered seat tilt. Height range 14" to 25" tall from seat pan.</i> | FOC |

Enter height preference here.

| | | | | | | |
|-----------------------|-------|-------|-------|-------|-------|-------|
| Back Cane Height. | 14"H. | 15"H. | 16"H. | 17"H. | 18"H. | 19"H. |
| <i>From seat pan.</i> | 20"H. | 21"H. | 22"H. | 23"H. | 24"H. | 25"H. |

STANDARD BACK SUPPORT POSITIONING

| | |
|---|----------|
| Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i> | Standard |
|---|----------|

ARM SUPPORTS

| | |
|---|--|
| HD Single Post Tool Height Adjustable Removable. | Standard |
| Bariatric Single Post Tool Height Adjustable Removable. | FOC |
| HD & Bariatric Height Package. <i>From seat pan.</i> | 8"-12" FOC 10"-14" Standard |
| Arm Support Pads | 10" Standard. <i>Not available with Bariatric.</i> FOC 10" Waterfall. <i>Not available with Bariatric.</i> FOC |
| | 14" Standard. <i>Not available with Bariatric.</i> Standard 14" Waterfall. <i>Not available with Bariatric.</i> FOC |
| | Bariatric. <i>Bariatric arm supports only.</i> FOC |

CENTRE MOUNT LEG SUPPORTS

| | | |
|--|--|------------------------------------|
| Base Mounted Flip-up Foot Platform. <i>Mounted to chassis.</i> | Standard | |
| Lower Leg Extensions. <i>From seat pan.</i> | 10"-16" Standard 13"-19" FOC | |
| Footplates. | 10"W x 8"D Split Footplates. FOC 10"W x 11"D Split Footplates. Standard | |
| Calf Support. | Individual Calf Pads. <i>6"W x 6"H.</i> FOC | Calf Panel. <i>18"W x 8"H.</i> FOC |
| | Individual Large Calf Pads. <i>8"W x 8"H tapering to 6"W.</i> Standard | |

INDIVIDUALLY MOUNTED LEG SUPPORTS

| | | | | | | |
|--|---------------------------------------|----------|-------------|----------|-------------|----------|
| 70° HD Drop-in Swing-away Legrests. | | | | | | FOC |
| 70° HD Lower Extensions. <i>From seat pan.</i> | 8.5"-12". | FOC | 13"-18". | Standard | 16"-21". | FOC |
| 70° HD Footplates. | 6"W x 8"D. | Standard | 6"W x 11"D. | FOC | 8"W x 11"D. | FOC |
| HD Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20". Min. seat width 15"</i> | | | | | | FOC |
| HD Manual Elevating Lower Extensions. <i>From seat pan.</i> | 8.5"-12". | FOC | 13"-18". | Standard | 16"-21". | FOC |
| HD Manual Elevating Footplates. | Flip-up Extruded. <i>7.25"W x 6"D</i> | | | | | FOC |
| | Flip-up Extruded. <i>8.5"W x 6"D</i> | | | | | Standard |
| | 6"W x 8"D. | FOC | 6"W x 11"D. | FOC | 8"W x 11"D. | FOC |

STEALTH HEAD SUPPORTS

| | |
|---|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M/CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M/CP14CCFPM.</i> | FOC |
| UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M/ACP150M.</i> | FOC |

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

| | | | | | | |
|--|------|-----|------|-----|------|----------|
| Auto Buckle Lap Belt, Push Button Release. | | | | | | Standard |
| Auto Buckle Lap Belt Length. | 38". | FOC | 50". | FOC | 60". | Standard |
| | 70". | FOC | 80". | FOC | 90". | FOC |

EXPANDABLE CONTROLS

| | | | | | | |
|---|--|--|--|--|--|----------|
| Q-Logic 3e Controls. <i>2 power option ready.</i> | | | | | | Standard |
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i> | | | | | Standard |
| | Attendant Control. <i>Hand control mounted at rear.</i> | | | | | FOC |

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

| | | | | | | |
|---|------|--|-------|--|--|----------|
| Fixed. | | | | | | Standard |
| Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i> | | | | | | FOC |
| Mounting Position. | Left | | Right | | | |

SPEED PROGRAM & LIGHTS PACKAGE

| | | | | | | |
|---|--|--|--|--|--|----------|
| 4 Mph Program with Lights and indicators. | | | | | | Standard |
|---|--|--|--|--|--|----------|

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 721025 - 13/12/23 V1



Crash Tested



V1: 14/06/2022