



MOTABILITY ID 1001071

R-TRAK WITH TB FLEX SEATING / POWER TILT / PRO BACK

DEALERS DETAILS

Dealers Name _____
 Order Date _____
 Address _____
 _____ Postcode _____
 Phone _____ Email _____
 PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____
 _____ Postcode _____
 Quantum Specialist _____
 Ⓢ Clients Weight _____

BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 3" x 9" Flat Free Front Castors



- SRS Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
3" x 9" Pneumatic Black Tyres.	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush. <i>Extended lead-time.</i>	FOC
Midnight Blue	FOC	Night Train. <i>Extended lead-time.</i>	FOC
Electric Blue	FOC	Evolution Orange. <i>Extended lead-time.</i>	FOC
Blue Jean	FOC	Military Green. <i>Extended lead-time.</i>	FOC
Deep Purple	FOC	Grey Street. <i>Extended lead-time.</i>	FOC
Red Corvette	FOC	Raspberry Beret. <i>Extended lead-time.</i>	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

15"W.	FOC	18"W.	STD	21"W.	FOC
16"W.	FOC	19"W.	FOC	22"W.	FOC
17"W.	FOC	20"W.	FOC		

SEAT PACKAGES - SEAT DEPTH

15"D.	FOC	18"D.	Standard
16"D.	FOC	19"D.	FOC
17"D.	FOC	20"D.	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>				Standard
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18"	<i>Changes Tilt range to -5° to 25°</i>	18.5"
				<i>Changes Tilt range to -10° to 20°</i>

STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
--	----------

STANDARD BACK SUPPORT

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	FOC
---	-----

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	FOC			
Active, Coolcore Cover. <i>Solid back, 12 cm contour</i>	Standard			
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	FOC			
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	FOC			
Back Height	50cm. <i>Actual height.</i>	Standard	57cm. <i>Actual height.</i>	FOC

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard			
TB Flex Removable, Quick Release. <i>Adj. range 9" to 14" from seat pan.</i>	FOC			
Single Post Quick Height Adjustable Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14"</i>	FOC			
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16"</i>	FOC			
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14"</i>	Standard					
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>	FOC					
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". <i>Manually Adj. & powered elevating only</i>	FOC	13"-19". <i>Manually Adj. & powered elevating only</i>	Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>					FOC
	10"W x 8"D Split Footplates.	FOC	10"W x 11"D Split Footplates.		Standard	
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>		FOC	

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Adjustment range 12" to 15".</i>	Standard
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	FOC
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i>	FOC

Enter seat pan to heel measurement.

If no measurement is supplied the standard factory length option will be fitted.

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>2 power option ready.</i>	Standard	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

HAND CONTROL MOUNTS

Fixed.	Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>	FOC

Mounting Position.	Left	Right
--------------------	------	-------

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	FOC
6 Mph Program with Lights and indicators.	Standard
LED Front Castor and Rear Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Front Castor and Rear Fender Lights.	FOC

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 1001071 - 06/10/23 V1



V1: 14/06/2022