

## R-TRAK WITH TB3 SEATING / POWER TILT / POWER RECLINE

### DEALERS DETAILS

Dealers Name \_\_\_\_\_  
 Order Date \_\_\_\_\_  
 Address \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Phone \_\_\_\_\_ Email \_\_\_\_\_  
 PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Quantum Specialist \_\_\_\_\_  
 Ⓢ Clients Weight \_\_\_\_\_

### BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Front Castors



- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

12"-16"W.
14"-16"W.
16"-20"W.
20"-22"W.
Set width at:

## SEAT PACKAGES - SEAT DEPTH

12"-16"D.
14"-16"D.
16"-20"D.
20"-22"D.
Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	Standard		
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.75"	18.75"	19.75"

## BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
--	----------

Include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
--	----------

## TB3 BACK SUPPORT POSITIONING

TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	Standard
---	----------

## BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB3 Back Support System.

Backrest Support Width Package.	12"-14"W.	14"-16"W.	16"-20"W.	
	18"-22"W.	20"-24"W.		
Back Support Width.	12"W.	14"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.		
Back Support Height. <i>From seat pan.</i>	19"H.	20"H.	21"H.	22"H.
	23"H.	24"H.	25"H.	26"H.
	27"H.	28"H.	29"H.	

## ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i>		Standard
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>		FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>		FOC
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>		FOC
Height Adjustable, Removable, Heavy Duty.		FOC
HD Removable Height Package. <i>From seat pan.</i>	8"-12".	FOC 10"-14". Standard
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>		FOC
Arm Support Pads.	10" Standard	Standard
	14" Standard	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	FOC
	10" Waterfall.	FOC
	14" Waterfall.	FOC

## CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>		Standard		
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>		FOC		
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC 10"-16". <i>Manually adjustable only.</i>	FOC 13"-19". <i>Manually adjustable only.</i>	Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>		FOC	
	10"W x 8"D Split Footplates.		FOC	
	10"W x 11"D Split Footplates.		Standard	
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>	FOC

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>		FOC
90° Paediatric Swing-Away Leg Supports. <i>4"W x 6"D angle &amp; depth Adj. footplates. Adjustment range 7" to 11".</i>		FOC
Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>		FOC

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° TO 90°	60° to 90° Tool Adjustable Swing-Away Leg Supports. <i>Flip-up extruded footplates 8.5"W x 6"D. Select lowers below. Minimum seat width 15"</i>		FOC
	Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC 13.5"-16.5" Standard 15.5"-19.5" FOC

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>		FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>		Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>		FOC

## PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.

Standard

## EXPANDABLE CONTROLS

Q-Logic 3e hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls. *2 power option ready.*

Standard

Control Method.

Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

## HAND CONTROL MOUNTS

Fixed.

Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.*

FOC

Mounting Position.

Left

Right

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.

Standard

6 Mph Program with Lights and indicators.

FOC

LED Front Castor and Rear Fender Lights. *Includes single switch button mounted to side rail. Not road legal.*

Standard

Omit LED Front Castor and Rear Fender Lights.

FOC

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**MOTABILITYID 1000814 - 31/06/23 V1**



**Crash Tested**

