

FUSION E EXTRA DEEP CONTOUR 6 MPH

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

Height Adjustable Removable Armrests
 70° Swing-Away Legrests
 Adjustable Seat Width & Depth Package
 TRU-Comfort 2 Seat Cushion
 Stealth Headrest
 Fold-forward Back with Stealth Pro Back Support
 14" Pneumatic Drive Wheel
 10" Pneumatic Castors



Rear Suspension System
 Black Frame
 Yellow Wheel
 Accent Inserts
 VR2 Control System
 2 x 50 ah Batteries
 Max User Weight:
 23 st 6lbs. (150 kg)
 Stealth Pro Back Support

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

SEAT PACKAGE

18"W x 20"D. *Seat width can be adjusted between 16-23". Seat depth can be adjusted between 16-21".* Standard

SEAT CUSHIONS

2" Foam Seat Cushions. *Supplied with 1 x 16"W x 18"D & 1 x 20" x 22".* Standard

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.* Standard

Enter TC2 cushion size here

WIDTH

DEPTH

Refer to cushion order forms for size availability.

BACK SUPPORT POSITIONING

Angle Adjustable, Fold-Forward. *3° - 33° angle adjustment in 4 steps of 10°.* Standard

STEALTH PRO BACK SUPPORTS

Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				Standard
Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC
Back Width	16"W.		16.5"W.	17"W.
	18"W.		19"W.	20"W.
	21"W.		22"W.	

SEAT POSITIONING

Static Seat. <i>18.5" front seat to floor and 17" rear seat to floor height from seat pan.</i>	Standard
--	----------

ARM SUPPORTS

Height Adjustable, Removable. <i>Includes 15"L x 3"W arm pads. Adjustment range 10.5" to 14.5" from seat pan.</i>	Standard
---	----------

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 13" to 16" from seat pan.</i>	Standard
--	----------

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

NON-EXPANDABLE CONTROLS

VR2 Hand Control.	Standard
-------------------	----------

HAND CONTROL MOUNTS

Right Side with Fixed Position Mount. <i>Can be moved to left side.</i>	Standard
---	----------

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

SPEED PACKAGE

6 Mph with Lights and Indicators.	Standard
-----------------------------------	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 1000397 - 05/12/24 V1

