

OUTBACK WITH CONTOURED SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

£11,500

- 5 Colour Options
- Flip-up Arm Supports
- Centre Mount Foot Platform
- Contoured Seating
- 15" Pneumatic Drive Tyres
- 4 Motor Independent Drive

- SRS Suspension System
- Q-Logic 3e Control System
- 2 x 95 ah Batteries
- 6 Mph (9.65 Km/h) 4-Pole Motors
- Max User Weight: 25 st (159 kg)



Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

COLOURS

Cadet Matte	Standard
Tan Matte	FOC
Desert Camo	£150
Black Camo	£150
White Camo	£150



Please download our free PQ-365 app to view colours and other product features in more detail.

CONTOURED SEAT PACKAGES

16"W x 16"D	FOC	20"W x 16"D*	FOC
16"W x 18"D*	FOC	20"W x 18"D	FOC
16"W x 20"D*	FOC	20"W x 20"D	FOC
18"W x 16"D	FOC	20"W x 22"D*	FOC
18"W x 18"D	Standard	22"W x 18"D*	FOC
18"W x 20"D	FOC	22"W x 20"D	FOC
18"W x 22"D*	FOC	22"W x 22"D	FOC

*Extended lead time. All seat packages above include a solid seat pan, contoured seat cushion, and a standard head support.

STEALTH SEAT CUSHIONS

No Seat Cushion Required.	FOC
Contoured Seat Cushion. <i>136 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£165
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£182
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£209
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£231
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£231
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£292
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£292

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.	Standard
Powered Seat Lift. <i>10" lift and up to 3.2 Mph drive speed while elevated. Uses 1 power option.</i>	£1950
Seat to Floor Height. <i>From seat pan. Only the lowest seat to floor height can be selected with powered seat lift.</i>	19.5" 20.5" 21.5"

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes 14"L contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>	£100
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	£150
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>	£150
HD Height Adjustable, Removable. <i>Includes HD arm pads. Do not select arm pads below.</i>	8"-12"H. 10"-14"H. £220
Bariatric Removable Height Adjustable. <i>Includes Bariatric arm pads. Do not select arm pads below.</i>	£250
Arm Support Pads.	10" Standard STD 10" Waterfall. <i>Only with single post or paediatric arms.</i> FOC
	14" Standard FOC 14" Waterfall. <i>Only with single post or paediatric arms.</i> FOC
	10" Falcon Gel £40 Falcon Gel Mounting Plates. <i>Required for Single Post & Paediatric arm supports.</i> £40
	14" Falcon Gel £40

CENTRE MOUNTED LEG SUPPORTS

Basic Centre Mount Foot Platform	Standard
Lower Leg Extensions. <i>8"-10". Only available from seat pan. with 16" seat width.</i>	FOC 10"-14". Standard 14"-19". FOC
Footplates.	10"W x 9"D Split Footplates. Standard 10"W x 12"D Split Footplates. FOC
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i> Standard
	Individual Large Calf Pads. <i>6"W x 8"H. Only with long lower ext.</i> FOC
Calf Panel. <i>18"W x 8"H.</i>	£43
Footplate Accessories.	Footplate Connector. £10

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

	70° HD Drop-in Swing-away Legrests. <i>Includes heel loops.</i>						£360
70°	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12"	FOC	13"-18"	Standard	16"-21"	FOC
	Footplate Size.	6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
	Accessories.	Wheel Bumpers					£10
LIMB SUPPORT	Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>		LEFT	£300	Limb Support	£300	RIGHT

The 8"W footplates may cause interference with seat widths less than 19".
Each limb support includes a 60-90° Swing-away leg support hanger for mounting.

PELVIC BELTS

	Auto Buckle Lap Belt, Push Button Release.		Standard	
Standard Pelvic Belts.	38"	FOC	70"	FOC
	50"	Standard	80"	FOC
	60"	FOC	90"	FOC
Stealth Positioning 2-Point Padded Pelvic Belt, Single Pull, Push Button Release	Medium			£84
	Large			£84
	Extra Large			£84
Retractable Lap Belt	Push Button Release.			£169

STEALTH THIGH & KNEE SUPPORTS

UniLink Thigh Supports.	LEFT	£185	UniLink Hip Guide Long Removable.	£185	RIGHT
		£185	Unilink Hip Guide Short Removable.	£185	
		£51	4x4 Pelvic Thigh Pad.	£51	
		£51	4x8 Pelvic Thigh Pad.	£51	
		£51	4x12 Pelvic Thigh Pad.	£51	
Knee Adductors.	LEFT	£163	Fixed Cane Mounted.	£163	RIGHT
		£179	Swing-Away Track Mounted.	£179	
		£38	3" Round Knee Pad.	£38	
		£38	4" Round Knee Pad.	£38	

STEALTH ELBOW SUPPORTS

UniLink Elbow Supports. <i>Only available with Single Post, Paediatric, and HD arm supports.</i>	LEFT	£133	Swing-Away, Elbow Support Bracket. <i>ARMSES-L ----- ARMSES-R</i>	£133	RIGHT
		£38	3"W x 3"H Elbow Pad. <i>ARMSES-3X3</i>	£38	
		£38	4"Wx4"H Elbow Pad. <i>ARMSES-4X4</i>	£38	

STEALTH FOOT SUPPORTS

Stealth Foot Boxes (Pair).	Small. <i>5"Wx8"Dx7"H.</i>	£361
	Medium. <i>6"Wx10"Dx9"H.</i>	£361
	Large. <i>7"Wx12"Dx11"H.</i>	£361
Stealth Shoe Holders (Pair).	Small Shoe Holder, 3/4" Webbing. <i>IFSH4211</i>	£81
	Medium Shoe Holder, 3/4" Webbing. <i>IFSH4212</i>	£81
	Large Shoe Holder, 3/4" Webbing. <i>IFSH4213</i>	£81

EXPANDABLE CONTROLS Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls. <i>No power options.</i>	Standard
Q-Logic 3e Controls. <i>3 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£350
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£550
Control Method. Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard

HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away.						£75
Advanced Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						£100
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£25	Extended I-Handle	£25
	Mushroom	£25	Medium Ball	£25	3" Goal Post	£50
	Chin Cup	£25	Soft Ball	£25	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.

Standard

ACCESSORIES Supplied uninstalled.

Cup Holder.	£20
Mobile Phone Holder.	£40
Clothes Guards (Pair).	£100
Crutch/Cane Holder.	£170
Glove Box.	£30
Oxygen Holder.	£50
Transfer Handles.	£125
Stealth Tablet Holder.	£342
Large Black Weather Cover.	£125
Rear View Mirror.	£15
Quantum Back Pack.	£110
Rear Basket.	£40
Saddle Bag.	£25
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if requested.</i>	£379

NOTES & SPECIAL REQUESTS

OUTBACK WITH CONTOURED SEATING



BASE PRICE

£11,500

COST OPTIONS

TOTAL

5



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 01/04/25 V1

