

OUTBACK WITH TB3 SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

£15,950

- 5 Colour Options
- Height Adjustable Arm Supports
- Centre Mount Foot Platform
- TRU-Balance 3 Seating
- 15" Pneumatic Drive Tyres
- 4 Motor Independent Drive



- SRS Suspension System
- Q-Logic 3e Control System
- 2 x 95 ah Batteries
- 6 Mph (9.65 Km/h) 4-Pole Motors
- Max User Weight: 21st 6 lbs (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

COLOURS

Cadet Matte	Standard
Tan Matte	FOC
Desert Camo	£150
Black Camo	£150
White Camo	£150



Please download our free PQ-365 app to view colours and other product features in more detail.

SEAT PACKAGES - SEAT WIDTH

16"-20"W.
20"-22"W.
Set width at:

SEAT PACKAGES - SEAT DEPTH

16"-18"D.
19"-22"D.
Set depth at:

STEALTH SEAT CUSHIONS

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	Standard
No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£126
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£126
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£148
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£148
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations.
Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. 0° to 50° COG tilt. Uses 1 power option.			Standard
10" Powered HD Seat Lift & Tilt. 10" lift & 50° tilt. Up to 3.2 Mph drive speed while elevated. Uses 2 power options.			£ 1950
Powered Tilt Seat to Floor Height. From seat pan.	19"	20"	21"
Powered HD Lift & Tilt Seat to Floor Height.. From seat pan.	19.5"	20.5"	21.5"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a pair of plastic lateral wedge inserts

TB3 BACK SUPPORT

Active Contour Back Support with Padded Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
Stealth Adapter Interface Plate. <i>SSAPT4. For mounting third party back supports to TB3 back system. Not required for Q-Back 2.</i>	£358
ADI Q-Back 2. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts directly to TB3 system.</i>	POA

Enter ADI Q-Back 2 part number here for quotation.

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>85° to 105° angle adjustment via gas strut within back cover.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not suitable for occupied transport. Not available with TB3 flip-back arm supports or powered recline.</i>	£150
TRU-Balance 3 Powered Recline. <i>85° to 175° with sheer reduction. Uses 1 power option.</i>	£700

BACK CANE ONLY SUPPORT

Active Contour Back Support with Padded Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to canes with cane clamps.</i>	POA

Enter ADI Back Support part number here for quotation

Back Cane Only Support Positioning.	Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>			Standard
	Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>			£200
	2" Endo Cane Mounts. <i>1" each side. Reduces back width by 2". Minimum seat width 12".</i>			£40
	4" Endo Cane Mounts. <i>2" each side. Reduces back width by 4". Minimum seat width 14".</i>			£40
Back Cane Height.	16"H.	17"H.	20"H.	
<i>From seat pan.</i>	21"H.	24"H.	25"H.	

BACK SUPPORT DIMENSIONS

Back width must be within +/- 4" of seat width when TB3 Back Support has been selected, and must match the seat width when the Back Canes Only option has been selected. The only exception is when the 2" or 4" Endo mounts have been selected, in which case the back can be 2" or 4" narrower than the seat width.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.	18"-22"W.
		20"-24"W.			
Back Support Width.	12"W.	14"W.	15"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.	22"W.
Back Support Height. <i>From seat pan.</i>	14"H.	15"H.	16"H.	17"H.	18"H.
	19"H.	20"H.	21"H.	22"H.	23"H.
	24"H.	25"H.	26"H.	27"H.	28"H.
	29"H.	30"H.	31"H.	32"H.	

19" minimum backrest height when powered recline is selected.

BACK WIDTH STATIC SEAT OR POWER TILT

	12"W	14"W	15"W	16"W	17"W	18"W	19"W	20"W	21"W	22"W
14"H	T									
15"H	T									
16"H	S	S	T	S	T	S	T	S	T	S
17"H	S	S	T	S	T	S	T	S	T	S
18"H		S	T	S	T	S	T	S	T	S
19"H		S	Q	T	S	Q	T	S	Q	T
20"H		S	Q	T	S	Q	T	S	Q	T
21"H		S	Q	T	S	Q	T	S	Q	T
22"H		S	Q	T	S	Q	T	S	Q	T
23"H		S	T	S	T	S	Q	T	S	Q
24"H		S	T	S	T	S	Q	T	S	Q
25"H		S	T	S	T	S	Q	T	S	Q
26"H		S	T	S	T	S	Q	T	S	Q
27"H		S	T	S	T	S	Q	T	S	Q
28"H		S	T	S	T	S	Q	T	S	Q
29"H		S	T	S	T	S	Q	T	S	Q

BACK WIDTH POWER RECLINE

	12"W	14"W	15"W	16"W	17"W	18"W	19"W	20"W	21"W	22"W
19"H	S	S	T	S	T	S	T	S	T	S
20"H	S	S	T	S	T	S	T	S	T	S
21"H		S	T	S	T	S	T	S	T	S
22"H		S	Q	T	S	Q	T	S	Q	T
23"H		S	Q	T	S	Q	T	S	Q	T
24"H		S	Q	T	S	Q	T	S	Q	T
25"H		S	Q	T	S	Q	T	S	Q	T
26"H		S	T	S	T	S	Q	T	S	Q
27"H		S	T	S	T	S	Q	T	S	Q
28"H		S	T	S	T	S	Q	T	S	Q
29"H		S	T	S	T	S	Q	T	S	Q
30"H		S	T	S	T	S	Q	T	S	Q
31"H		S	T	S	T	S	Q	T	S	Q
32"H		S	T	S	T	S	Q	T	S	Q

Back Type.

Tru-Comfort 2

Sport Back

ADI/Q-Back 2

ARM SUPPORTS

TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>								Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i>								£279
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>								£100
HD Height Adjustable, Removable. <i>Includes HD arm pads. Do not select arm pads below.</i>						8"-12" H.	10"-14" H.	£220
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>								£350
Down Posts for TB3 Flip-back. <i>Only with TB3 Flip-back arm supports. Minimum seat/back width 16".</i>								£100
Arm Support Pads.	10" Standard	<i>Not available with HD arms.</i>	STD	10" Falcon Gel	£40	10" Waterfall.	<i>Not available with TB Flex or HD arms.</i>	FOC
	14" Standard	<i>Not available with HD arms.</i>	FOC	14" Falcon Gel	£40	14" Waterfall.	<i>Not available with TB Flex or HD arms.</i>	FOC
	17" Standard.	<i>TB3 only. Minimum seat depth 16".</i>	£50	Falcon Gel Mounting Plates. <i>Required for 2-post, Single Post & Paediatric arm supports.</i>				£40

BASIC CENTRE MOUNTED LEG SUPPORTS

Basic Centre Mount Foot Platform.						Standard
Lower Leg Extensions. <i>From seat pan.</i>	8"-10".	FOC	10"-14".	Standard	14"-19".	FOC
Split Footplates.	10"W x 9"D.		Standard	10"W x 12"D.	<i>Tapered footplates. Not available with 8"-10" lower extensions.</i>	FOC
Calf Support.	Individual Calf Pads. 6"W x 6"H.	Standard	}	Calf Panel. 18"W x 8"H.		£43
	Individual Large Calf Pads. 6"W x 8"H. <i>Only available with 14"-19" lower extensions.</i>	FOC				
Footplate Accessories.	Footplate Connector.					£10

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

POWERED CENTRE MOUNTED LEG SUPPORTS

Power Elevating Articulating Centre Mount Foot Platform. <i>Uses 1 power option.</i>							£950
Lower Leg Extensions. <i>From seat pan.</i>	11"-17".		FOC	14"-20".		Standard	
Split Footplates.	10"W x 6"D.	FOC	10"W x 8"D.	FOC	10"W x 11"D.	Standard	
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	}	Small Calf Panel. <i>18"W x 8"H.</i>		£43	
	Individual Large Calf Pads. <i>6"W x 8"H.</i>	FOC		Large Calf Panel. <i>24"W x 8"H.</i> <i>Only available with seat widths 20"-22".</i>		£85	
	Reverse Dartex Calf Pad Covers.	£42		Reverse Dartex Calf Panel Cover.		£62	
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>					£16	
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>					£16	
	Footplate Connector.	£10	Outer Edge Supports.		£10		
	Wheel Bumpers.	£10	Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.				

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

Power Elevating Articulating Swing-Away Leg Supports. <i>Includes calf pads and heel loops. Uses 2 power options.</i>							£850
POWER ELEVATING ARTICULATING	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-14.5"		FOC	12.5"-21"		Standard
	Footplate Type.	Angle & Depth Adjustable. <i>Select size below.</i>					Standard
	Footplate Size.	4"W x 6"D	FOC	5"W x 6"D	FOC	5"W x 8"D	FOC
		6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
	Accessories.	Wheel Bumpers					Standard
		Gel Padded Knee Adductors					FOC
Leg Support Operation Method.	Independent			Combined			
70° HD Drop-in Swing-away Leg Support. <i>Includes heel loops.</i>							£360
70°	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12"	FOC	13"-18"	Standard	16"-21"	FOC
	Footplate Size.	6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
	Accessories.	Wheel Bumpers					£10
LIMB SUPPORT	Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>		LEFT	£300	Limb Support	£300	RIGHT

The 8"W footplates may cause interference with seat widths less than 19". Certain configurations may result in reduced ground clearance below footplates. Each limb support includes a 60-90° Swing-away leg support hanger for mounting. Contact Quantum Technical for information.

STEALTH TRUNK SOLUTIONS

UniLink Thoracic Lateral Supports. <i>Supplied as a pair.</i>	Lateral Thoracic Swing-away Hardware. £444		Additional 1.5" Link (each). £54	
	3" H x 5" L Curved Lateral Pad. £84		Additional 2" Link (each). £54	
	4" H x 6" L Curved Lateral Pad. £84		Additional 3" Link (each). £54	
4-Point Harness.	Extra Small £143	Classic STD*	Dynamic STD*	Vinyl Face FOC
	Small £143			Zipper Face** FOC
	Medium £143	Contoured FOC	Structured FOC	
	Large £143			
	Extra Large £143			
Chest Strap.	Extra Small £91	Small £91	Side Release Standard	
	Medium £91	Large £91	Velcro Release £20	
	Extra Large £91			

*STD = Standard option **Classic Only

STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	LEFT	£185	UniLink Hip Guide Long Removable.	£185	RIGHT	
		£185	Unilink Hip Guide Short Removable.	£185		
		£51	4x4 Pelvic Thigh Pad.	£51		
		£51	4x8 Pelvic Thigh Pad.	£51		
		£51	4x12 Pelvic Thigh Pad.	£51		
Knee Adductors.	LEFT	£163	Fixed Cane Mounted.	£163	RIGHT	
		£170	Swing-Away Track Mounted.	£179		
		£38	3" Round Knee Pad.	£38		
		£38	4" Round Knee Pad.	£38		
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.				Standard	
	38".	FOC	70".	FOC		
	50".	Standard	80".	FOC		
	60".	FOC	90".	FOC		
Retractable Lap Belt.	Push Button Release.				£169	
Positioning 2-Point Padded Pelvic Belt.	Extra Small	£84	Single Pull	Standard	Side Release	Standard
	Small	£84	Rear Pull	FOC	Push Button	FOC
	Medium	£84	Dual Pull	FOC		
	Large	£84				
	Extra Large	£84				
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£143	Small	£143	Push Button	Standard
	Medium	£143	Large	£143	Side Release	£21
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>				£189	
	Small Wedge. <i>MTMD-12M.</i>		£108	Medium Wedge. <i>MTMD-13M.</i>		£108
	Large Wedge. <i>MTMD-14M.</i>		£108	Adjustable Wedge. <i>MTMD-42M.</i>		£209

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£230
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M/ ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92

Please visit stealthproducts.com for alternative head supports.

STEALTH ELBOW SOLUTIONS

UniLink Elbow Supports.	LEFT	£133	Swing-Away, Elbow Support Bracket. <i>ARMS-ES-L ----- ARMS-ES-R</i>	£133	RIGHT
		£38	3"W x 3"H Elbow Pad. <i>ARMSES-3X3</i>	£38	
		£38	4"Wx4"H Elbow Pad. <i>ARMSES-4X4</i>	£38	

STEALTH FOOT SOLUTIONS

Stealth Foot Boxes (Pair).	Small. <i>5"Wx8"Dx7"H.</i>	£361
	Medium. <i>6"Wx10"Dx9"H.</i>	£361
	Large. <i>7"Wx12"Dx11"H.</i>	£361
Stealth Shoe Holders (Pair).	Small Shoe Holder, 3/4" Webbing. <i>IFSH4211</i>	£81
	Medium Shoe Holder, 3/4" Webbing. <i>IFSH4212</i>	£81
	Large Shoe Holder, 3/4" Webbing. <i>IFSH4213</i>	£81

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>3 power option ready.</i>	Standard
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£375
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount. Not available with 3e controls.</i>	£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.
*Mounted to opposite armrest unless specified.
Supplied FOC when iLevel is selected*

£300

Push Button iAccess
Toggle Lever iAccess

If iAccess is selected choose either Push Button or Toggle lever option.

Control Method. Hand Control. *Select position in Hand Control Mounts below.*

Standard

HAND CONTROL MOUNTS

Fixed.					Standard
Swing-away.					£75
Advanced Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>					£100
Mounting Position.	Left	FOC	Right	Standard	
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£25	Extended I-Handle £25
	Mushroom	£25	Medium Ball	£25	3" Goal Post £50
	Chin Cup	£25	Soft Ball	£25	4" Goal Post £50

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.

Standard

ACCESSORIES Supplied uninstalled.

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail</i>	£110
Cup Holder.	£20
Mobile Phone Holder.	£40
Clothes Guards (Pair).	£100
Crutch/Cane Holder.	£170
Glove Box.	£30
Oxygen Holder.	£50
Transfer Handles.	£125
Stealth Tablet Holder.	£342
Large Black Weather Cover.	£125
Back-up Camera.	£300
Rear View Mirror.	£15
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if requested.</i>	£379

NOTES & SPECIAL REQUESTS

PRICE

OUTBACK WITH TB3 SEATING



BASE PRICE

£15,950

COST OPTIONS

TOTAL

7



Quantum®

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

RETAIL UK - 01/04/25 V1

