

Technical Bulletin

Bulletin Ref.:

UK-Q1-6-25

Date Issued:

20 February 2025

BULLETIN DETAILS

Product Group:

Powerchairs

Effective Date:

20 February 2025

Model(s) Affected:

Edge 3 Stretto with Contoured Seating

Description of Changes/ Diviations:

This technical bulletin outlines the updates to the Retail order form for the Quantum Edge 3 Stretto with contoured seating.

The following changes have been made:

Edge 3 Stretto with Contoured Seating – Retail

- **Page 3 – Non-Expandable Controls**
 - NE+ replaced with QLNES control system.
- **Page 3 – Non-Expandable Controls**
 - Note added to VR2 stating this option is not available with fender lights
- **Page 3 – Hand Control Mounts**
 - Swing-away model mount change
- **Page 3 – Speed Program & Lights Package**
 - Note added to LED Fender Lights stating this option is not available with VR2 controls.
- **Page 3 – ISO Branding**
 - Updated ISO logo to reflect current branding guidelines.

Documentation Affected: ☐ Owner's Manual ☒ Order Form(s) ☐ Spec Sheets ☒ Other ☐ N/A

Link to new version of affected document: <https://quantumrehab.co.uk/powerchairs/edge-3-stretto/>

For more detailed information, please contact your local designated Area Sales Manager.

NB: This Bulletin is for informational purposes only and DOES NOT replace the Specifications.

Pride Mobility Products Ltd.
32 Wedgwood Road, Bicester,
Oxfordshire, OX26 4UL

E: sales@quantumrehab.co.uk
T: +44 (0)1869 324600