

EDGE 3 STRETTO - CONTOURED SEAT / iLEVEL

DEALERS DETAILS

| | |
|--------------|-----------|
| Dealers Name | |
| Order Date | |
| Address | |
| | Postcode |
| Phone | Email |
| PO No | Order Ref |

DELIVERY DETAILS

| | |
|---|----------|
| Dealers Address | |
| | |
| | Postcode |
| Quantum Specialist | |
| Ⓢ Clients Weight (required for suspension settings) | |

BASE MODEL

14 Colour Options
 Flip-up Armrests
 70° Swing-Away Legrests
 Contoured Seat with Fold-Forward Backrest
 12.5" Flat Free Drive Tyres
 2" x 5" Flat Free Front Castors



2" x 6" Flat Free Rear Castors
 SRS Suspension System
 QLNE Control System
 2 x 62 ah Batteries
 4 Mph (6 Km/h) Motors
 Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|--|----------|
| 3" x 12.5" Flat Free Black Tyres. | Standard |
| 3" x 12.5" Pneumatic Black Tyres. | FOC |
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | FOC |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC |

COLOURS

| | | | |
|---------------|----------|-------------------|-----|
| Back in Black | Standard | Pretty in Pink | FOC |
| Green Machine | FOC | Lemon Crush* | FOC |
| Midnight Blue | FOC | Night Train* | FOC |
| Electric Blue | FOC | Evolution Orange* | FOC |
| Blue Jean | FOC | Military Green* | FOC |
| Deep Purple | FOC | Grey Street* | FOC |
| Red Corvette | FOC | Raspberry Beret* | FOC |

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



CONTOURED SEAT PACKAGES

| 16"W x 16"D | FOC | 18"W x 18"D | Standard | 20"W x 20"D | FOC |
|---|-----|-------------|----------|-------------|-----|
| Seat package includes standard headrest and auto buckle lap belt. | | | | | |

SEAT CUSHIONS

| | |
|--|----------|
| Contoured Seat Cushion. <i>136 Kg weight limit.</i> | Standard |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i> | FOC |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i> | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC |
| Seat cushion dimensions will match seat package selected above. | |

SEAT POSITIONING

| | |
|---|-------------|
| iLevel Powered Seat Lift. <i>12" lift. Uses 1 power option.</i> | Standard |
| iLevel Powered Lift Seat to Floor Height. <i>From seat pan.</i> | 18" 19" 20" |

ARM SUPPORTS

| | |
|--|--|
| Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i> | Standard |
| TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i> | FOC |
| Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i> | FOC |
| Arm Support Pads | 10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC |
| | 14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC |

LEG SUPPORTS

| | |
|--|------------------------------------|
| 70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i> | Standard |
| Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20". Min. seat width 15".</i> | FOC |
| Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range 10" to 16". 5"W x 11"L Footplates.</i> | FOC |
| Seat pan to heel measurement | (Leave blank for standard factory) |

NON-EXPANDABLE CONTROLS

| | |
|--|----------|
| QLNES Controls. <i>2 power option ready.</i> | Standard |
|--|----------|

HAND CONTROL MOUNTS

| | |
|--------------------|-------------------------|
| Fixed. | Standard |
| Basic swing-away. | FOC |
| Mounting Position. | Left FOC Right Standard |

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. *Lights and indicators not included.*

Standard

NOTES & SPECIAL REQUESTS

3



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
Motability ID 1000398 - 28/02/25 V1

