

R-TRAK WITH TB FLEX SEATING / POWER TILT

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

14 Colour Options
Height Adjustable Removable Armrests
Centre Mount Foot Platform
Angle Adjustable Fold-Forward Backrest
14" Flat Free Drive Tyres
3" x 9" Flat Free Front Castors



SRS Suspension System
QLNES Control System
2 x 70 ah Batteries
4 Mph (6 Km/h) Motors
Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
3" x 9" Pneumatic Black Tyres.	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

15"W.	FOC	19"W.	FOC
16"W.	FOC	20"W.	FOC
17"W.	FOC	21"W.	FOC
18"W.	Standard	22"W.	FOC

SEAT PACKAGES - SEAT DEPTH

15"D.	FOC	19"D.	FOC
16"D.	FOC	20"D.	FOC
17"D.	FOC		
18"D.	Standard		

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations.
Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Powered Seat Tilt. 0° to 30°. Uses 1 power option.				Standard	
Powered Tilt Seat to Floor Height. From seat pan.	17.5"	18"	Changes Tilt range to -5° to 25°	18.5"	Changes Tilt range to -10° to 20°

STANDARD BACK SUPPORT POSTIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
--	----------

STANDARD BACK SUPPORT

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	FOC
---	-----

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adj. range 9" to 14" from seat pan.</i>				Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				FOC
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not available with TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not available with TB Flex.</i>	FOC

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>				Standard
Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	FOC	13"-19".	Standard
Footplates.	10"W x 8"D Split Footplates.			FOC
	10"W x 11"D Split Footplates.			Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>			Standard

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Adjustment range 12" to 15".</i>	Standard
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20". Minimum seat width 16".</i>	FOC
Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i>	
Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.	

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

NON-EXPANDABLE CONTROLS

QLNES Controls. <i>2 power option ready.</i>	Standard
Control Method. Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
Attendant Control. <i>Hand control mounted at rear.</i>	FOC

HAND CONTROL MOUNTS

Fixed.	Standard
Basic swing-away.	FOC
Mounting Position. Left Right	

SPEED PROGRAM & LIGHTS PACKAGE

2-Pole Motors & 4 Mph Program with Lights and indicators.	FOC
2-Pole Motors & 6 Mph Program with Lights and indicators.	Standard

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 1001054 - 04/03/25 V1

