

## R-TRAK WITH TB FLEX SEATING / POWER TILT / PRO BACK

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

### BASE MODEL

14 Colour Options  
Height Adjustable Removable Armrests  
Centre Mount Foot Platform  
Angle Adjustable Fold-Forward Backrest  
14" Flat Free Drive Tyres  
3" x 9" Flat Free Front Castors



SRS Suspension System  
QLNES Control System  
2 x 70 ah Batteries  
4 Mph (6 Km/h) Motors  
Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
3" x 9" Pneumatic Black Tyres.	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush. <i>Extended lead-time.</i>	FOC
Midnight Blue	FOC	Night Train. <i>Extended lead-time.</i>	FOC
Electric Blue	FOC	Evolution Orange. <i>Extended lead-time.</i>	FOC
Blue Jean	FOC	Military Green. <i>Extended lead-time.</i>	FOC
Deep Purple	FOC	Grey Street. <i>Extended lead-time.</i>	FOC
Red Corvette	FOC	Raspberry Beret. <i>Extended lead-time.</i>	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

15"W.	FOC	18"W.	STD	21"W.	FOC
16"W.	FOC	19"W.	FOC	22"W.	FOC
17"W.	FOC	20"W.	FOC		

## SEAT PACKAGES - SEAT DEPTH

15"D.	FOC	18"D.	Standard
16"D.	FOC	19"D.	FOC
17"D.	FOC	20"D.	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations.  
Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Powered Seat Tilt. 0° to 30°. Uses 1 power option.				Standard	
Powered Tilt Seat to Floor Height. From seat pan.	17.5"	18"	Changes Tilt range to -5° to 25°	18.5"	Changes Tilt range to -10° to 20°

## STANDARD BACK SUPPORT POSTIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
--	----------

## STANDARD BACK SUPPORT

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	FOC
---	-----

## STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>				FOC
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>				Standard
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				FOC
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>				FOC
Back Height	50cm. <i>Actual height.</i>	Standard	57cm. <i>Actual height.</i>	FOC

## ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adj. range 9" to 14" from seat pan.</i>				Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				FOC
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not available with TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not available with TB Flex.</i>	FOC

## CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>				Standard
Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	FOC	13"-19".	Standard
Footplates.	10"W x 8"D Split Footplates.			FOC
	10"W x 11"D Split Footplates.			Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>			Standard

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. *Adjustment range 12" to 15".*

Standard

Manual Elevating Swing-Away Legrests. *Adjustment range 15" to 20". Minimum seat width 16".*

FOC

Enter seat pan to heel measurement.

*If no measurement is supplied the standard factory length option will be fitted.*

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops.  
Elevating legrests supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. *HMO475-17M / CP08CCFPM.*

FOC

UniLink 10" Comfort Plus. *HMO475-17M / CP10CCFPM.*

Standard

UniLink 14" Comfort Plus. *HMO475-17M / CP14CCFPM.*

FOC

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.

Standard

## NON-EXPANDABLE CONTROLS

QLNES Controls. *2 power option ready.*

Standard

Control  
Method.

Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

## HAND CONTROL MOUNTS

Fixed.

Standard

Basic swing-away.

FOC

Mounting Position.

Left

Right

## SPEED PROGRAM & LIGHTS PACKAGE

2-Pole Motors & 4 Mph Program with Lights and indicators.

FOC

2-Pole Motors & 6 Mph Program with Lights and indicators.

Standard

## NOTES & SPECIAL REQUESTS



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**Motability ID 1001071 - 04/03/25 V1**

