

4FRONT 2 HD WITH TB3 SEAT / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight (required for suspension settings)	

BASE MODEL

14 Colour Options
Height Adj. Removable Arm Supports
Centre Mount Foot Platform
TB3 Seating
14" Flat Free Drive Tyres
2" x 9" Flat Free Rear Castors



SRS Suspension System
Smart Traction Control
Q-Logic 3e Control System
2 x 70 ah Batteries
4 Mph (6 Km/h) 4-Pole Motors
Max User Weight: 32 st 14lbs. (204 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

20"-24"W.

Set width at:

SEAT PACKAGES - SEAT DEPTH

20"-24"D.

Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>294 Kg weight limit.</i>	Standard
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>294 Kg weight limit.</i>	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 3 HD Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	Standard
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.75" 18.75" 19.75"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Includes a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.	

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
---	----------

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 HD Powered Recline. <i>85° to 175° with sheer reduction. Uses 1 power option.</i>	Standard
---	----------

BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width

Backrest Support Width Package.		16"-20"W.	20"-24"W.	24"-28"W.
		28"-32"W		
Back Support Width.	16"W.	17"W.	18"W.	19"W.
	20"W.	21"W.	22"W.	24"W.
	26"W.	28"W.		
Back Support Height.	19"H.	20"H.	21"H.	22"H.
	23"H.	24"H.	25"H.	26"H.
	<i>From seat pan.</i> 27"H.	28"H.		

19" minimum backrest height when powered recline is selected.

ARM SUPPORTS

Height Adjustable, Removable, Heavy Duty.					Standard
HD Removable Height Package. <i>From seat pan.</i>	8"-12"	FOC	10"-14"		Standard
TB3 Flip-back with Down Posts. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>					FOC
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>		FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>		FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	£50			

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform.					Standard
Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	Standard	13"-19".		FOC
Footplates.	10"W x 8"D Split Footplates.	FOC	10"W x 11"D Split Footplates.		Standard
Calf Support.	Individual Calf Pads. <i>8"W x 6"H.</i>				Standard
	Calf Panel. <i>18"W x 8"H.</i>				FOC

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

70° Heavy Duty, Drop-in, Swing-away Leg Supports. <i>Select lowers and footplates below.</i>							FOC
70° HD	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12"	FOC	13"-18"	Standard	16"-21"	FOC
	Footplate Type.	Angle & Depth Adjustable. <i>Select size below.</i>					Standard
	Footplate Size.	6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
Heavy Duty, Manual Elevating, Swing-away Leg Supports. <i>Flip-up Extruded footplates 8.5"W x 6"D.</i>							FOC
ELEVATING HD	Lower Leg Extensions. <i>From seat pan.</i>	11"-14"	FOC	14"-17"	Standard	12.5"-16.5"	FOC
		15.5"-19.5"	FOC				

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit stealthproducts.com for alternative head supports.

PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>2 power option ready.</i>		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.		Standard
Basic Swing-away.		FOC
Mounting Position.	Left FOC Right	Standard

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	Standard
---	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 706574 - 01/04/25 V1

