

4FRONT 2 - TB3 SEAT / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
	Postcode
Quantum Specialist	
Ⓢ Clients Weight (required for suspension settings)	

BASE MODEL

14 Colour Options
Height Adj. Removable Arm Supports
Centre Mount Foot Platform
TB3 Seating
14" Flat Free Drive Tyres
2" x 9" Flat Free Rear Castors



SRS Suspension System
Smart Traction Control
Q-Logic 3e Control System
2 x 70 ah Batteries
6 Mph (9.65 Km/h) 4-Pole Motors
Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-16"W.
16"-20"W.
20"-22"W.
Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-16D.
16"-20"D.
20"-22"D.
Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>			Standard
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.75"	18.75"	19.75"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts. Select TB3 back</i>	Standard
--	----------

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
--	----------

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	Standard
---	----------

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.
		18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.		
Back Support Height. <i>From seat pan.</i>	19"H.	20"H.	21"H.	22"H.
	23"H.	24"H.	25"H. <i>Not available with 24"W back.</i>	26"H. <i>Not available with 24"W back.</i>
	27"H. <i>Not available with 24"W back.</i>	28"H. <i>Not available with 24"W back.</i>	29"H. <i>Not available with 24"W back.</i>	

19" minimum backrest height when powered recline is selected.

ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i>				Standard
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>				FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan.</i>				FOC
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				FOC
Arm Support Pads.	10" Standard	Standard	10" Waterfall.	FOC
	14" Standard	FOC	14" Waterfall.	FOC

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform.						Standard
Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>						FOC
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16".	Standard	13"-19".	FOC
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>					FOC
	10"W x 8"D Split Footplates.		FOC	10"W x 11"D Split Footplates.		Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>					Standard
Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.						

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-away. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>				FOC
Manual Elevating Swing-Away. <i>Flip-up footplates. Adjustment range 15" to 20". Min. seat width 15".</i>				FOC

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-away Leg Supports. <i>Select lowers and footplates below. Minimum seat width 15".</i>							FOC
60° TO 90°	Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	13.5"-16.5"	Standard	15.5"-19.5"	FOC
	Footplate Type.	Angle & Depth Adjustable. <i>6"W x 8"D.</i>					Standard

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>				FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>				Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>				FOC

PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.				Standard
--	--	--	--	----------

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>2 power option ready.</i>				£300
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>			Standard
	Attendant Control. <i>Hand control mounted at rear.</i>			FOC

HAND CONTROL MOUNTS

Fixed.					Standard
Basic Swing-away.					FOC
Mounting Position.	Left	FOC	Right	Standard	

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.	Standard
4 Mph Program with Lights and indicators.	FOC
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
MOTABILITY ID 715986 - 01/04/25 V1

