

## ASPEN WITH TB FLEX SEAT / POWER LIFT/ POWER TILT / POWER RECLINE

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

### BASE MODEL

4 Colour Options  
Height Adj. Removable Arm Supports  
70° Swing-Away Supports  
Angle Adj. Fold-Forward Back Support  
14" Flat Free Drive Tyres  
3" x 9" Flat Free Front Castors  
ATX Suspension System



Powered Lift  
Powered Tilt  
Powered Recline  
Q-Logic 3e Control System  
2 x 62 ah Batteries  
4 Mph (6 Km/h) Motors  
Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

### FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
--------------------------------	----------

### COLOURS

Back in Black	Standard
Electric Blue	FOC
Silver Bullet	FOC
Red Corvette	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

## SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified.

## SEAT POSITIONING

Powered Seat Lift & Tilt. <i>10" lift &amp; 50° tilt. Uses 2 power options.</i>				Standard
Seat to Floor Height. <i>From seat pan.</i>	17.5"	18"	18.5"	

## STANDARD BACK SUPPORT POSITIONING

Powered Back Recline. <i>Adjustable range 95° to 125°. Uses 1 power option.</i>	Standard
---	----------

## STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat pan.</i>	Standard
---	----------

## ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>				Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				FOC
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with TB Flex.</i>	FOC

## LEG SUPPORTS

70° Swing-away. <i>Adjustment range 12" to 15".</i>	Standard
60° to 90° Tool Adjustable Swing-away. <i>Adjustment range 11" to 16.5".</i>	FOC
Manual Elevating Swing-away. <i>Adjustment range 15" to 20".</i>	FOC
Seat pan to heel measurement	(Leave blank for standard factory)

Adjustment ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit [stealthproducts.com](https://stealthproducts.com) for alternative head supports.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.

Standard (STD)

## NON-EXPANDABLE CONTROLS

Q-Logic 3e Controls. *3 power option ready.*

Standard

Control Method.

Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

## HAND CONTROL MOUNTS

Fixed.

Standard

Basic Swing-away.

FOC

Mounting Position.

Left

FOC

Right

Standard

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. *Lights and indicators not included.*

Standard

6 Mph Program. *With Lights and indicators.*

FOC

## ACCESSORIES

Kerb Climber.

FOC

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**MOTABILITY ID 1001366 - 01/04/25 V1**

