

ASPEN WITH TB FLEX SEAT / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

BASE MODEL

4 Colour Options
 Height Adjustable Removable Armrests
 Height Adj. Removable Arm Supports
 70° Swing-Away Leg Supports
 Angle Adj. Fold-Forward Back Support
 3" x 9" Flat Free Front Castors
 ATX Suspension System



Powered Tilt
 Powered Recline
 QLNES Control System
 2 x 62 ah Batteries
 4 Mph (6 Km/h) Motors
 Max User Weight:
 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
--------------------------------	----------

COLOURS

Back in Black	Standard
Electric Blue	FOC
Silver Bullet	FOC
Red Corvette	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified.

SEAT POSITIONING

Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>				Standard
Seat to Floor Height. <i>From seat pan.</i>	17.5"	18" <i>Changes tilt range to -5° to 25°</i>	18.5" <i>Changes tilt range to -5° to 25°</i>	

STANDARD BACK SUPPORT POSITIONING

Powered Back Recline. <i>Adjustable range 95° to 125°. Uses 1 power option.</i>	Standard
---	----------

STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
---	----------

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>				Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				FOC
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>		Standard
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>		FOC
Manual Elevating Swing-away Legrests. <i>Adjustment range 15" to 20".</i>		FOC
Seat pan to heel measurement	(Leave blank for standard factory)	Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.

Standard

NON-EXPANDABLE CONTROLS

QLNES Controls. *2 power option ready.*

Standard

Control Method. Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

HAND CONTROL MOUNTS

Fixed.

Standard

Basic Swing-away.

FOC

Mounting Position.

Left

FOC

Right

Standard

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. *Lights and indicators not included.*

Standard

4 Mph Program with Lights and indicators.

FOC

6 Mph Program with Lights and indicators.

FOC

ACCESSORIES

Kerb Climber.

FOC

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
MOTABILITY ID 721040 - 01/04/24 V1

