

## EDGE 3 STRETTO - CONTOURED SEAT / 6 MPH

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓢ Clients Weight (required for suspension settings) \_\_\_\_\_

### BASE MODEL

- 14 Colour Options
- Flip-up Arm Supports
- 70° Swing-Away Supports
- Contoured Seat / Fold-Forward Back Support
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- QLNES Control System
- 2 x 62 ah Batteries
- 6 Mph (9.5 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## CONTOURED SEAT PACKAGES

16"W x 16"D	FOC	18"W x 18"D	Standard	20"W x 20"D	FOC
-------------	-----	-------------	----------	-------------	-----

Seat package includes standard head support and auto buckle lap belt.

## SEAT CUSHIONS

Contoured Seat Cushion. <i>136 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Seat cushion dimensions will match seat package selected above.

## SEAT POSITIONING

Static Seat.	Standard
Static Seat to Floor Height. <i>From seat pan.</i>	16.5"      17.5"      18.5"

## ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC
Arm Support Pads	10" Standard      Standard      10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	14" Standard      FOC      14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC

## LEG SUPPORTS

70° Swing-away. <i>Adjustment range 12" to 15".</i>	Standard
Manual Elevating Swing-Away. <i>Inc. calf pads. Adjustment range 15" to 20". Min. seat width 15".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range 10" to 19". 5"W x 11"L Footplates.</i>	FOC
Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i>	

## NON-EXPANDABLE CONTROLS

QLNES Controls. <i>With no power options.</i>	Standard
---	----------

## HAND CONTROL MOUNTS

Fixed.	Standard
Basic swing-away.	FOC
Mounting Position.	Left      FOC      Right      Standard

## SPEED PROGRAM & LIGHTS PACKAGE

---

6 Mph Program with Lights and indicators.

Standard

---

## NOTES & SPECIAL REQUESTS

---



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**Motability ID 706910 - 01/04/25 V1**

