

EDGE 3 STRETTO WITH TB FLEX SEAT

DEALERS DETAILS

| | |
|--------------|-----------|
| Dealers Name | |
| Order Date | |
| Address | |
| Postcode | |
| Phone | Email |
| PO No | Order Ref |

DELIVERY DETAILS

| | |
|---|--|
| Dealers Address | |
| Postcode | |
| Quantum Specialist | |
| Ⓢ Clients Weight (required for suspension settings) | |

BASE MODEL

14 Colour Options
Height Adj. Removable Arm Supports
70° Swing-Away Leg Supports
Angle Adj. Fold-Forward Back Support
12.5" Flat Free Drive Tyres
2" x 5" Flat Free Front Castors



2" x 6" Flat Free Rear Castors
SRS Suspension System
QLNES Control System
2 x 62 ah Batteries
4 Mph (6 Km/h) Motors
Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|--|----------|
| 3" x 12.5" Flat Free Black Tyres. | Standard |
| 3" x 12.5" Pneumatic Black Tyres. | FOC |
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | FOC |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC |

COLOURS

| | | | |
|---------------|----------|-------------------|-----|
| Back in Black | Standard | Pretty in Pink | FOC |
| Green Machine | FOC | Lemon Crush* | FOC |
| Midnight Blue | FOC | Night Train* | FOC |
| Electric Blue | FOC | Evolution Orange* | FOC |
| Blue Jean | FOC | Military Green* | FOC |
| Deep Purple | FOC | Grey Street* | FOC |
| Red Corvette | FOC | Raspberry Beret* | FOC |

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



STANDARD SEAT DIMENSIONS

| | | | | | | | |
|-----------|-----|-----------|----------|-----------|-----|-----------|----------|
| 15" Width | FOC | 16" Width | FOC | 15" Depth | FOC | 16" Depth | FOC |
| 17" Width | FOC | 18" Width | Standard | 17" Depth | FOC | 18" Depth | Standard |
| 19" Width | FOC | 20" Width | FOC | 19" Depth | FOC | 20" Depth | FOC |

STEALTH SEAT CUSHIONS

| | |
|--|----------|
| No Seat Cushion Required | FOC |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i> | Standard |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i> | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

| | |
|--|-------------------|
| Static Seat. | Standard |
| Static Seat to Floor Height. <i>From seat pan.</i> | 16.5" 17.5" 18.5" |

STANDARD BACK SUPPORT POSITIONING *Compatible with standard seat dimensions.*

| | |
|---|----------|
| Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°.</i> | Standard |
|---|----------|

STANDARD BACK SUPPORTS

| | |
|---|----------|
| TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i> | Standard |
|---|----------|

ARM SUPPORTS

| | |
|---|--|
| TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i> | Standard |
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i> | FOC |
| Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i> | FOC |
| Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i> | FOC |
| Arm Support Pads | 10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC |
| | 14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC |

LEG SUPPORTS

| | |
|---|----------|
| 70° Swing-away. <i>Adjustment range 12" to 15".</i> | Standard |
| 60° to 90° Tool Adjustable Swing-away. <i>Adjustment range 11" to 16.5".</i> | FOC |
| Manual Elevating Swing-Away. <i>Inc. calf pads. Adjustment range 15" to 20".</i> | FOC |
| Manually Adjustable Centre Mount Foot Platform. <i>Extension range range 10" to 19". 5"W x 11"L Footplates.</i> | FOC |
| Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i> | |

Adjustment ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops.
Elevating leg supports supplied with calf supports.

STEALTH HEAD SUPPORTS

| | |
|--|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | FOC |

PELVIC SOLUTIONS

| | |
|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--|----------|

NON-EXPANDABLE CONTROLS

| | | |
|-----------------|--|----------|
| QLNES Controls. | | Standard |
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard |
| | Attendant Control. <i>Hand control mounted at rear.</i> | FOC |

HAND CONTROL MOUNTS

| | | | | |
|--------------------|------|----------|-------|----------|
| | | | | |
| Fixed. | | Standard | | |
| Basic swing-away. | | FOC | | |
| Mounting Position. | Left | FOC | Right | Standard |

SPEED PROGRAM & LIGHTS PACKAGE

| | |
|---|----------|
| 4 Mph Program. <i>Lights and indicators not included.</i> | Standard |
|---|----------|

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 707226 - 01/04/25 V1

