

## EDGE 3 STRETTO TB FLEX SEAT / POWER LIFT / POWER TILT

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight (required for suspension settings)	

### BASE MODEL

14 Colour Options  
 Height Adj. Removable Arm Supports  
 70° Swing-Away Leg Supports  
 12.5" Flat Free Drive Tyres  
 2" x 5" Flat Free Front Castors  
 2" x 6" Flat Free Rear Castors  
 SRS Suspension System



Powered Seat Lift  
 Powered Tilt  
 QLNES Control System  
 2 x 40 ah Batteries  
 4 Mph (6 Km/h) Motors  
 Max User Weight:  
 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## STANDARD SEAT DIMENSIONS

15" Width	FOC	16" Width	FOC	15" Depth	FOC	16" Depth	FOC
17" Width	FOC	18" Width	Standard	17" Depth	FOC	18" Depth	Standard
19" Width	FOC	20" Width	FOC	19" Depth	FOC	20" Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Powered Seat Lift & Tilt. <i>10" lift &amp; 50° tilt. Uses 2 power options.</i>	Standard
Powered Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	18" 19" 20"

## STANDARD BACK SUPPORT POSITIONING *Compatible with standard seat dimensions.*

Manual Angle Settable. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
--	----------

## STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
---	----------

## ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i>	FOC
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>	FOC
Arm Support Pads	10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC

## LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
60° to 90° Tool Adjustable Swing-away. <i>Adjustment range 11" to 16.5".</i>	FOC
Manual Elevating Swing-Away. <i>Inc. calf pads. Adjustment range 15" to 20".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Extension range range 10" to 19". 5"W x 11"L Footplates.</i>	FOC
Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i>	

Adjustment ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops.  
Elevating leg supports supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

## NON-EXPANDABLE CONTROLS

QLNES Controls. <i>2 power option ready.</i>		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

## HAND CONTROL MOUNTS

Fixed.		Standard			
Basic swing-away.		FOC			
●●●●●	Mounting Position.	Left	FOC	Right	Standard

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
---	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**Motability ID 706573 - 01/04/25 V1**

