

EDGE 3 STRETTO WITH TB FLEX SEAT / PRO BACK

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight (required for suspension settings)	

BASE MODEL

14 Colour Options
Height Adj. Removable Arm Supports
70° Swing-Away Leg Supports
Angle Adj. Fold-Forward Back Support
12.5" Flat Free Drive Tyres
2" x 5" Flat Free Front Castors



2" x 6" Flat Free Rear Castors
SRS Suspension System
QLNES Control System
2 x 62 ah Batteries
4 Mph (6 Km/h) Motors
Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



STANDARD SEAT DIMENSIONS

15" Width	FOC	16" Width	FOC	15" Depth	FOC	16" Depth	FOC
17" Width	FOC	18" Width	Standard	17" Depth	FOC	18" Depth	Standard
19" Width	FOC	20" Width	FOC	19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations.
Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.	Standard
Static Seat to Floor Height. <i>From seat pan.</i>	16.5" 17.5" 18.5"

STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°.</i>	Standard
---	----------

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	FOC
Active, Coolcore Cover. <i>Solid back, 12 cm contour</i>	Standard
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	FOC
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	FOC
Back Height	50cm. <i>Actual height.</i> Standard 57cm. <i>Actual height.</i> FOC

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i>	FOC
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>	FOC
Arm Support Pads	10" Standard Standard 10" Waterfall. <i>Not with TB Flex.</i> FOC
	14" Standard FOC 14" Waterfall. <i>Not with TB Flex.</i> FOC

LEG SUPPORTS

70° Swing-away. <i>Adjustment range 12" to 15".</i>	Standard
60° to 90° Tool Adjustable Swing-away. <i>Adjustment range 11" to 16.5".</i>	FOC
Manual Elevating Swing-Away. <i>Inc. calf pads. Adjustment range 15" to 20".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Extension range 10" to 19". 5"W x 11"L Footplates.</i>	FOC
Enter seat pan to heel measurement.	<i>If no measurement is supplied the standard factory length option will be fitted.</i>

Adjustment ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops.
Elevating leg supports supplied with calf supports.

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

NON-EXPANDABLE CONTROLS

QLNES Controls.		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

HAND CONTROL MOUNTS

Fixed.				Standard
Basic swing-away.				FOC
Mounting Position.	Left	FOC	Right	Standard

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
---	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 1000838 - 01/04/25 V1

