

MOTABILITY ID 710381

EDGE 3 WITH TB3 SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 14 Colour Options
- Height Adj. Removable Arm Supports
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- Q-Logic 3e Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|--|----------|
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | Standard |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC |

COLOURS

| | | | |
|---------------|----------|-------------------|-----|
| Back in Black | Standard | Pretty in Pink | FOC |
| Green Machine | FOC | Lemon Crush* | FOC |
| Midnight Blue | FOC | Night Train* | FOC |
| Electric Blue | FOC | Evolution Orange* | FOC |
| Blue Jean | FOC | Military Green* | FOC |
| Deep Purple | FOC | Grey Street* | FOC |
| Red Corvette | FOC | Raspberry Beret* | FOC |

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

| |
|---------------|
| 12"-16"W. |
| 14"-16"W. |
| 16"-20"W |
| 20"-22"W |
| Set width at: |

SEAT PACKAGES - SEAT DEPTH

| |
|---------------|
| 12"-16"D. |
| 14"-16"D. |
| 16"-20"D |
| 20"-22"D |
| Set depth at: |

STEALTH SEAT CUSHIONS

| | |
|--|----------|
| No Seat Cushion Required | FOC |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i> | Standard |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i> | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

| | |
|--|-----------------------------|
| Static Seat. | Standard |
| Static Seat to Floor Height. <i>From seat pan.</i> | 16.5" 17.5" 18.5" |

BACK SUPPORT SPECIFICATIONS

| | |
|---|----------|
| TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i> | Standard |
| Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i> | FOC |

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

| | |
|--|----------|
| Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i> | Standard |
|--|----------|

TB3 BACK SUPPORT POSTIONING

| | |
|---|----------|
| TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i> | Standard |
| TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports, Not suitable for occupied transport.</i> | FOC |

BACK CANE ONLY SUPPORT

| | |
|---|--|
| Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i> | Standard |
| Back Cane Only Support Positioning. Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i> | Standard |
| Back Cane Height. <i>From seat pan.</i> | 16"H. 17"H. 18"H. 19"H. 20"H. 21"H. 22"H. 23"H. 24"H. 25"H. |

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Ignore Back Support Dimensions if omit back support is selected in TB3 or Back Canes sections. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

| | | | | |
|---|------------------------------|------------------------------|-----------|-----------|
| Backrest Support Width Package. | | 12"-14"W. | 14"-16"W. | 16"-20"W. |
| | | 18"-22"W. | 20"-24"W. | |
| Back Support Width. | 12"W. | 14"W. | 16"W. | 17"W. |
| | 18"W. | 19"W. | 20"W. | 21"W. |
| | 22"W. | 24"W. | | |
| Back Support Height. <i>From seat pan.</i> | 14"H. <i>12"W back only.</i> | 15"H. <i>12"W back only.</i> | 16"H. | 17"H. |
| | 18"H. | 19"H. | 20"H. | 21"H. |
| | 22"H. | 23"H. | 24"H. | 25"H. |
| | 26"H. | 27"H. | 28"H. | 29"H. |
| | 30"H. | 31"H. | 32"H. | |

ARM SUPPORTS

| | | |
|--|---|-----|
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i> | Standard | |
| Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i> | FOC | |
| Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i> | FOC | |
| TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i> | FOC | |
| Arm Support Pads. | 10" Standard Standard 10" Waterfall. <i>Not with TB3 Flip-back.</i> | FOC |
| | 14" Standard FOC 14" Waterfall. <i>Not with TB3 Flip-back.</i> | FOC |

CENTRE MOUNTED LEG SUPPORTS

| | | | | |
|--|--|---|---|---|
| Base Mounted Flip-up Foot Platform. <i>Mounted to chassis.</i> | FOC | | | |
| Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i> | FOC | | | |
| Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i> | Standard | | | |
| Lower Leg Extensions. <i>From seat pan.</i> | 8"-12". <i>Paediatric centre mount only.</i> FOC | 10"-16". <i>Manually adjustable only.</i> FOC | 13"-19". <i>Manually adjustable only.</i> STD | |
| Foot Supports. | 10"W x 6"D Split Foot Supports. <i>Paediatric centre mount only.</i> | | | FOC |
| | 10"W x 8"D Split Foot Supports. | | FOC | 10"W x 11"D Split Foot Supports. Standard |
| Calf Support. | Individual Calf Pads. <i>6"W x 6"H.</i> Standard | Calf Panel. <i>18"W x 8"H.</i> | FOC | |

INDIVIDUALLY MOUNTED LEG SUPPORTS

| | |
|---|-----|
| 70° Swing-Away. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i> | FOC |
| Manual Elevating Swing-Away. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i> | FOC |

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

| | | | | | | |
|------------|---|---------|-----|-------------|----------|-------------|
| 60° TO 90° | 60° to 90° Tool Adjustable Swing-Away. <i>Angle & depth adjustable footplates 6"W x 8"D. Select lowers below. Minimum seat width 15".</i> | | | | | FOC |
| | Lower Leg Extensions. <i>From seat pan.</i> | 11"-13" | FOC | 13.5"-16.5" | Standard | 15.5"-19.5" |

STEALTH HEAD SUPPORTS

| | |
|--|----------|
| UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i> | FOC |
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | FOC |

Please visit stealthproducts.com for alternative head supports.

PELVIC BELTS

| | |
|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--|----------|

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth.

| | | |
|--|--|----------|
| Q-Logic 3e Controls. <i>With no power options.</i> | Standard | |
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard |
| | Attendant Control. <i>Hand control mounted at rear.</i> | FOC |

HAND CONTROL MOUNTS

| | |
|--------------------|------------|
| Fixed. | Standard |
| Basic Swing-away. | FOC |
| Mounting Position. | Left Right |

SPEED PROGRAM & LIGHTS PACKAGE

| | |
|---|----------|
| 4 Mph Program with Lights and indicators. | Standard |
| 6 Mph Program with Lights and indicators. | FOC |
| LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i> | Standard |
| Omit LED Fender Lights. | FOC |

ACCESSORIES

| | |
|---|-----|
| Prepare Chassis for Dahl Docking System. <i>Dahl system not included.</i> | FOC |
|---|-----|

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 710381 - 01/04/25 V1

