

## FUSION E ACTIVE CONTOUR 6 MPH

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

### BASE MODEL

Height Adj. Removable Arm Supports  
 70° Swing-Away Leg Supports  
 Adjustable Seat Width & Depth Package  
 TRU-Comfort 2 Seat Cushion  
 Stealth Headrest  
 Fold-forward Back with Stealth Pro Back Support  
 14" Pneumatic Drive Wheel  
 10" Pneumatic Castors



Rear Suspension System  
 Black Frame  
 Yellow Wheel  
 Accent Inserts  
 VR2 Control System  
 2 x 50 ah Batteries  
 Max User Weight:  
 23 st 6lbs. (150 kg)  
 Stealth Pro Back Support

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### SEAT PACKAGE

18"W x 20"D. *Seat width can be adjusted between 16-23". Seat depth can be adjusted between 16-21".*

Standard

### SEAT CUSHIONS

2" Foam Seat Cushions. *Supplied with 1 x 16"W x 18"D & 1 x 20" x 22".*

Standard

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.*

Standard

Enter TC2 cushion size here

**WIDTH**

**DEPTH**

*Refer to cushion order forms for size availability.*

### BACK SUPPORT POSITIONING

Angle Adjustable, Fold-Forward. *3° - 33° angle adjustment in 4 steps of 10°.*

Standard

## STEALTH PRO BACK SUPPORTS

Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>				Standard
Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>				FOC
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				FOC
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>				FOC
Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC
Back Width	16"W.		16.5"W.	17"W.
	18"W.		19"W.	20"W.
	21"W.		22"W.	

## SEAT POSITIONING

Static Seat. <i>18.5" front seat to floor and 17" rear seat to floor height from seat pan.</i>	Standard
------------------------------------------------------------------------------------------------	----------

## ARM SUPPORTS

Height Adjustable, Removable. <i>Includes 15"L x 3"W arm pads. Adjustment range 10.5" to 14.5" from seat pan.</i>	Standard
-------------------------------------------------------------------------------------------------------------------	----------

## LEG SUPPORTS

70° Swing-away Leg Supports. <i>Adjustment range 13" to 16" from seat pan.</i>	Standard
--------------------------------------------------------------------------------	----------

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--------------------------------------------	----------

## NON-EXPANDABLE CONTROLS

VR2 Hand Control.	Standard
-------------------	----------

## HAND CONTROL MOUNTS

Right Side with Fixed Position Mount. <i>Can be moved to left side.</i>	Standard
-------------------------------------------------------------------------	----------

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

## SPEED PACKAGE

6 Mph with Lights and Indicators.	Standard
-----------------------------------	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**Motability ID 706355 - 01/04/25 V1**

