

KOZMO PORTABLE

DEALERS DETAILS

Dealers Name	
<hr/>	
Order Date	
<hr/>	
Address	
<hr/>	
	Postcode
<hr/>	<hr/>
Phone	Email
<hr/>	<hr/>
PO No	Order Ref
<hr/>	<hr/>

DELIVERY DETAILS

Dealers Address	
<hr/>	
<hr/>	
Postcode	
<hr/>	
Quantum Specialist	
<hr/>	
Clients Weight	
<hr/>	

BASE MODEL

- 6 Colour Options
- Height Adj. Removable Arm Supports
- Flip-up Foot Platform
- Padded Fold-Forward Back Support
- 9" Flat Free Drive Tyres
- 6" Flat Free Front Castors
- Dismantles into 4 Pieces



- 12" - 18" Seat Width
- 12" - 18" Seat Depth
- GC3 Control System
- 2 x 18 ah Batteries
- 3.6 Mph (5.79 Km/h) Motors
- Max User Weight:
17 st 8lbs. (113 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

COLOURS

Ruby Red	Standard	Rose Quartz	FOC
Pearl White	FOC	Citrine Yellow	FOC
Sapphire Blue	FOC		
Amber Orange	FOC		

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

12-16" Width	FOC
16-18" Width	FOC

Set width at:

SEAT DEPTH

12-16" Depth	FOC
16-18" Depth	FOC

Set depth at:

SEAT HEIGHT - SHORT SEAT POST

16"H <i>Seat pan to floor.</i>	FOC
17"H <i>Seat pan to floor.</i>	FOC

SEAT HEIGHT - STANDARD SEAT POST

19"H <i>Seat pan to floor.</i>	FOC
20"H <i>Seat pan to floor.</i>	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

BACK SUPPORTS

Visco Padded Back Support.	Standard
Back Height. <i>From Seat pan to top of back.</i>	16" 17" 18" 19" <i>Not available with 17" seat width.</i>

ARM SUPPORTS

Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>	Standard
Flip-up Armrests Mounted to Back Canes. <i>Range of adjustment 2".</i>	FOC
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>	FOC
Arm Support Pads	10" Standard Standard 10" Waterfall. FOC
	14" Standard FOC 14" Waterfall. FOC

LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis.</i>	FOC
70° Swing-away. <i>Adjustment range 12" to 15". Min. seat width 16".</i>	FOC
Manual Elevating Swing-away. <i>Adjustment range 15" to 20". Min. seat width 16".</i>	FOC
Paediatric Fixed Foot Platform. <i>Mounted to chassis.</i>	FOC

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

NON-EXPANDABLE CONTROLS

GC3 Controls.	Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i> Standard

HAND CONTROL MOUNTS

Fixed.	Standard
Advanced Swing-away. <i>Chain driven. Swings in & out with height, tension, and pivot adjustment.</i>	FOC
Mounting Position.	Left FOC Right Standard

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 713666 - 01/04/25 V1

