



MOTABILITY ID 1002475

OUTBACK WITH CONTOURED SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 2 Colour Options
- Flip-up Arm Supports
- Centre Mount Foot Platform
- Contoured Seating
- 15" Pneumatic Drive Tyres
- 4 Motor Independent Drive



- SRS Suspension System
- Q-Logic 3e Control System
- 2 x 95 ah Batteries
- 6 Mph (9.65 Km/h) 4-Pole Motors
- Max User Weight: 25 st (159 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

COLOURS

Cadet Matte	Standard
Tan Matte	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



CONTOURED SEAT PACKAGES

16"W x 16"D	FOC	20"W x 16"D*	FOC
16"W x 18"D*	FOC	20"W x 18"D	FOC
16"W x 20"D*	FOC	20"W x 20"D	FOC
18"W x 16"D	FOC	20"W x 22"D*	FOC
18"W x 18"D	Standard	22"W x 18"D*	FOC
18"W x 20"D	FOC	22"W x 20"D	FOC
18"W x 22"D*	FOC	22"W x 22"D	FOC

*Extended lead time. All seat packages above include a contoured seat cushion and a standard head support.

STEALTH SEAT CUSHIONS

No Seat Cushion Required.	FOC
Contoured Seat Cushion. <i>136 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.	Standard
Seat to Floor Height. <i>From seat pan.</i>	19.5" 20.5" 21.5"

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes 14"L contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75" - 9.75" from seat pan. 10" pads only.</i>	FOC
Arm Support Pads.	10" Standard STD 10" Waterfall. <i>Only with single post or paediatric arms.</i> FOC
	14" Standard FOC 14" Waterfall. <i>Only with single post or paediatric arms.</i> FOC

CENTRE MOUNTED LEG SUPPORTS

Basic Centre Mount Foot Platform	Standard
Lower Leg Extensions. <i>From seat pan.</i>	8"-10". <i>Only available with 16" seat width.</i> STD 10"-14". Standard 14"-19". FOC
Foot Supports	10"W x 9"D Split Foot Supports. Standard 10"W x 12"D Split Foot Supports. FOC
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i> Standard
	Individual Large Calf Pads. <i>6"W x 8"H. Only with long lower ext.</i> FOC

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

70° HD Drop-in Swing-away. <i>Includes heel loops.</i>	FOC	
70°	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12" FOC 13"-18" Standard 16"-21" FOC
	Foot Support Size.	6"W x 8"D Standard 6"W x 11"D FOC 8"W x 11"D FOC
	Accessories.	Wheel Bumpers FOC

PELVIC BELTS

Standard Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.	Standard
	38".	FOC 70". FOC
	50".	Standard 80". FOC
	60".	FOC 90". FOC

EXPANDABLE CONTROLS Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls.	Standard
Control Method. Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard

HAND CONTROL MOUNTS

Fixed.	Standard
Basic Swing-away.	FOC
Mounting Position. Left FOC Right Standard	

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.	Standard
---	----------

NOTES & SPECIAL REQUESTS

**OUTBACK WITH
CONTOURED SEATING**



BASE PRICE

COST OPTIONS

TOTAL

3



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 1002475 - 01/04/25 V1

