

Q1450 WITH CONTOURED SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

5 Colour Options
Height Adj. Removable Arm Supports
Flip-up Centre Mount Foot Platform
Angle Adjustable Back Support
Seat Cushion
3" x 14" Flat Free Drive Tyres
2" x 9" Flat Free Rear Castors



1.5" Anti-tip Wheels
Swing-arm Suspension System
VR2 Control System
2 x 70 ah Batteries
4 Mph (6 Km/h) Motors
Max User Weight:
42 st 8lbs. (272 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Grey Tyres.

Standard

REAR CASTORS

2" x 9" Flat Free Grey Tyres.

Standard

COLOURS

Back in Black	Standard
Electric Blue	FOC
Red Corvette	FOC
Evolution Orange	FOC
Pretty in Pink	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEATING OPTIONS

Pinchless Hinge High Back Captain's Seat. <i>Includes a head support and contoured seat cushion.</i>	22"W-20"D.	FOC	22"W-22"D.	FOC	22"W-24"D.	Standard
	24"W-20"D.	FOC	24"W-22"D.	FOC	24"W-24"D.	FOC
HD High Back Captain's Seat. <i>Includes a head support, lumbar cut-out, static hinges, contoured seat cushion.</i>			26"W-20"D.	FOC	26"W-22"D	FOC
Pinchless Hinge High Back Contoured Seat. <i>Includes a head support and solid seat pan.</i>					22"W-20"D	FOC

STEALTH SEAT CUSHIONS

Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	STD
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may only ordered with the Contoured Pinchless Hinge High Back seat and the solid seat pan option.

SEAT POSITIONING

Static Seat.				Standard
Static Seat to Floor Height. <i>From seat pan.</i>		17.5"	18.5"	19.5"

ARM SUPPORTS

Bariatric Single Post Tool Height Adjustable Removable. <i>Includes Bariatric full length arm pads.</i>					Standard
Single Post HD Height Adjustable Removable.					FOC
HD & Bariatric Height Package. <i>From seat pan.</i>		8"-12".	FOC	10"-14".	Standard
Single Post HD Arm Support Pads.	10" Standard. <i>Not available with Bariatric.</i>				FOC
	14" Standard. <i>Not available with Bariatric.</i>				Standard
	10" Waterfall. <i>Not available with Bariatric.</i>				FOC
	14" Waterfall. <i>Not available with Bariatric.</i>				FOC
Flip-up Arm Supports <i>Includes arm pads.</i>	Full Length, Flip-up, Height Adjustable, Removable. <i>Required with pinchless hinge high back. Not available with HD seat</i>				FOC
	Desk Length, Flip-up, Height Adjustable. <i>Required with pinchless hinge high back. Not available with HD seat.</i>				FOC

CENTRE MOUNTED LEG SUPPORTS

17"W x 13"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis.</i>	Standard
17"W x 16"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis.</i>	FOC

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away.						FOC
70° HD Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
70° HD Foot Supports.	6"W x 8"D.	Standard	6"W x 11"D.	FOC	8"W x 11"D.	FOC
Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to foot supports. Individual leg supports are supplied with heel loops.						
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT	

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard
Auto Buckle Lap Belt Length.	60".	Standard	70".	FOC	80".	FOC
	90".	FOC				

NON-EXPANDABLE CONTROLS

VR2 Controls. <i>With no power options.</i>						Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>					Standard

HAND CONTROL MOUNTS

Fixed.						Standard
Basic Swing-away.						FOC
Mounting Position.	Left		Right			

SPEED PROGRAM

4 Mph Program						Standard
---------------	--	--	--	--	--	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 713664 - 01/04/25 V1

