

## Q6 EDGE 2.0 HD WITH TB3 SEATING

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

### BASE MODEL

13 Colour Options  
Height Adj. Removable Arm Supports  
Centre Mount Foot Platform  
TB3 Seating  
14" Flat Free Drive Tyres  
2" x 6" Flat Free Front Castors



2" x 6" Flat Free Rear Castors  
SRS Suspension System  
Q-Logic 3e Control System  
2 x 62 ah Batteries  
4 Mph (6 Km/h) 4-Pole Motors  
Max User Weight: 28 st 5lbs. (181.0 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Deep Purple	FOC	Military Green*	FOC
Red Corvette	FOC	Grey Street*	FOC

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

20"-24"W

Set width at:

## SEAT PACKAGES - SEAT DEPTH

20"-24"D

Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	STD
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>181.4 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>181.4 Kg weight limit.</i>	FOC
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	FOC
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>226.7 Kg weight limit.</i>	FOC
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Static Seat.	Standard
Static Seat to Floor Height. <i>From seat pan.</i>	16" 17" 18"

## BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
--	----------

## TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports or powered recline. Not suitable for occupied transport.</i>	FOC

## BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>						Standard
Back Cane Only Support Positioning.		Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>				Standard
Back Cane Height. <i>From seat pan.</i>	16"H.	17"H.	18"H.	19"H.	20"H.	
	21"H.	22"H.	23"H.	24"H.	25"H.	

## BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.		16"-20"W.		20"-24"W.
Back Support Width.	16"W.	17"W.	18"W.	19"W.
	20"W.	21"W.	22"W.	24"W.
Back Support Height. <i>From seat pan.</i>	16"H.	17"H.	18"H.	19"H.
	20"H.	21"H.	22"H.	23"H.
	24"H.	25"H.		

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

Height Adjustable, Removable, Heavy Duty.				Standard
HD Removable Height Package. <i>From seat pan.</i>	8"-12"H.		FOC	
	10"-14"H.		Standard	
TB3 Flip-back with Down Posts. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				FOC
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	£50		

## CENTRE MOUNTED LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis.</i>				Standard
Manually Adjustable Centre Mount Foot Platform.				FOC
Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	Standard	13"-19".	FOC
Foot Supports.	10"W x 8"D Split Foot Supports.	FOC	10"W x 11"D Split Foot Supports.	Standard
Calf Support.	Individual Calf Pads. <i>8"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>	FOC

Certain configurations may result in reduced ground clearance below foot supports. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

70° Heavy Duty, Drop-in, Swing-away. <i>Select lowers and footplates below.</i>							FOC
70° HD	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12"	FOC	13"-18"	Standard	16"-21"	FOC
	Foot Support Type.	Angle & Depth Adjustable. <i>Select size below.</i>					Standard
	Foot Support Size.	6"W x 8"D	STD	6"W x 11"D	FOC	8"W x 11"D	FOC

## INDIVIDUALLY MOUNTED ELEVATING LEG SUPPORTS

Heavy Duty, Manual Elevating, Swing-away. <i>Flip-up Extruded foot supports 8.5"W x 6"D.</i>							FOC
ELEVATING HD	Lower Leg Extensions.	11"-14"	FOC	14"-17"	Standard	12.5"-16.5"	FOC
	<i>From seat pan.</i>	15.5"-19.5"	FOC				

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to foot supports. Individual foot supports supplied with heel loops. Elevating leg supports supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	FOC
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	Standard
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

## EXPANDABLE CONTROLS

Q-Logic 3e. <i>No power options.</i>	Standard
Control Method. Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
Attendant Control. <i>Hand control mounted at rear.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.			Standard	
Basic Swing-away.			FOC	
Mounting Position.	Left	FOC	Right	Standard

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	Standard
---	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
 Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**Motability ID 706571 - 01/04/25 V1**

