

Q6 EDGE 2.0 WITH TB FLEX SEATING

DEALERS DETAILS

| | |
|--------------|-----------|
| Dealers Name | |
| Order Date | |
| Address | |
| | Postcode |
| Phone | Email |
| PO No | Order Ref |

DELIVERY DETAILS

| | |
|--------------------|----------|
| Dealers Address | |
| | |
| | Postcode |
| Quantum Specialist | |
| Clients Weight | |

BASE MODEL

- 12 Colour Options
- Height Adj. Removable Arm Supports
- 70° Swing-Away Leg Supports
- Angle Adj. Fold-Forward Back Support
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- QLNES Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|---------------------------------|----------|
| 3" x 14" Flat Free Black Tyres. | Standard |
| 3" x 14" Pneumatic Black Tyres. | FOC |

COLOURS

| | | | |
|----------------|----------|-------------------|-----|
| Back in Black | Standard | Red Corvette | FOC |
| Green Machine | FOC | Lemon Crush* | FOC |
| Midnight Blue | FOC | Night Train* | FOC |
| Electric Blue | FOC | Evolution Orange* | FOC |
| Pretty in Pink | FOC | Military Green* | FOC |
| Deep Purple | FOC | Grey Street* | FOC |

* Extended lead-time

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

| | | | |
|-----------|-----|-----------|----------|
| 15" Width | FOC | 16" Width | FOC |
| 17" Width | FOC | 18" Width | Standard |
| 19" Width | FOC | 20" Width | FOC |

SEAT DEPTH

| | | | |
|-----------|-----|-----------|----------|
| 15" Depth | FOC | 16" Depth | FOC |
| 17" Depth | FOC | 18" Depth | Standard |
| 19" Depth | FOC | 20" Depth | FOC |

STEALTH SEAT CUSHIONS

| | |
|--|----------|
| No Seat Cushion Required | FOC |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i> | Standard |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

| | |
|--|-----------------------|
| Static Seat. | Standard |
| Static Seat to Floor Height. <i>From seat pan.</i> | 16" 17" 18" |

STANDARD BACK SUPPORT POSITIONING

| | |
|---|----------|
| Manual Angle Settable, Fold-Forward. <i>Adjustment range 95° to 125°.</i> | Standard |
|---|----------|

STANDARD BACK SUPPORTS

| | |
|--|----------|
| TRU-Comfort 2. <i>Plastic shell, 4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i> | Standard |
|--|----------|

ARM SUPPORTS

| | |
|---|---|
| TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i> | Standard |
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16"</i> | FOC |
| Arm Support Pads | 10" Standard Standard 10" Waterfall. FOC |
| | 14" Standard FOC 14" Waterfall. FOC |

LEG SUPPORTS

| | |
|---|----------|
| 70° Swing-away. <i>Adjustment range 12" to 15"</i> | Standard |
| Base Mounted Foot Platform. <i>Mounted to chassis.</i> | FOC |
| Manual Elevating Swing-Away. <i>Adjustment range 15" to 20"</i> | FOC |
| Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range range 10" to 19" W x 11" L Footplates.</i> | FOC |
| Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i> | |

Adjustment ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

STEALTH HEAD SUPPORTS

| | |
|--|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | FOC |

PELVIC SOLUTIONS

| | |
|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--|----------|

NON-EXPANDABLE CONTROLS

| | | |
|-----------------|--|----------|
| QLNES Controls. | | Standard |
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard |
| | Attendant Control. <i>Hand control mounted at rear.</i> | FOC |

HAND CONTROL MOUNTS

| | | |
|--------------------|------|----------|
| Fixed. | | Standard |
| Basic swing away. | | FOC |
| Mounting Position. | Left | Right |

SPEED PROGRAM & LIGHTS PACKAGE

| | |
|---|----------|
| 4 Mph Program. <i>Lights and indicators not included.</i> | Standard |
|---|----------|

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
Motability ID 1002463 - 01/04/25

