

## Q6 EDGE 2.0 WITH TB FLEX SEATING / PRO BACK

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

### BASE MODEL

12 Colour Options  
Height Adj. Removable Arm Supports  
70° Swing-Away Leg Supports  
Angle Adj. Fold-Forward Back Support  
14" Flat Free Drive Tyres  
2" x 6" Flat Free Front Castors



2" x 6" Flat Free Rear Castors  
ATX Suspension System  
QLNES Control System  
2 x 62 ah Batteries  
4 Mph (6 Km/h) Motors  
Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

### COLOURS

Back in Black	Standard	Red Corvette	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Pretty in Pink	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC

\* Extended lead-time

Please download our free PQ-365 app to view colours and other product features in more detail.



### SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

### SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

## SEAT POSITIONING

Static Seat.	Standard
Static Seat to Floor Height. <i>From seat pan.</i>	16" 17" 18"

## STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adjustment range 95° to 125°.</i>	Standard
---	----------

## STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	FOC
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>	Standard
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	FOC
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	FOC
Back Height 50 cm. <i>Actual height.</i> Standard 57 cm. <i>Actual height.</i>	FOC

The back support width will match the seat width selected in seat dimensions.

## ARM SUPPORTS

TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>	Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	FOC
Arm Support Pads 10" Standard Standard 10" Waterfall.	FOC
14" Standard FOC 14" Waterfall.	FOC

## LEG SUPPORTS

70° Swing-away. <i>Adjustment range 12" to 15".</i>	Standard
Base Mounted Foot Platform. <i>Mounted to chassis.</i>	FOC
Manual Elevating Swing-Away. <i>Adjustment range 15" to 20".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range range 10" to 19"W x 11"L Footplates.</i>	FOC
Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i>	

Adjustment ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

## NON-EXPANDABLE CONTROLS

QLNES Controls.		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

## HAND CONTROL MOUNTS

Fixed.		Standard
Basic swing away.		FOC
Mounting Position.	Left	Right

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
---	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**Motability ID 1002468 - 01/04/25**

