

MOTABILITY ID 721029

Q6 EDGE 2.0 WITH TB3 SEAT / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 12 Colour Options
- Height Adj. Removable Arm Supports
- Centre Mount Foot Platform
- TB3 Seating
- Auto Buckle Lap Belt
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- QLNES Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) 4-Pole Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Red Corvette	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Pretty in Pink	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-16"W.
14"-16"W.
16"-20"W.
20"-22"W.
Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-16"D.
14"-16"D.
16"-20"D.
20"-22"D.
Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required.	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	Standard
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5" 18.5" 19.5"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts. Select TB3 back options below.</i>	Standard
---	----------

TB3 BACK SUPPORT POSITIONING

TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	Standard
---	----------

BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB3 Back Support System.

Backrest Support Width Package.	12"-14"W.	14"-16"W.	16"-20"W.
	18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	17"W.
	18"W.	19"W.	21"W.
	22"W.	24"W.	
Back Support Height.	19"H.	20"H.	22"H.
	23"H.	24"H.	26"H. <i>Not available with 24"W back.</i>
	<i>From seat pan.</i> 27"H. <i>Not available with 24"W back.</i>	28"H. <i>Not available with 24"W back.</i>	29"H. <i>Not available with 24"W back.</i>

19" minimum backrest height when powered recline is selected.

ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>						Standard
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>						FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>						FOC
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>						FOC
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with TB3 Flip-back</i>			FOC
	14" Standard	FOC	14" Waterfall. <i>Not with TB3 Flip-back</i>			FOC

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>						Standard
Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>						FOC
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16".	Standard	13"-19".	FOC
Foot Supports.	10"W x 6"D Split Foot Supports. <i>Paediatric centre mount only.</i>					FOC
	10"W x 8"D Split Foot Supports.	FOC	10"W x 11"D Split Foot Supports.			Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>					Standard

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>						FOC
Manual Elevating Swing-Away. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>						FOC

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-Away. <i>Angle & depth adjustable footplates 6"W x 8"D. Select lowers below. Minimum seat width 15"</i>						FOC
Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	13.5"-16.5"	Standard	15.5"-19.5"	FOC

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>						FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>						Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>						FOC

Please visit stealthproducts.com for alternative head supports.

NON-EXPANDABLE CONTROLS

QLNES Controls. <i>2 power option ready.</i>						Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>					Standard
	Attendant Control. <i>Hand control mounted at rear.</i>					FOC

HAND CONTROL MOUNTS

Fixed.			Standard
Basic swing-away.			FOC
Mounting Position.	Left	Right	

SPEED PROGRAM & LIGHTS PACKAGE

2-Pole Motors & 4 Mph Program with Lights and indicators.	Standard
2-Pole Motors & 6 Mph Program with Lights and indicators.	FOC

ACCESSORIES

Prepare Chassis for Dahl Docking System. <i>Dahl system not included.</i>	FOC
---	-----

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
MOTABILITY ID 721029 - 01/04/25 V1

