

## Q6 EDGE 2.0 WITH TB3 SEAT / iLEVEL / POWER TILT / POWER RECLINE / POWER FOOTBOARD

### DEALERS DETAILS

|              |           |
|--------------|-----------|
| Dealers Name |           |
| Order Date   |           |
| Address      |           |
| Postcode     |           |
| Phone        | Email     |
| PO No        | Order Ref |

### DELIVERY DETAILS

|                    |  |
|--------------------|--|
| Dealers Address    |  |
| Postcode           |  |
| Quantum Specialist |  |
| Ⓢ Clients Weight   |  |

### BASE MODEL

12 Colour Options  
Height Adj. Removable Arm Supports  
Centre Mount Foot Platform  
TB3 Seating  
Auto Buckle Lap Belt  
14" Flat Free Drive Tyres  
2" x 6" Flat Free Front Castors



2" x 6" Flat Free Rear Castors  
ATX Suspension System  
Q-Logic 3e Control System  
2 x 62 ah Batteries  
4 Mph (6 Km/h) 4-Pole Motors  
Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

|                                                                            |          |
|----------------------------------------------------------------------------|----------|
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | Standard |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC      |

### COLOURS

|                |          |                   |     |
|----------------|----------|-------------------|-----|
| Back in Black  | Standard | Red Corvette      | FOC |
| Green Machine  | FOC      | Lemon Crush*      | FOC |
| Midnight Blue  | FOC      | Night Train*      | FOC |
| Electric Blue  | FOC      | Evolution Orange* | FOC |
| Pretty in Pink | FOC      | Military Green*   | FOC |
| Deep Purple    | FOC      | Grey Street*      | FOC |

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

|               |
|---------------|
| 14"-16"W.     |
| 16"-20"W      |
| 20"-22"W      |
| Set width at: |

## SEAT PACKAGES - SEAT DEPTH

|               |
|---------------|
| 14"-16"D.     |
| 16"-20"D      |
| 20"-22"D      |
| Set depth at: |

Minimum seat size of 14"W x 14"D required for powered centre mount foot platform when iLevel is selected.

## STEALTH SEAT CUSHIONS

|                                                                                                              |          |
|--------------------------------------------------------------------------------------------------------------|----------|
| No Seat Cushion Required                                                                                     | FOC      |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>                  | Standard |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>                          | FOC      |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC      |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

|                                                                                                                       |       |       |       |          |
|-----------------------------------------------------------------------------------------------------------------------|-------|-------|-------|----------|
| iLevel Powered Seat Lift & Tilt. 12" lift & 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options. |       |       |       | Standard |
| Powered iLevel Lift & Tilt Seat to Floor Height. From seat pan.                                                       | 17.5" | 18.5" | 19.5" |          |

## BACKREST SPECIFICATIONS

|                                                                          |          |
|--------------------------------------------------------------------------|----------|
| TRU-Balance 3 Back Support System. <i>Select TB3 back options below.</i> | Standard |
|--------------------------------------------------------------------------|----------|

Include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

|                                                                                                                               |          |
|-------------------------------------------------------------------------------------------------------------------------------|----------|
| Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i> | Standard |
|-------------------------------------------------------------------------------------------------------------------------------|----------|

## TB3 BACK SUPPORT POSTIONING

|                                                                                       |          |
|---------------------------------------------------------------------------------------|----------|
| TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i> | Standard |
|---------------------------------------------------------------------------------------|----------|

## BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB3 Back Support System.

| Backrest Support Width Package.            |       | 12"-14"W. | 14"-16"W. | 16"-20"W. |
|--------------------------------------------|-------|-----------|-----------|-----------|
|                                            |       | 18"-22"W. | 20"-24"W. |           |
| Back Support Width.                        | 12"W. | 14"W.     | 16"W.     | 17"W.     |
|                                            | 18"W. | 19"W.     | 20"W.     | 21"W.     |
|                                            | 22"W. | 24"W.     |           |           |
| Back Support Height. <i>From seat pan.</i> | 19"H. | 20"H.     | 21"H.     | 22"H.     |
|                                            | 23"H. | 24"H.     | 25"H.     | 26"H.     |
|                                            | 27"H. | 28"H.     | 29"H.     | 30"H.     |
|                                            | 31"H. | 32"H.     |           |           |

## ARM SUPPORTS

|                                                                                                              |              |          |                                              |          |
|--------------------------------------------------------------------------------------------------------------|--------------|----------|----------------------------------------------|----------|
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>           |              |          |                                              | Standard |
| Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i> |              |          |                                              | FOC      |
| Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>        |              |          |                                              | FOC      |
| TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>             |              |          |                                              | FOC      |
| Arm Support Pads.                                                                                            | 10" Standard | Standard | 10" Waterfall. <i>Not with TB3 Flip-back</i> | FOC      |
|                                                                                                              | 14" Standard | FOC      | 14" Waterfall. <i>Not with TB3 Flip-back</i> | FOC      |

## CENTRE MOUNTED LEG SUPPORTS

|                                                                                                              |                                         |     |          |          |
|--------------------------------------------------------------------------------------------------------------|-----------------------------------------|-----|----------|----------|
| Power Elevating Articulating Centre Mount Foot Platform. <i>Minimum seat width 14". Uses 1 power option.</i> |                                         |     |          | Standard |
| Lower Leg Extensions.                                                                                        | 10"-16".<br><i>From seat pan.</i>       | FOC | 13"-19". | Standard |
| Foot Supports.                                                                                               | 10"W x 8"D Split Foot Supports.         |     |          | FOC      |
|                                                                                                              | 10"W x 11"D Split Foot Supports.        |     |          | Standard |
| Calf Support.                                                                                                | Individual Calf Pads. <i>6"W x 6"H.</i> |     |          | Standard |
|                                                                                                              | Calf Panel. <i>18"W x 8"H.</i>          |     |          | Standard |

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## STEALTH HEAD SUPPORTS

|                                                          |          |
|----------------------------------------------------------|----------|
| UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>  | FOC      |
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | FOC      |

Please visit [stealthproducts.com](https://stealthproducts.com) for alternative head supports.

## EXPANDABLE CONTROLS

|                                                                                              |                                                                    |   |                                                           |                                                                          |
|----------------------------------------------------------------------------------------------|--------------------------------------------------------------------|---|-----------------------------------------------------------|--------------------------------------------------------------------------|
| Q-Logic 3e Controls. <i>5 power option ready.</i>                                            |                                                                    |   |                                                           | Standard                                                                 |
| iAccess OLED 4 Button Switch Module.<br><i>Mounted to opposite armrest unless specified.</i> | FOC                                                                | } | <u>Push Button iAccess</u><br><u>Toggle Lever iAccess</u> | If iAccess is selected choose either Push Button or Toggle lever option. |
|                                                                                              |                                                                    |   |                                                           |                                                                          |
| Control Method.                                                                              | Hand Control. <i>Select position in Hand Control Mounts below.</i> |   |                                                           | Standard                                                                 |
|                                                                                              | Attendant Control. <i>Hand control mounted at rear.</i>            |   |                                                           | FOC                                                                      |

Q-Logic 3e hand controls are supplied with built-in Bluetooth.

## HAND CONTROL MOUNTS

|                    |            |
|--------------------|------------|
| Fixed.             | Standard   |
| Basic swing-away.  | FOC        |
| Mounting Position. | Left Right |

## SPEED PROGRAM & LIGHTS PACKAGE

|                                                           |          |
|-----------------------------------------------------------|----------|
| 2-Pole Motors & 4 Mph Program with Lights and indicators. | Standard |
| 2-Pole Motors & 6 Mph Program with Lights and indicators. | FOC      |

## ACCESSORIES

|                                                                           |     |
|---------------------------------------------------------------------------|-----|
| Prepare Chassis for Dahl Docking System. <i>Dahl system not included.</i> | FOC |
|---------------------------------------------------------------------------|-----|

## NOTES & SPECIAL REQUESTS



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**MOTABILITY ID 706911 - 01/04/25 V1**

