



MOTABILITY ID 1001591

Q6 EDGE HD WITH TB3 SEATING / POWER TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Clients Weight _____

BASE MODEL

- 10 Colour Options
- Height Adj. Removable Arm Supports
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 32 st 2lbs. (204 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature


DRIVE WHEELS

3" x 14" Flat Free Charcoal Tyres.	Standard
3" x 14" Pneumatic Charcoal Tyres.	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Night Train	FOC	Evolution Orange	FOC
Electric Blue	FOC	Grey Street	FOC
Midnight Blue	FOC	Green Machine	FOC
Deep Purple	FOC		
Red Corvette	FOC		

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

20"-24"W _____

24"-28"W _____

28"-32"W _____

Set width at: _____

SEAT PACKAGES - SEAT DEPTH

20"-25"D _____

Set depth at: _____

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	STD
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

HD Powered Seat Tilt. <i>0° to 45°. Uses 1 power option. Only available with R-Net controls.</i>	Standard
HD Tilt Seat to Floor Height . <i>From seat pan.</i>	18.5" 19.5" 20.5"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
--	----------

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	Standard
---	----------

BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard	
Back Cane Only Support Positioning.	Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
Back Cane Height.	16"H. 17"H. 18"H. 19"H. 20"H.	
<i>From seat pan.</i>	21"H. 22"H. 23"H. 24"H. 25"H.	

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.	16"-20"W.	20"-24"W.
	14"-28"W.	28"-32"W.
Back Support Width.	18"W. 19"W. 20"W. 21"W.	
	22"W. 24"W. 26"W. 28"W.	
	30"W. 32"W.	
Back Support Height.	16"H. 17"H. 18"H. 19"H.	
	20"H. 21"H. 22"H. 23"H.	
<i>From seat pan.</i>	24"H. 25"H. 26"H. 27"H.	
	28"H.	

ARM SUPPORTS

Height Adjustable, Removable, Heavy Duty.				Standard
HD Removable Height Package. <i>From seat pan.</i>	8"-12".	FOC	10"-14".	Standard
TB3 Flip-back with Down Posts. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				FOC
Arm Support Pads.	10" Standard			Standard
	14" Standard			FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>			FOC
	10" Waterfall.			FOC
	14" Waterfall.			FOC

CENTRE MOUNTED LEG SUPPORTS

Adjustable Centre Mount Foot Platform. <i>Mounted to seat.</i>				Standard
Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	Standard	13"-19".	FOC
Foot Supports.	10"W x 8"D Split Foot Supports.	FOC	10"W x 11"D Split Foot Supports.	Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	FOC	Calf Panel. <i>18"W x 8"H.</i>	FOC
	Individual Large Calf Pads. <i>8"W x 8"H tapering to 6"W.</i>	Standard		

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

70° Heavy Duty, Drop-in, Swing-away Leg Supports. <i>Select lowers and footplates below.</i>						FOC	
70° HD	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12"	FOC	13"-18"	Standard	16"-21"	FOC
	Foot Support Type.	Angle & Depth Adjustable. <i>Select size below.</i>				Standard	
	Foot Support Size.	6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
Heavy Duty, Manual Elevating, Swing-away Leg Supports.						FOC	
ELEVATING HD	Lower Leg Extensions. <i>From seat pan.</i>	11"-14"	FOC	14"-17"	Standard	12.5"-16.5"	FOC
		15.5"-19.5"	FOC				
	Foot Support Type.	Flip-up Extruded. <i>8.5"W x 6"D.</i>				Standard	
		Angle & Depth Adjustable. <i>Select size below.</i>				FOC	
	Foot Support Size. <i>Angle adj. only.</i>	6"W x 8"D				Standard	
	6"W x 11"D				FOC		
	8"W x 11"D				FOC		

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	FOC

STEALTH PELVIC & THIGH SOLUTIONS

Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.	Standard
---------------	--	----------

EXPANDABLE CONTROLS

Q-Logic 3e Controls.	<i>3 power option ready.</i>	Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

HAND CONTROL MOUNTS

Fixed.					Standard
Basic Swing-away.					FOC
Mounting Position.	Left	FOC	Right	Standard	

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	Standard
---	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 1001591 - 01/04/25 V1

