

## R-TRAK WITH TB FLEX SEATING / POWER TILT / PRO BACK

### DEALERS DETAILS

|              |           |
|--------------|-----------|
| Dealers Name |           |
| Order Date   |           |
| Address      |           |
| Postcode     |           |
| Phone        | Email     |
| PO No        | Order Ref |

### DELIVERY DETAILS

|                    |  |
|--------------------|--|
| Dealers Address    |  |
| Postcode           |  |
| Quantum Specialist |  |
| Ⓢ Clients Weight   |  |

### BASE MODEL

14 Colour Options  
Height Adj. Removable Arm Supports  
Centre Mount Foot Platform  
Angle Adj. Fold-Forward Back Support  
14" Flat Free Drive Tyres  
3" x 9" Flat Free Front Castors



SRS Suspension System  
QLNES Control System  
2 x 70 ah Batteries  
4 Mph (6 Km/h) Motors  
Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

|                                                                            |          |
|----------------------------------------------------------------------------|----------|
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | Standard |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC      |

### FRONT CASTORS

|                                |          |
|--------------------------------|----------|
| 3" x 9" Flat Free Black Tyres. | Standard |
| 3" x 9" Pneumatic Black Tyres. | FOC      |

### COLOURS

|               |          |                                              |     |
|---------------|----------|----------------------------------------------|-----|
| Back in Black | Standard | Pretty in Pink                               | FOC |
| Green Machine | FOC      | Lemon Crush. <i>Extended lead-time.</i>      | FOC |
| Midnight Blue | FOC      | Night Train. <i>Extended lead-time.</i>      | FOC |
| Electric Blue | FOC      | Evolution Orange. <i>Extended lead-time.</i> | FOC |
| Blue Jean     | FOC      | Military Green. <i>Extended lead-time.</i>   | FOC |
| Deep Purple   | FOC      | Grey Street. <i>Extended lead-time.</i>      | FOC |
| Red Corvette  | FOC      | Raspberry Beret. <i>Extended lead-time.</i>  | FOC |

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

|       |     |       |     |       |     |
|-------|-----|-------|-----|-------|-----|
| 15"W. | FOC | 18"W. | STD | 21"W. | FOC |
| 16"W. | FOC | 19"W. | FOC | 22"W. | FOC |
| 17"W. | FOC | 20"W. | FOC |       |     |

## SEAT PACKAGES - SEAT DEPTH

|       |     |       |          |
|-------|-----|-------|----------|
| 15"D. | FOC | 18"D. | Standard |
| 16"D. | FOC | 19"D. | FOC      |
| 17"D. | FOC | 20"D. | FOC      |

## STEALTH SEAT CUSHIONS

|                                                                                                              |          |
|--------------------------------------------------------------------------------------------------------------|----------|
| No Seat Cushion Required                                                                                     | FOC      |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>                          | Standard |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>                  | FOC      |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC      |

Please note that seat cushions may not be available in all width and depth combinations.  
Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

| Powered Seat Tilt. 0° to 30°. Uses 1 power option. |       |     |                                  | Standard |                                   |
|----------------------------------------------------|-------|-----|----------------------------------|----------|-----------------------------------|
| Powered Tilt Seat to Floor Height. From seat pan.  | 17.5" | 18" | Changes Tilt range to -5° to 25° | 18.5"    | Changes Tilt range to -10° to 20° |

## STANDARD BACK SUPPORT POSTIONING

|                                                                     |          |
|---------------------------------------------------------------------|----------|
| Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°.</i> | Standard |
|---------------------------------------------------------------------|----------|

## STANDARD BACK SUPPORT

|                                                                                           |          |
|-------------------------------------------------------------------------------------------|----------|
| TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i> | Standard |
|-------------------------------------------------------------------------------------------|----------|

## STEALTH PRO BACK SUPPORTS

|                                                                               |                             |          |                             |          |
|-------------------------------------------------------------------------------|-----------------------------|----------|-----------------------------|----------|
| Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>                    |                             |          |                             | FOC      |
| Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>                     |                             |          |                             | Standard |
| Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>          |                             |          |                             | FOC      |
| Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i> |                             |          |                             | FOC      |
| Back Height                                                                   | 50cm. <i>Actual height.</i> | Standard | 57cm. <i>Actual height.</i> | FOC      |

## ARM SUPPORTS

|                                                                                                    |              |          |                                                   |          |
|----------------------------------------------------------------------------------------------------|--------------|----------|---------------------------------------------------|----------|
| TB Flex Removable, Quick Release. <i>Adj. range 9" to 14" from seat pan.</i>                       |              |          |                                                   | Standard |
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i> |              |          |                                                   | FOC      |
| Arm Support Pads.                                                                                  | 10" Standard | Standard | 10" Waterfall. <i>Not available with TB Flex.</i> | FOC      |
|                                                                                                    | 14" Standard | FOC      | 14" Waterfall. <i>Not available with TB Flex.</i> | FOC      |

## CENTRE MOUNTED LEG SUPPORTS

|                                                 |                                         |     |          |          |
|-------------------------------------------------|-----------------------------------------|-----|----------|----------|
| Manually Adjustable Centre Mount Foot Platform. |                                         |     |          | Standard |
| Lower Leg Extensions.<br><i>From seat pan.</i>  | 10"-16".                                | FOC | 13"-19". | Standard |
| Foot Supports.                                  | 10"W x 8"D Split Foot Supports.         |     |          | FOC      |
|                                                 | 10"W x 11"D Split Foot Supports.        |     |          | Standard |
| Calf Support.                                   | Individual Calf Pads. <i>6"W x 6"H.</i> |     |          | Standard |

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED LEG SUPPORTS

|                                                                                                                             |          |
|-----------------------------------------------------------------------------------------------------------------------------|----------|
| 70° Swing-Away. <i>Adjustment range 12" to 15".</i>                                                                         | Standard |
| Manual Elevating Swing-Away. <i>Adjustment range 15" to 20". Minimum seat width 16".</i>                                    | FOC      |
| Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i> |          |

Adjustment ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops.  
Elevating leg supports supplied with calf supports.

## STEALTH HEAD SUPPORTS

|                                                          |          |
|----------------------------------------------------------|----------|
| UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>  | FOC      |
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | FOC      |

## PELVIC SOLUTIONS

|                                            |          |
|--------------------------------------------|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--------------------------------------------|----------|

## NON-EXPANDABLE CONTROLS

|                                              |                                                                    |          |
|----------------------------------------------|--------------------------------------------------------------------|----------|
| QLNES Controls. <i>2 power option ready.</i> |                                                                    | Standard |
| Control Method.                              | Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard |
|                                              | Attendant Control. <i>Hand control mounted at rear.</i>            | FOC      |

## HAND CONTROL MOUNTS

|                    |      |          |
|--------------------|------|----------|
| Fixed.             |      | Standard |
| Basic swing-away.  |      | FOC      |
| Mounting Position. | Left | Right    |

## SPEED PROGRAM & LIGHTS PACKAGE

|                                                           |          |
|-----------------------------------------------------------|----------|
| 2-Pole Motors & 4 Mph Program with Lights and indicators. | FOC      |
| 2-Pole Motors & 6 Mph Program with Lights and indicators. | Standard |

## NOTES & SPECIAL REQUESTS



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**Motability ID 1001071 - 01/04/25 V1**

