

R-TRAK WITH TB3 SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 14 Colour Options
- Height Adj. Removable Arm Supports
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Front Castors



- SRS Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-15"W.
16"-19"W.
20"-22"W.
Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-15"D.
16"-19"D.
20"-22"D.
Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.	Standard
Static Seat to Floor Height. <i>From seat pan.</i>	15" 16" 17"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
--	----------

TB3 BACK SUPPORT POSITIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports. Not suitable for occupied transport.</i>	FOC

BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i>	Standard
Back Cane Only Support Positioning. Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	FOC
Back Cane Height. <i>From seat pan.</i>	16"H. 17"H. 20"H. 21"H. 24"H. 25"H.

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.
		18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.		
Back Support Height. <i>From seat pan.</i>	14"H. <i>12"W back only.</i>	15"H. <i>12"W back only.</i>	16"H.	17"H.
	18"H.	19"H.	20"H.	21"H.
	22"H.	23"H.	24"H.	25"H.
	26"H.	27"H.	28"H.	29"H.

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i>				Standard	
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>				FOC	
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>				FOC	
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				FOC	
Height Adjustable, Removable, Heavy Duty.				FOC	
HD Removable Height Package. <i>From seat pan.</i>		8"-12".	FOC	10"-14".	Standard
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment. Not available with back canes.</i>				FOC	
Arm Support Pads.	10" Standard			Standard	
	14" Standard			FOC	
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>			FOC	
	10" Waterfall.			FOC	
	14" Waterfall.			FOC	

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>				Standard		
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>				FOC		
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". <i>Manually adjustable only.</i>	FOC	13"-19". <i>Manually adjustable only.</i>	Standard
Foot Supports.	10"W x 6"D Split Foot Supports. <i>Paediatric centre mount only.</i>				FOC	
	10"W x 8"D Split Foot Supports.				FOC	
	10"W x 11"D Split Foot Supports.				Standard	
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>	FOC		

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>				FOC
Manual Elevating Swing-Away. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>				FOC

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-Away Leg Supports. *Angle & depth adjustable footplates 6"W x 8"D. Select lowers below. Minimum seat width 15"* FOC

60° TO 90°

Lower Leg Extensions. *From seat pan.* 11"-13" FOC 13.5"-16.5" Standard 15.5"-19.5" FOC

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. *HMO475-17M / CP08CCFPM.* FOC

UniLink 10" Comfort Plus. *HMO475-17M / CP10CCFPM.* Standard

UniLink 14" Comfort Plus. *HMO475-17M / CP14CCFPM.* FOC

PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release. Standard

EXPANDABLE CONTROLS

Q-Logic 3e hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls. *No power options.* Standard

Control Method. Hand Control. *Select position in Hand Control Mounts below.* Standard

Attendant Control. *Hand control mounted at rear.* FOC

HAND CONTROL MOUNTS

Fixed. Standard

Basic Swing-away. FOC

Mounting Position. Left Right

SPEED PROGRAM & LIGHTS PACKAGE

2-Pole Motors & 4 Mph Program with Lights and indicators. Standard

2-Pole Motors & 6 Mph Program with Lights and indicators. FOC

LED Front Castor and Rear Fender Lights. *Includes single switch button mounted to side rail. Not road legal.* Standard

Omit LED Front Castor and Rear Fender Lights. FOC

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITYID 1000793 - 01/04/25 V1

