

Technical Bulletin

Bulletin Ref.:

UK-Q2-3-25

Date Issued:

23 May 2025

BULLETIN DETAILS

Product Group: Footplates

Model(s) Affected: Powerchairs

Effective Date:

1 May 2025

Description of Changes/ Diviations:

Quantum Rehab has begun rolling out a new style of footplates across our powerchair range.

Key Improvements

1. Weight Reduction

- The new footplates are 0.6 kg lighter per footplate, reducing overall chair weight.

2. Enhanced Durability

- The previous stick-on coating, which could peel over time, has been replaced with a durable etched coating, more resilient and permanently bonded.



OLD STYLE FOOTPLATES



NEW STYLE FOOTPLATES



The transition is already in effect. As existing stocks of the old-style footplates run out by size, orders and production are automatically fulfilled with the new version.

Documentation Affected: ☐ Owner's Manual ☒ Order Form(s) ☐ Spec Sheets ☒ Other ☐ N/A

Link to new version of affected document:

For more detailed information, please contact your local designated Area Sales Manager.

NB: This Bulletin is for informational purposes only and DOES NOT replace the Specifications.

Pride Mobility Products Ltd.
32 Wedgwood Road, Bicester,
Oxfordshire, OX26 4UL
E: sales@quantumrehab.co.uk
T: +44 (0)1869 324600