



Price Guide

1st April – 30th June 2025

Content

How does the scheme work?	4
Edge 3 Stretto with Contoured Seat	6
Edge 3 Stretto with TB Flex Seat	8
Edge 3 Stretto with TRU-Balance 3 Seat	10
Edge 3 with TRU-Balance 3 Seat	12
4Front 2 with TRU-Balance 3 Seat	14
R-Trak with TB Flex Seat	16
R-Trak with TRU-Balance 3 Seat	18

Outback with Contoured Seat	20
Q6 Edge 2.0 with TB Flex Seat	22
Q6 Edge 2.0 with TRU-Balance 3 Seat	24
Q6 Edge 2.0 HD with TRU-Balance 3 Seat	26
Q6 Edge HD with Synergy Seat	28
Q6 Edge HD with TRU-Balance 3 Seat	30
4Front 2 HD with TRU-Balance 3 Seat	32
Q1450 with Contoured Seat	34
Q1450 with Synergy Seat	36
Aspen with TB Flex Seat	38
Kozmo Portable	40
Fusion e	42



Featuring the R-Trak with
Tru-balance 3 seating



How Does the Motability Scheme work?



Use your government-funded mobility allowance to lease a powered wheelchair.



A brand-new product + Full support package:



How to Get Started

How it works



Exchange all or part of your allowance

Lease duration:



3 years (standard)
5 years (custom-made)

What's Included



Insurance



Servicing & Maintenance



Full Motability Assist
Breakdown Cover



Yearly Tax
(for Class 3 products)



Tyre & Battery Replacements

Getting Started

1

Check Eligibility
Use the Motability Scheme Tool online.

2

Arrange a Demonstration
At home or at a local showroom.

3

Choose & Confirm
Select your powerchair and consult a dealer.

4

Place Your Order

5

Home Delivery
Enjoy hassle-free delivery straight to your door!

EDGE 3 Stretto

CONTOURED
SEATING


Weight capacity
136 kg


Battery range
Up to 15 miles


Crash tested




Fold forward
backrest



iLevel Stability
Technology



Narrow base
20.47"
(520 mm)



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The Edge 3 Stretto is Quantum's narrowest power chair, perfect for navigating tight spaces. Paired with the Contoured seat, it offers a supportive design with a vinyl, easy-to-clean surface that promotes optimal posture and pressure relief for all-day comfort.

Specifications

*All measurements provided are subject to a variation of +/- 1/8"

Standard

Turning Radius	21" (533 mm)
Length	36" (917 mm)
Width	20.47" (520 mm) with 12.5" wheels 21.75" (552 mm) with 14" wheels
Ground Clearance	2.9" (74 mm)
Battery Size	2 x 62 Ah
Seat Dimensions	16" W x 16" D / 18" W x 18" D / 20" W x 20" D
Motor Package	2 - Pole, Mid - Wheel
Suspension	Smooth Ride Suspension

Optional

Powered Positioning	12" iLevel Powered Seat Lift
Speeds	4 mph (6.44 kmh) 6 mph (9.65 kmh)
Electronics	Q-Logic NES / VR2

Motability Options

ID Reference

PIP weekly costs

Edge 3 Stretto - Contoured seat / standard	ID 706811	£ 34.48
Edge 3 Stretto - Contoured seat / 6 mph	ID 706910	£ 40.87
Edge 3 Stretto - Contoured seat	ID 706914	£ 38.93
Edge 3 Stretto - Contoured seat / iLevel	ID 1000398	£ 50.43



Weight capacity
136 kg

Battery range
Up to 15 miles

Crash tested




SRS
Smooth Ride
Suspension



Safe Seat
Elevation



Narrow base at
20.47"
(520 mm)



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The Edge 3 Stretto is Quantum's narrowest power chair, perfect for navigating tight spaces. Paired with the TB Flex seat, it offers flexibility and comfort. Ideal for daily activities and exercise, it adapts to your movements, enhancing performance and reducing strain.

Specifications

*All measurements provided are subject to a variation of $\pm 1/8"$

Standard	Turning Radius	21" (533 mm)
	Length	36" (917 mm)
	Width	20.47" (520 mm) with 12.5" wheels 21.75" (552 mm) with 14" wheels
	Ground Clearance	2.9" (74 mm)
	Battery Size	2 x 40 Ah or 2 x 62 Ah ^{^^}
	Seat Dimensions	15" - 20" W / 15" - 20" D
	Motor Package	2 - Pole, Mid - Wheel
	Suspension	Smooth Ride Suspension
	Speeds	4 mph (6.44 kmh)
Optional		
Powered Positioning		Tilt / Recline / Lift
Electronics		Q-Logic NES / Q-Logic 3e

Motability Options

ID Reference

PIP weekly costs

Edge 3 Stretto with TB Flex seat ^{^^}	ID 707226	£ 48.54
Edge 3 Stretto with TB Flex Seat/ Pro Back ^{^^}	ID 1000838	£ 50.06
Edge 3 Stretto with TB Flex /Power Tilt	ID 707225	£ 50.46
Edge 3 Stretto with TB Flex /Power Tilt / Pro Back	ID 1000815	£ 52.69
Edge 3 Stretto with TB Flex /Power Tilt/ Power Recline	ID 707223	£ 52.23
Edge 3 Stretto with TB Flex/ Power tilt/Recline	ID 706572	£ 59.26
Edge 3 Stretto with TB Flex/ Power Lift /Power Tilt ^{^^}	ID 706573	£ 56.52

EDGE 3 Stretto

TRU-BALANCE³
POWER POSITIONING SYSTEMS


Weight capacity
136 kg


Battery range
Up to 15 miles


Crash tested




LED Fender
Lights



iLevel Stability
Technology



Narrow base at
20.47"
(520 mm)



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The Edge 3 Stretto is Quantum's narrowest power chair, perfect for navigating tight spaces. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications

*All measurements provided are subject to a variation of +/- 1/8"

Standard	Turning Radius	21" (533 mm)
	Length	36" (917 mm)
	Width	20.47" (520 mm) with 12.5" wheels 21.75" (552 mm) with 14" wheels
	Ground Clearance	2.9" (74 mm)
	Battery Size	2 x 62 Ah
	Seat Dimensions	12" - 20" W / 12" - 20" D
	Motor Package	2 - Pole, Mid - Wheel
	Suspension	Smooth Ride Suspension
Optional		
Powered Positioning		Tilt / Recline / iLevel
Speeds		4 mph (6.44 kmh) 6 mph (9.65 kmh)
Electronics		QLNES / Q-Logic 3e / Q-Logic 3 EX

Motability Options

ID Reference

PIP weekly costs

Edge 3 Stretto with TB3 Seat	ID 707384	£ 52.65
Edge 3 Stretto with TB3 Seat/Powered Tilt	ID 706913	£ 58.69
Edge 3 Stretto with TB3 Seat/Powered Tilt/ Powered Recline	ID 707383	£ 59.83
Edge 3 Stretto with TB3 Seat/iLevel/Powered Tilt Powered Recline	ID 706909	£ 62.50
Edge 3 Stretto with TB3 Seat/iLevel/Powered Tilt/Powered Recline	ID 707382	£ 62.55

EDGE 3

TRU-BALANCE³
POWER POSITIONING SYSTEMS


Weight capacity
136 kg


Battery range
Up to 15 miles


Crash tested




**SRS
Suspension**



**iLevel Stability
Technology**



**Colour Wheel
Accents**



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.




The Edge 3 is a versatile power chair delivering advanced performance, comfort, and customisation for diverse user needs. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications		*All measurements provided are subject to a variation of +/- 1/8"
Standard	Turning Radius	20.5" (521 mm)
	Length	32.23" (819 mm)
	Width	24.4" (620 mm) with 14" wheels
	Ground Clearance	2.7" (69 mm)
	Battery Size	2 x 62 Ah
	Seat Dimensions	12" - 22" W / 12" - 22" D
	Motor Package	2 - Pole, Mid - Wheel
	Suspension	Smooth Ride Suspension
	Electronics	Q-Logic 3e
Optional		
Powered Positioning		Tilt / Recline / iLevel
Speeds		4 mph (6.44 kmh) 6 mph (9.65 kmh)

Motability Options	ID Reference	PIP weekly costs
Edge 3 with TB3 Seat	ID 710381	£ 54.57
Edge 3 with TB3 Seat/Power Tilt	ID 710380	£ 57.52
Edge 3 with TB3 Seat/Power Tilt/Power Recline	ID 710379	£ 60.58
Edge 3 with TB3 Seat/iLevel/Power Tilt	ID 1001593	£ 63.14

4FRONT²

TRU-BALANCE³
POWER POSITIONING SYSTEMS



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.

The 4Front 2 is a front-wheel-drive power chair built for outdoor performance, with advanced features and all-day comfort. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications		*All measurements provided are subject to a variation of + 1/8"
Standard	Turning Radius	25.9" (658 mm)
	Length	39.9" (1013 mm)
	Width	24" (610 mm) with 14" wheels
	Ground Clearance	3.15" (80 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	12" - 22" W / 12" - 22" D
	Motor Package	4 - Pole, Front - wheel
	Suspension	Smooth Ride Suspension
	Electronics	Q-Logic 3e
Optional		
Powered Positioning		Tilt / Recline / Lift
Speeds		4 mph (6.44 kmh) 6 mph (9.65 kmh)

Motability Options	ID Reference	PIP weekly costs
4Front 2 with TB3 Seat	ID 714946	£ 59.01
4 Front 2 with TB3 Seat/Power Tilt	ID 706577	£ 62.25
4Front 2 with TB3 Seat/Power Tilt/ Power Recline	ID 715986	£ 62.57
4Front 2 with TB3 Seat/Power Lift/Power Tilt	ID 706352	£ 62.58
4Front 2 with TB3 Seat/Power Lift/Power Tilt/ Power Recline	ID 1000835	£ 62.70

QUANTUM® R-TRAK

TB-flex

Weight capacity
136 kg

Battery range
Up to 20 miles

Crash tested



Front/Rear LED-
Fender Lights



Safe Seat
Elevation



Colour Wheel
Accents



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The R-Trak combines the stability of a traditional rear-wheel drive power chair with comfort and style. Paired with the TB Flex seat, it offers flexibility and comfort. Ideal for daily activities and exercise, it adapts to your movements, enhancing performance and reducing strain.

Specifications

*All measurements provided are subject to a variation of +/- 1/8"

Standard	Turning Radius	25.7" (653 mm)
	Length	34.8" (884 mm)
	Width	24.2" (615 mm) with 14" wheels
	Ground Clearance	3.97" (101 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	15" - 20" W / 15" - 20" D
	Motor Package	2 - Pole, Rear - wheel
	Suspension	Smooth Ride Suspension
	Electronics	Q-Logic 3e
Optional		
Powered Positioning		Tilt / Recline / Lift
Speeds		4 mph (6.44 kmh) 6 mph (9.65 kmh)

Motability Options

ID Reference

PIP weekly costs

R-Trak with TB-Flex Seat / Power Tilt	ID 1001054	£ 52.06
R-Trak with TB-Flex Seat / Power Tilt / Power Recline	ID 1001592	£ 53.72
R-Trak with TB-Flex Seat / Power Tilt / Pro Back	ID 1001071	£ 54.14

QUANTUM®
R-TRAK

TRU-BALANCE³
POWER POSITIONING SYSTEMS


Weight capacity
136 kg


Battery range
Up to 20 miles


Crash tested



Front/Rear LED-
Fender Lights



Safe Seat
Elevation



Colour Wheel
Accents



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The R-Trak combines the stability of a traditional rear-wheel drive power chair with comfort and style. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications

*All measurements provided are subject to a variation of $\pm 1/8"$

Standard	Turning Radius	25.7" (653 mm)
	Length	34.8" (884 mm)
	Width	24.2" (615 mm) with 14" wheels
	Ground Clearance	3.97" (101 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	12" - 22" W / 12" - 22" D
	Motor Package	2 - Pole, Rear - wheel
	Suspension	Smooth Ride Suspension
	Electronics	Q-Logic 3e
Optional		
	Powered Positioning	Tilt / Recline / Lift
	Speeds	4 mph (6.44 kmh) 6 mph (9.65 kmh)

Motability Options

ID Reference

PIP weekly costs

R-Trak with TB3 Seat	ID 1000793	£ 55.11
R-Trak with TB3 Seat/Power Tilt	ID 1000800	£ 58.23
R-Trak with TB3 Seat/Power Tilt/ Power Recline	ID 1000814	£ 60.74
R-Trak with TB3 Seat/Power Lift/ Power Recline	ID 1000834	£ 62.70



**CONTOURED
SEATING**


Weight capacity
159 kg


Battery range
Up to 16 miles


Non-Crash
tested



Fold forward
backrest



All Terrain



15" Pneumatic
Drive Tyres



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



Explore the great outdoors like never before with the Quantum® Outback, a powerful 4x4 all-terrain power chair built for adventurers and nature lovers. Paired with the Contoured seat, it offers a supportive design with a vinyl, easy-to-clean surface that promotes optimal posture and pressure relief for all-day comfort.

Specifications		*All measurements provided are subject to a variation of +/- 1/8"
Standard	Turning Radius	49.2" (1249.7 mm)
	Length	42" (1067 mm)
	Width	28.82" (732.03 mm)
	Ground Clearance	4.92" (125 mm)
	Battery Size	2 x 95 Ah
	Seat Dimensions	16" - 22" W / 16" - 22" D
	Motor Package	600 W 4 motor, independent
	Suspension	Smooth Ride Suspension
	Powered Positioning	N/A
	Speeds	6 mph (9.65 kmh)
	Electronics	Q-Logic 3e

Motability Options	ID Reference	PIP weekly costs
Outback with Contoured Seat	ID 1002475	£ 63.20

Q6 edge[®]
2.0

TB-flex


Weight capacity
136 kg


Battery range
Up to 15 miles



Crash tested



ATX Suspension



QLNES Controls
with Magnetic
Charger Port



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The revolutionary Q6 Edge 2.0 Power Chair has been engineered to meet the performance needs and demands of the most active user. Paired with the TB Flex seat, it offers flexibility and comfort. Ideal for daily activities and exercise, it adapts to your movements, enhancing performance and reducing strain.

Specifications

*All measurements provided are subject to a variation of +/- 1/8"

Standard	Turning Radius	20.5" (521 mm)
	Length	35.5" (902 mm)
	Width	24.4" (620 mm) with 14" wheels
	Ground Clearance	2.7" (69 mm)
	Battery Size	2 x 62 Ah
	Seat Dimensions	15" - 20" W / 15" - 20" D
	Motor Package	2 - Pole, Mid - wheel
	Suspension	ATX Suspension
	Speeds	4 mph (6.44 kmh)
	Electronics	Q-Logic NES
Optional		
Powered Positioning		Tilt / Recline / Lift

Motability Options

ID Reference

PIP weekly costs

Q6 Edge 2.0 with TB Flex Seat	ID 1002463	£ 42.43
Q6 Edge 2.0 with TB Flex Seat / Pro Back	ID 1002468	£ 42.47
Q6 Edge 2.0 TB Flex Seat / Power Tilt	ID 1002464	£ 47.50
Q6 Edge 2.0 with TB Flex Seat / Power Tilt / Pro Back	ID 1002469	£ 48.45
Q6 Edge 2.0 with TB Flex Seat / Power Tilt / Power Recline	ID 1002467	£ 49.26
Q6 Edge 2.0 with TB Flex Seat / Power Tilt / Power Recline / Pro Back	ID 1002474	£ 50.21

Q6 edge^{2.0}

TRU BALANCE³
POWER POSITIONING SYSTEMS


Weight capacity
136 kg


Battery range
Up to 15 miles


Crash tested




ATX Suspension



iLevel Stability
Technology



QLNES Controls
with Magnetic
Charger Port



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.

The revolutionary Q6 Edge 2.0 power chair has been engineered to meet the performance needs and demands of the most active user. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications

*All measurements provided are subject to a variation of $\pm 1/8"$

Standard	Turning Radius	20.5" (521 mm)
	Length	35.5" (902 mm)
	Width	24.4" (620 mm) with 14" wheels
	Ground Clearance	2.7" (69 mm)
	Battery Size	2 x 62 Ah
	Seat Dimensions	12" - 22" W / 12" - 22" D
	Motor Package	2 - Pole, Mid - Wheel
	Suspension	ATX Suspension
Optional		
	Powered Positioning	Tilt / Recline / Lift / Power Foot board
	Speeds	4 mph (6.44 kmh) 6 mph (9.65 kmh)
	Electronics	Q-Logic NES / Q-Logic 3

Motability Options

ID Reference

PIP weekly costs

Q6 Edge 2.0 with TB3 Seat	ID 721038	£ 48.58
Q6 Edge 2.0 with TB3 Seat/Power Tilt	ID 721032	£ 49.78
Q6 Edge 2.0 with TB3 Seat/Power Tilt /ADI Back	ID 714914	£ 55.12
Q6 Edge 2.0 with TB3 Seat/Power Tilt / Power Recline	ID 721029	£ 56.53
Q6 Edge 2.0 with TB3 Seat/iLevel/Tilt	ID 721026	£ 58.38
Q6 Edge 2.0 with TB3 Seat/iLevel/Power Tilt / Power Recline	ID 718709	£ 60.70
Q6 Edge 2.0 with TB3 Seat/iLevel/Power Tilt/ Power Recline/Power Footboard	ID 706911	£ 62.50

Q6 edge[®] 2.0HD

TRU-BALANCE³
POWER POSITIONING SYSTEMS


Weight capacity
181 kg


Battery range
Up to 15 miles


Crash tested




**ATX
Suspension**



**iLevel Stability
Technology**



4-Pole Motors



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.

The revolutionary Q6 Edge 2.0 HD power chair has been engineered to meet the needs and demands of the most active user with added heavy duty capabilities. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications		*All measurements provided are subject to a variation of $\pm 1/8"$
Standard	Turning Radius	20.5" (521 mm)
	Length	35.5" (902 mm)
	Width	24.4" (620 mm) with 14" wheels
	Ground Clearance	2.7" (69 mm)
	Battery Size	2 x 62 Ah
	Seat Dimensions	20" - 24" W / 20" - 24" D
	Motor Package	4 - Pole, Mid - wheel
	Suspension	ATX Suspension
	Speeds	4 mph (6.44 kmh)
	Electronics	Q-Logic 3e
Optional		
Powered Positioning		Tilt / iLevel

Motability Options	ID Reference	PIP weekly costs
Q6 Edge 2.0 HD with TB3 Seat	ID 706571	£ 51.56
Q6 Edge 2.0 HD with TB3 Seat / Power Tilt	ID 706570	£ 57.91
Q6 Edge 2.0 HD with TB3 Seat / iLevel / Power Tilt	ID 1001594	£ 63.13

Q6 edge[®] HD SYNERGY[®]


Weight capacity
226 kg


Battery range
Up to 20 miles



Crash tested




Q-Logic 3e
Control System


ATX
Suspension


4-Pole Motors



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.

The Q6 Edge HD is a heavy-duty power chair designed for users who need a higher weight capacity without compromising on performance or comfort. Paired with Synergy seat, the Q6 Edge HD offers adjustable seat sizes and positioning for a comfortable, customisable fit. It supports various cushion and backrest combinations for optimal pressure relief and postural support

Specifications		*All measurements provided are subject to a variation of $\pm 1/8"$
Standard	Turning Radius	19.75" (502 mm)
	Length	35.6" (904 mm)
	Width	26.25" (667 mm) with 14" wheels
	Ground Clearance	3" (76 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	20" - 28" W / 15" - 28" D
	Motor Package	4 - Pole, Mid - wheel
	Suspension	ATX Suspension
	Speeds	4 mph (6.44 kmh)
Optional		
Powered Positioning		Tilt
Electronics		Q-Logic 3e / R-Net

Motability Options	ID Reference	PIP weekly costs
Q6 Edge HD with Synergy Seat	ID 721035	£ 52.73
Q6 Edge HD with Synergy Seat / Power Tilt	ID 721025	£ 60.62

Q6 edge[®] HD

TRU-BALANCE³
POWER POSITIONING SYSTEMS


Weight capacity
204 kg


Battery range
Up to 20 miles


Crash tested




Q-Logic 3e
Control System



ATX Suspension



4-Pole Motors



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.

The Q6 Edge HD is a heavy-duty power chair designed for users who need a higher weight capacity without compromising on performance or comfort. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications

*All measurements provided are subject to a variation of $\pm 1/8"$

Standard	Turning Radius	20.5" (521 mm)
	Length	35.5" (902 mm)
	Width	24.4" (620 mm) with 14" wheels
	Ground Clearance	2.7" (69 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	20" - 32" W / 20" - 25" D
	Motor Package	4 - Pole, Mid - wheel
	Suspension	ATX Suspension
	Speeds	4 mph (6.44 kmh)
	Electronics	Q-Logic 3e
Optional		
Powered Positioning		Tilt / Recline

Motability Options

	ID Reference	PIP weekly costs
Q6 Edge HD with TB3 Seat	ID 714945	£ 57.56
Q6 Edge HD with TB3 Seat / Power Tilt	ID 714944	£ 62.64
Q6 Edge HD with TB3 Seat / Power Tilt / Power Recline	ID 1001591	£ 63.12

QUANTUM®
4FRONT² HD

TRU-BALANCE³
POWER POSITIONING SYSTEMS


Weight capacity
204 kg


Battery range
Up to 20 miles


Crash tested




Front/ Rear LED
Fender lights



Smart Traction
Control



SRS Smooth Ride
Suspension



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The 4Front 2 HD has all the benefits of the 4 Front 2 but with Heavy-Duty capabilities. Paired with the TRU-Balance 3 seat, it offers a customisable seat and positioning system that ensures optimal alignment and pressure distribution. It offers maximum comfort and support for various postures and activities.

Specifications

*All measurements provided are subject to a variation of $\pm 1/8"$


Standard	Turning Radius	24.5" (622 mm)
	Length	39.7" (1008 mm)
	Width	24" (610 mm) with 14" wheels
	Ground Clearance	2.6" (66 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	20" - 24" W / 20" - 24" D
	Motor Package	4 - Pole, Front - wheel
	Suspension	Smooth Ride Suspension
	Speeds	4 mph (6.44 kmh)
	Electronics	Q-Logic 3e
Optional		
Powered Positioning		Tilt / Recline

Motability Options

	ID Reference	PIP weekly costs
4Front 2 HD with TB3 Seat	ID 706576	£ 60.28
4Front 2 HD with TB3 Seat / Power Tilt	ID 706575	£ 61.37
4Front 2 HD with TB3 Seat / Power Tilt / Power Recline	ID 706574	£ 62.50

QUANTUM 1450

CONTOURED
SEATING


Weight capacity
272 kg


Battery range
Up to 15 miles


Non-Crash
tested




High Traction
Drive Tyres



Oversized Rear
Castor Wheels



Enhanced Weight
Capacity




*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.

The Q1450 power chair front-wheel drive bariatric power chair has an innovative design that gives remarkable outdoor performance while allowing for confined space manoeuvrability. Paired with the Contoured seat, it offers a supportive design with a vinyl, easy-to-clean surface that promotes optimal posture and pressure relief for all-day comfort.

Specifications		*All measurements provided are subject to a variation of + 1/8"
Standard	Turning Radius	26.5" (673 mm)
	Length	38.25" (972 mm)
	Width	29.25" (743 mm)
	Ground Clearance	1.5" (38 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	22" - 26" W / 20" - 24" D
	Motor Package	4 - Pole, Front - wheel
	Suspension	Swing - Arm Suspension
	Powered Positioning	N/A
	Speeds	4 mph (6.44 kmh)
	Electronics	VR2

Motability Options	ID Reference	PIP weekly costs
Q1450 with Contoured Seat	ID 713664	£ 62.60

QUANTUM 1450 SYNERGY


Weight capacity
272 kg


Battery range
Up to 15 miles


Non-Crash
tested




High Traction
Drive Tyres



Oversized Rear
Castor Wheels



Enhanced Weight
Capacity



*Colour shades shown here are for reference only and may vary from actual product.
Please visit our website for the full colour range.



The Q1450 power chair front-wheel drive bariatric power chair has an innovative design that gives remarkable outdoor performance while allowing for confined space manoeuvrability. Paired with Synergy seat, the Q6 Edge HD offers adjustable seat sizes and positioning for a comfortable, customisable fit. It supports various cushion and backrest combinations for optimal pressure relief and postural support.

Specifications		*All measurements provided are subject to a variation of $\pm 1/8"$
Standard	Turning Radius	26.5" (673 mm)
	Length	38.25" (972 mm)
	Width	29.25" (743 mm)
	Ground Clearance	1.5" (38 mm)
	Battery Size	2 x 70 Ah
	Seat Dimensions	20" - 32" W / 20" - 28" D
	Motor Package	4 - Pole, Front - wheel
	Suspension	Swing - Arm Suspension
	Powered Positioning	N/A
	Speeds	4 mph (6.44 kmh)
	Electronics	VR2 Controller

Motability Options	ID Reference	PIP weekly costs
Q1450 with Synergy Seat	ID 713665	£ 62.60

ASPEN

TB-Flex

Weight capacity
136 kg

Battery range
Up to 15 miles


Crash tested



Rear - Wheel
Drive



Kerb Climber



*Colour shades shown here are for reference only and may vary from actual product. Please visit our website for the full colour range.



The Aspen rear-wheel drive power chair boasts a compact frame, making it more manoeuvrable than other rear-wheel drive power chairs. Paired with the TB Flex seat, it offers flexibility and comfort. Ideal for daily activities and exercise, it adapts to your movements, enhancing performance and reducing strain.

Specifications

*All measurements provided are subject to a variation of $\pm 1/8"$

Standard	Turning Radius	25.25" (642 mm)
	Length	33.41" (849 mm)
	Width	22.79" (579 mm)
	Ground Clearance	3.58" (91 mm)
	Battery Size	2 x 62 Ah
	Seat Dimensions	15" - 20" W / 15" - 20" D
	Motor Package	2 - Pole, Rear - wheel
	Suspension	Independent Rear Suspension
Optional		
	Powered Positioning	Tilt / Recline / Lift
	Speeds	4 mph (6.44 kmh) 6 mph (9.65 kmh)
	Electronics	QLNES / Q-Logic 3e

Motability Options

ID Reference

PIP weekly costs


Aspen with TB Flex Seat	ID 721049	£ 42.90
Aspen with TB Flex Seat / Power Tilt	ID 721044	£ 47.06
Aspen with TB Flex Seat / Power Tilt / Power Recline	ID 721040	£ 47.66
Aspen with TB Flex Seat / Power Lift / Power Tilt / Power Recline	ID 1001366	£ 54.13



Disassembles
into 4 pieces



Compact Design



*Colour shades shown here are for reference only
and may vary from actual product.
Please visit our website for the full colour range.

The Kozmo Portable combines manoeuvrability with a compact design and an exceptionally tight turning circle of 646 mm, making it ideal for tight spaces. The Kozmo Portable is packed with fantastic features and adjustable seat, and it's easy to disassemble for transport or storage.

Specifications		*All measurements provided are subject to a variation of $\pm 1/8"$
Standard	Turning Radius	25.45" (646 mm)
	Length	33.09" (840 mm)
	Width	22.6" (574 mm)
	Ground Clearance	1.9" (48 mm)
	Battery Size	2 x 18 Ah
	Seat Dimensions	12" - 18" W / 12" - 18" D
	Motor Package	Rear - wheel
	Suspension	Limited
	Powered Positioning	N/A
	Speeds	3.6 Mph (5.79 Km/h)
	Electronics	GC3

Motability Options	ID Reference	PIP weekly costs
Kozmo Portable	ID 713666	£ 31.90



10" castor wheels
for improved outdoor
performance



Enhanced
Weight Capacity



Back in
Black

*Colour shades shown here are for reference only
and may vary from actual product.
Please visit our website for the full colour range.

The Fusion e is an affordable rear-wheel drive power wheelchair that provides great performance and a smooth riding experience. The adaptable Fusion e power chair is also highly adjustable, allowing users to tailor it to their unique needs and preferences.

Specifications		*All measurements provided are subject to a variation of +/- 1/8"
Standard	Turning Radius	25.25" (642 mm)
	Length	41.33" (1050 mm)
	Width	24.8" (630 mm)
	Ground Clearance	3.14" (80 mm)
	Battery Size	2 x 50 Ah
	Seat Dimensions	16" - 23" W / 16" - 21" D
	Motor Package	2 - Pole, Rear - wheel
	Suspension	Independent Rear Suspension
	Powered Positioning	N/A
	Electronics	VR2
Optional		
Speeds		4 mph (6.44 kmh) 6 mph (9.65 kmh)

Motability Options	ID Reference	PIP weekly costs
Fusion e Pro Back 4 mph	ID 706356	£ 37.12
Fusion e Pro Back 6 mph	ID 706355	£ 37.82



01869 324600

quantumrehab.co.uk

sales@quantumrehab.co.uk



UK
CA | CE
V1 - 15/04/2025