



# Dyno SE UK Retail Order Form



June 2021

Stealth UK Customer Service: 01869324600

Email: sales@quantumrehab.co.uk

Date: \_\_\_\_\_  QUOTE  ORDER - PO# \_\_\_\_\_

Company Name: \_\_\_\_\_

Ship To: \_\_\_\_\_

Contact Name: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Billing Address: \_\_\_\_\_

Post Code: \_\_\_\_\_

Post Code: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Mark for: \_\_\_\_\_



Model shown with accessories.

## DYNO SE™ INCLUDES:

- Frame Suspension
- Hip Angle Adjustment
- QR 16" Rear Wheels
- 10" Front Wheel with Swivel Lock
- Height Adjustable Push Handle
- Attendant Wheel Lock
- Under Seat Storage Bag
- Rear Pocket
- Sunshade

## Specifications

Max Occupant Weight	80 lbs.	36 kg
Push Handle Height	30-43 in	76-109 cm
Weight	26.5 lbs.	12 kg
Overall Width	26.5 in	67 cm
Front Wheel Diameter	10 in	25 cm
Rear Wheel Diameter	16 in	41 cm
Folded without rear wheels	32x17x25 in	82x43x64 cm
Folded with wheels	39x17x26.5 in	99x43x67 cm
Reversing Width	41 in	104 cm
Tilt Positions	25-30°	
Maximum Safe Slope	5°	
Knee Active Range	10°	
Seat-to-back Active Range	30°	

MEASUREMENT WORKSHEET	
A. Seat Depth	_____
B. Trunk Depth	_____
C. Trunk Width	_____
D. Hip Width	_____
E. Seat to Top of Shoulder	_____
F. Lower Leg Length	_____
G. Seat to Axilla	_____

## Adjustability

Seat Depth	8.5-13.5 in	22-34 cm
Seat Width	14 in	36 cm
Lower Leg Length	7.5-14 in	19-36 cm
Back Height	11.5-21 in	29-53 cm
Hip Angle Adjustment	85°, 90°, 105°	
Chest Width Range	7-14 in	18-35 cm

## DYNO™ Chair

### DYNO SE™

<input type="checkbox"/>	3001-1002	E1238	DYNO SE™	£1,499.00
--------------------------	-----------	-------	----------	-----------

### Upholstery Color

<input type="checkbox"/>	3001-0101	BIG SKY BLUE - Blue/Black	N/C
<input type="checkbox"/>	3001-0102	ROSSO RED - Red/Black	N/C
<input type="checkbox"/>	3001-0103	ICE SILVER - Silver/Black	N/C
<input type="checkbox"/>	3001-0104	SNOWBERRY PINK - Pink/Black	N/C

*Upholstery Color applies to back, seat cushions, and sunshade*

## Transit Option & Transit Harness Selection

### Transit Option

*Requirements: Sure-Lok belt, 3001-4001 back rest and a head rest.*

<input type="checkbox"/>	3001-7007	Transit Option	£187.00
--------------------------	-----------	----------------	---------

*This option is not retrofittable.*

*The transit option includes 4 installed transportation anchors and belt mount studs*

### Transit Belt or Harness

<input type="checkbox"/>	3001-7008	K0108	Transit Sure-Lok Lap Belt	£68.00
--------------------------	-----------	-------	---------------------------	--------

## Features & Accessories

<input checked="" type="checkbox"/>	3001-5001	Height Adjustable Foot Plate	N/C
<input checked="" type="checkbox"/>	3001-7009	Attendant Wheel Lock	N/C
<input checked="" type="checkbox"/>	3001-7010	Under Seat Storage	N/C
<input checked="" type="checkbox"/>	3001-7012	Deluxe Sunshade	N/C

<input type="checkbox"/>	3001-7015	DYNO Grab Bar	£101.00
<input type="checkbox"/>	3001-7020	DYNO Rain Cover	£146.00
<input checked="" type="checkbox"/>	3001-2001	10 in (25 cm) Front Wheel with Pneu. Tire	N/C
<input checked="" type="checkbox"/>	3001-2002	16 in (41 cm) Rear Wheels with Pneu. Tires	N/C

## ACTIVATOR™ Dynamic Component

*If ACTIVATOR™ Dynamic Component is selected, select spring to right.*

### A. ACTIVATOR™ Dynamic Component

<input type="checkbox"/>	3001-7018	E2295	ACTIVATOR™ Dynamic Component	£295.00
--------------------------	-----------	-------	------------------------------	---------

### ACTIVATOR™ Springs (Only available with A)

<input type="checkbox"/>	3001-7003		Light Springs for ACTIVATOR™	N/C
<input type="checkbox"/>	3001-7004	Default	Medium Springs for ACTIVATOR™	N/C
<input type="checkbox"/>	3001-7005		Stiff Springs for ACTIVATOR™	N/C
<input type="checkbox"/>	3001-7006		Full ACTIVATOR™ Spring Kit	£38.00

## Back Rest & Trunk Supports

*Select back rest and review related options to the right.*

### A. Contoured Back

<input type="checkbox"/>	3001-4001	E2293	Contoured Configurable Back	N/C
--------------------------	-----------	-------	-----------------------------	-----

### Fixed Lateral Trunk Supports (Only available with A)

<input type="checkbox"/>	3003-6001	E0956	LEFT, Lateral Trunk Support 3x4 in (7.5x10 cm)	N/C
<input type="checkbox"/>	3003-6002	E0956	RIGHT, Lateral Trunk Support 3x4 in (7.5x10 cm)	N/C

### Swing Away Lateral Trunk Supports (Only available with A)

*Must select hardware & pad.*

<input type="checkbox"/>	3002-6111	E1028	LEFT, Leggero Swing In/Out Lateral	£168.00	
<input type="checkbox"/>	3002-6112	E1028	RIGHT, Leggero Swing In/Out Lateral	£168.00	
<input type="checkbox"/>	TWBTL-L-01	E1028	LEFT, Stealth Thin Lateral	£152.00	
<input type="checkbox"/>	TWBTL-R-01	E1028	RIGHT, Stealth Thin Lateral	£152.00	
L	R	3004-6035	E0956	Lateral Pad 3x5", Each	£34.00
L	R	3004-6045	E0956	Lateral Pad 4x5", Each	£34.00
L	R	3004-6046	E0956	Lateral Pad 4x6", Each	£34.00

### Seat Cushion

<input type="checkbox"/>	3001-3001	E2292	DYNO Cascading Contoured Seat Cushion	N/C
--------------------------	-----------	-------	---------------------------------------	-----

### Additional Positioning

#### DYNO Lateral/Pelvic Thigh Supports

<input type="checkbox"/>	3001-6001	E0956	LEFT, Lateral/Pelvic Thigh Support	N/C
Select Pad Thickness: <input type="checkbox"/> 1" <input type="checkbox"/> 2" <input type="checkbox"/> 3"				
<input type="checkbox"/>	3001-6002	E0956	RIGHT, Lateral/Pelvic Thigh Support	N/C
Select Pad Thickness: <input type="checkbox"/> 1" <input type="checkbox"/> 2" <input type="checkbox"/> 3"				

*Dyno hip width starts at 14". Adding 2 ea. x 2" pads reduces hip width to 14"-4"= 10".*

#### Abductor/Medial Knee Block

*Must select hardware & one size of abductor pad. Wedge sizes are width x depth x height.*

<input type="checkbox"/>	3005-7001	E1028	Leggero Flip Down Medial Hdw.	£180.00
<input type="checkbox"/>	MTMD-12	E0957	Small Wedge Pad (3x5x3")	£96.00
<input type="checkbox"/>	MTMD-13	E0957	Medium Wedge Pad (5x5x3")	£96.00

#### EVA Positioning Wedges for Seat and/or Back

<input type="checkbox"/>	3001-7011	K0108	Wedges QTY 4 and Lumbar QTY 1	£38.00
--------------------------	-----------	-------	-------------------------------	--------

### Positioning Tray

<input type="checkbox"/>	3001-7014	E1028	Tray Hardware, Adj./Removable	£157.00
--------------------------	-----------	-------	-------------------------------	---------

<input type="checkbox"/>	3001-7013	E0950	DYNO Clear Positioning Tray	£112.00
--------------------------	-----------	-------	-----------------------------	---------

*3001-7014 Tray hardware is required*

### Ankle Supports

<input type="checkbox"/>	3001-6007	E0951	LEFT, Extra Small (Pad 2x8")	£68.00
<input type="checkbox"/>	3001-6008	E0951	RIGHT, Extra Small (Pad 2x8")	£68.00
<input type="checkbox"/>	3001-6009	E0951	LEFT, Small (Pad 3x10")	£68.00
<input type="checkbox"/>	3001-6010	E0951	RIGHT, Small (Pad 3x10")	£68.00

### Stealth Products Trunk & Pelvic Positioning

#### SP™ Anterior Posture Supports

<input type="checkbox"/>	IFPS1101-01	E0960	SP™ Classic, Dynamic, XS	£100.00
<input type="checkbox"/>	IFPS1102-01	E0960	SP™ Classic, Dynamic, S	£100.00
<input type="checkbox"/>	IFPS1103-01	E0960	SP™ Classic Dynamic, M	£100.00
<input type="checkbox"/>	IFPS1201-01	E0960	SP™ Classic, Structured, XS	£100.00
<input type="checkbox"/>	IFPS1202-01	E0960	SP™ Classic, Structured, S	£100.00
<input type="checkbox"/>	IFPS1203-01	E0960	SP™ Classic, Structured, M	£100.00
<input type="checkbox"/>	IFPS2101-01	E0960	SP™ Contoured, Dynamic, XS	£100.00
<input type="checkbox"/>	IFPS2102-01	E0960	SP™ Contoured, Dynamic, S	£100.00
<input type="checkbox"/>	IFPS2103-01	E0960	SP™ Contoured, Dynamic, M	£100.00
<input type="checkbox"/>	IFPS2201-01	E0960	SP™ Contoured, Structured, XS	£100.00
<input type="checkbox"/>	IFPS2202-01	E0960	SP™ Contoured, Structured, S	£100.00
<input type="checkbox"/>	IFPS2203-01	E0960	SP™ Contoured, Structured, M	£100.00

#### SP™ Pelvic Positioning Belts

<input type="checkbox"/>	IFPB5141	E0978	4- Point Push Button 1", XS	£90.00
<input type="checkbox"/>	IFPB5242	E0978	4- Point Push Button 1.5", SM	£90.00
<input type="checkbox"/>	IFPB5111	E0978	2-Point, Push Button 1", XS	£76.00
<input type="checkbox"/>	IFPB5112	E0978	2-Point, Push Button 1", SM	£76.00

*Side release available on request.*

#### SP™ Buckle Covers

<input type="checkbox"/>	IFMDBC10		Buckle Cover Kit 1" Buckle	£21.00
<input type="checkbox"/>	IFMDBC15		Buckle Cover Kit 1.5" Buckle	£21.00

#### Posture Support Strap Risers

<input type="checkbox"/>	3001-6011	K0108	Strap Risers, Pair	£33.00
--------------------------	-----------	-------	--------------------	--------

### Leggero Head Positioning

<input type="checkbox"/>	3001-6005	E0955	8" Oval Headrest and Hardware	£146.00
<input type="checkbox"/>	3001-6006	E0955	RESPOND™ Cradle Headrest and Hardware	£146.00
<input type="checkbox"/>	3005-6099		Omit Head Positioning	N/C

Hardware is included with 8" Multi-Adjustable & RESPOND™ Headrests

### Stealth Products Head Positioning

Select **A or B**. Then continue selections to right to complete the package.

#### A. Link Headrest Hardware 1" Ball

<input type="checkbox"/>	TWB480-1-01	E1028	Link Hardware 1" Ball	£146.00
--------------------------	-------------	-------	-----------------------	---------

#### B. Link Headrest Hardware 17mm Ball/Flip Down

<input type="checkbox"/>	TWB480-17-01	E1028	Link Hdw. 17mm Ball/Flip Down	£201.00
--------------------------	--------------	-------	-------------------------------	---------

#### Comfort Plus Headrest (Only if A is selected)

<input type="checkbox"/>	CP250-SR	E0955	Comfort Plus 10"	£96.00
<input type="checkbox"/>	CP450-SR	E0955	Comfort Plus 8"	£96.00

#### i2i Head Positioning System (Only if B is selected)

<input type="checkbox"/>	i2i250	K0108	i2i Low Pro Nino XS, Pad Only	£213.00
<input type="checkbox"/>	i2i350	K0108	i2i Low Pro Nino SM, Pad Only	£213.00
<input type="checkbox"/>	i2i450	K0108	i2i Low Pro Adult SM, Pad Only	£213.00

#### Chin Prompts (Only if i2i Head Positioning is selected)

<input type="checkbox"/>	i2i460	K0108	i2i SM Chin Prompt	£49.00
<input type="checkbox"/>	i2i560	K0108	i2i MED Chin Prompt	£49.00

#### Combo Series Headrest

<input type="checkbox"/>	CB315	E0955	Combo Series, Small	£185.00
<input type="checkbox"/>	CB305	E0955	Combo Series, Medium	£185.00
<input type="checkbox"/>	CB325	E0955	Combo Series, Large	£185.00